



Bross **MEDICAL**




**Augeninstrumente
Ophthalmic Instruments
Oftalmik aletleri
instrumentos oftálmicos**



instruments made in Germany



Your partner for high class Ophthalmic instruments



Note

Dear interested customers,

we are delighted to introduce our Catalogue about our range of ophthalmic instruments.

Change of design and pattern reserved. © 2015 BROSS Medical GmbH
- All rights reserved! All rights of the manufacturer and the owner of this documentation are reserved. Copies of this documentation, including extracted forms, require the approval of BROSS Medical GmbH.

All pictures may vary from actual size. Descriptive, photographic and typographic errors are subject to correction.
We reserved the right to make alterations.

Bross Medical GmbH

-Headquarter-

Albstrasse 12
78532 Tuttlingen
GERMANY

-Administration-

Dreifaltigkeitsbergstrasse 12
78588 Denkingen
GERMANY

Phone: +49 (0) 7424 7037330
Fax: +49 (0) 7424 7037339
Mail: Ben@bross-medical.com
simone@bross-medical.com



Über uns

BROSS Medical GmbH ist spezialisiert auf die Herstellung und den Vertrieb Chirurgischer Instrumente für die Ophthalmologie, Orthopädie und Allgemein Chirurgie.

Hier bieten wir Ihnen viele Möglichkeiten wie z.B.: Instrumente „100% MADE IN GERMANY“ sowie dieselben Instrumente mit einem sehr guten Preis/Leistungsverhältnis welche nicht in Deutschland gefertigt werden. Diese werden aber nur von ausgesuchten Herstellern gefertigt welche nach Deutschen Standards produzieren und wir von BROSS Medical machen stets eine 100% Warenkontrolle so das Sie sicher sein können kein schlechtes Instrument verlässt unser Haus.

Wer steht hinter BROSS Medical

Benjamin Heiss (CEO) verfügt über mehr als 10 Jahre Erfahrung in der Medizintechnik. Zu Beginn beschäftigte er sich mit der Herstellung von Produkten für die Endoskopie wechselte aber nach seiner Ausbildung in das Familienunternehmen wo er sich für die nächsten 3.Jahre mit der Herstellung von Mikrochirurgischen Instrumenten auseinandersetzte. Danach wechselte er in den Vertrieb und das Marketing beendete eine Kaufmännische Weiterbildung sowie das Studium zum Betriebswirt GA. Hier betreute er ein Weltweites Netzwerk von Fachhändlern.

Simone Nadine Selim (CEO) verfügt über mehr als 15 Jahren Erfahrung in Finanzbuchhaltung. Nach Ihrer Ausbildung als Kauffrau im Groß- und Aussenhandel bei einem internationalem Grossunternehmen erweiterte sie ihr Spektrum in den Bereichen Einkauf und Export.

About Us

BROSS Medical GmbH is a certified manufacturer and distributor of surgical instruments. Our specialties include ophthalmology, orthopedics and general surgery.

Bross offers many options, such as: Instruments that are“100% MADE IN GERMANY“ and we offer the same instruments with an excellent price / performance ratio which are not manufactured in Germany but are made from high quality forgings (no Pakistan forgings). These are made by select manufacturers that produce to German standards. Most importantly at BROSS Medical we inspect to a 100% physical check so that you can be sure only BROSS high quality instruments leave our facility and arrive at your door.

Who is behind BROSS Medical:

Benjamin Heiss (CEO) has over 10 years experience in the medical industry. He Initially worked on the production of products for endoscopic surgery. However he switched after training in the family business to the production of micro-surgical instruments for the next 3 years. Next moving into an executive position in sales and marketing. There he supervised a worldwide network of dealers and ended up his commercial education by obtaining a degree in business administration, GA. Making Benjamin well prepared to manage and open BROSS for international business.

Simone Nadine Selim (CEO) has over 15 years experience in financial accounting. She also has extensive training as a management assistant in a large international firm specializing in wholesale and foreign trade. Since then Simone has expanded her expertise into the areas of purchasing and exporting.

Zertifikate

Zertifizierung nach DIN ISO 13485 und FDA

Vertrieb

BROSS MEDICAL GmbH
Bitte kontaktieren Sie uns wenn Sie Interesse
haben BROSS Medical GmbH Produkte in Ihrem
Land zu vertreiben!

Vielen Dank


Certificates

Certification for ISO 13485 and FDA

Distribution

BROSS MEDICAL GmbH
Please, contact us if you are interested in Distributing
BROSS Medical GmbH Products in your country

Thank you very much!





Cannula Sizes

Size in Gauge (G)	18	19	20	21	22	23	24	25	26	27	28	29	30
Outside Dimensions (mm) ISO/DIN 9626	1,2	1,1	0,9	0,8	0,7	0,6	0,55	0,5	0,45	0,4	0,36	0,33	0,3

Thread Sizes

Metric	USP	Diameter in mm
0,01	12 – 0	0,001 – 0,009
0,1	11 – 0	0,010 – 0,019
0,2	10 – 0	0,020 – 0,029
0,3	9 – 0	0,030 – 0,039
0,4	8 – 0	0,040 – 0,049
0,5	7 – 0	0,050 – 0,069
0,7	6 – 0	0,070 – 0,099
1	5 – 0	0,100 – 0,149
1,5	4 – 0	0,150 – 0,199
2	3 – 0	0,200 – 0,249
2,5 – 3	2 – 0	0,250 – 0,349
3,5	0	0,350 – 0,399
4	1	0,400 – 0,499
5	2	0,500 – 0,599
6	3	0,600 – 0,699
7	5	0,700 – 0,799
8	6	0,800 – 0,899
9	7	0,900 – 0,999

„BLACK-SERIES“

What does „BLACK-SERIES“ means ?




Materials we use



xx-xxx-01

Standard Material: STAINLESS STEEL, SILVER, PLATINUM



xx-xxx-02

TITANIUM



xx-xxx-03

BLACK SERIES




xx-xxx-06

Tungsten Carbide (TC)



xx-xxx-07

Super Cut Scissors





ZERTIFIKAT ◆ CERTIFICATE ◆ 認証證書 ◆ CERTIFICADO ◆ CERTIFICAT



CERTIFICATE

The certification body of Swiss TS Technical Services AG hereby confirms that the company

Bross Medical GmbH
D-78532 Tuttlingen



has introduced and applies a management system for

Manufacturing and distribution of surgical instruments

according to:

ISO 13485:2003 Quality management for medical devices

Registration number: 12-336-634
Initial certification: 30.11.2012
Valid until: 29.11.2015

Heinrich A. Bieler
Head of the certification body
Wallisellen, 30.11.2012
The certification body of Swiss TS Technical Services AG
A SVTI and TUV SUD company



SCESm013

ZERTIFIKAT ◆ CERTIFICATE ◆ 認証證書 ◆ CERTIFICADO ◆ CERTIFICAT



CE-Certificate

Full quality assurance system

The conformity assessment body for medical devices of Swiss TS Technical Services AG hereby certifies that:

Bross Medical GmbH
Albstrasse 12
D-78532 Tuttlingen



has introduced and applies a full quality assurance system according to annex II.3 of directive 93/42 EEC.

for the products

Monopolar and bipolar forceps

UMDNS: 11-502, rule 9, class IIb

The full quality assurance system covers design, manufacturing and final inspection of the products according to the requirements of annex II.3 of the directive 93/42 EEC. The quality assurance system is regularly controlled by Swiss TS Technical Services AG according to annex II.5.

Registration number: 12-380-718
Initial certification: 30.11.2012
Valid until: 29.11.2015

CE marking of conformity: **CE1253**

Heinrich A. Bieler
Head of the certification body

Wallisellen, 30.11.2012
Notified Body 1253
Swiss TS Technical Services AG
Richtstrasse 15
CH-8304 Wallisellen
Swiss TS, a company of SVTI and TUV SUD




SCESm013

CONTENT

KAT	English	Deutsch	Page
1.	Lid Speculums	Lidsperrer	10.-21.
2.	Diamond Knives	Diamantmesser	22.-27.
3.	Knives	Messer	28.-34.
4.	Gauge Instruments	Messinstrumente	35.-37.
5.	Hooks	Haken	38.-48.
6.	Spatulas & Curettes	Spatel & Löffel	49.-56.
7.	Probes & Trephines	Sonden & Trepane	57.-64.
8.	Markers	Markeure	65.-69.
9.	Cannulas	Kanülen	70.-85.
10.	Forceps	Pinzetten	86.-111.
11.	Scissors	Scheren	112.-128.
12.	Needle Holders	Nadelhalter	129.-133.
13.	Vitrectomy	Vitrektomie	134.-138.
14.	Miscellaneous Products	Sonstige Produkte	139.-147.
15.	Electro Surgery	Elektrochirurgie	148.-149.
95.	Sterilisation	Sterilisation	150.-163.
99.	Instruments Sets	Operationsbestecke	164.-184.


Barraquer-Colibri Eye Speculum
Barraquer-Colibri Lidsperrer

- 01-01-100-1** 5 mm - premature birth
- 01-01-100-3** BLACK SERIES
- 01-01-102-1** 10 mm - child
- 01-01-102-3** BLACK SERIES
- 01-01-104-1** 15 mm - adult
- 01-01-104-3** BLACK SERIES




Barraquer Eye Speculum sliding, 15mm blade
Barraquer Lidsperrer beweglich, 15mm Valve

- 01-01-106-1** open
- 01-01-107-1** closed




Thornton Eye Speculum 14 mm
Thornton Lidsperrer 14 mm

- 01-01-110-1** adult
- 01-01-110-3** BLACK SERIES




Kratz-Barraquer Eye Speculum
Kratz-Barraquer Lidsperrer

- 01-01-120-1** 10 mm - child
- 01-01-120-3** BLACK SERIES
- 01-01-122-1** 14 mm - adult
- 01-01-122-2** TITANIUM
- 01-01-122-3** BLACK SERIES



Barraquer Phaco Eye Speculum
Barraquer Phaco Lidsperrer

- 01-01-130-1** with irrigation
- 01-01-132-1** standard wire
- 01-01-132-2** TITANIUM
- 01-01-132-3** BLACK SERIES
- 01-01-134-1** strong wire - Ø 1 mm
- 01-01-134-3** BLACK SERIES




Bowman Eye Speculum 7 cm
Bowman Lidsperrer 7 cm

01-01-145-1



Barraquer Eye Speculum
Barraquer Lidsperrer

01-01-160-1 10 mm - solid blades child
01-01-160-3 BLACK SERIES
01-01-162-1 15 mm - solid blades adult
01-01-162-2 TITANIUM
01-01-162-3 BLACK SERIES




Nasal Temporal Approach Eye Speculum,
15 mm solid blades adult with spring on nasal side
Nasal Temporal Approach Lidsperrer
15mm geschlossene Valven Erwachsene

01-01-170-1
01-01-170-3 - BLACK SERIES



Nasal Temporal Approach Eye Speculum,
15 mm open wire adult with spring on nasal side
Nasal Temporal Approach Lidsperrer,
15mm Erwachsene

01-01-172-1
01-01-172-2 - TITANIUM
01-01-172-3 - BLACK SERIES




Pierse Eye Speculum,
3.5cm, 11mm valve
Pierse Lidsperrer,
3.5cm, 11mm valve

01-01-175-1
01-01-175-3 - BLACK SERIES



Corcelle Eye Speculum
Corcelle Lidsperrer

- 01-01-180-1** 8 mm - solid blades
- 01-01-180-3** BLACK SERIES
- 01-01-182-1** 12 mm - solid blades
- 01-01-182-3** BLACK SERIES



Lancaster Eye Speculum
Lancaster Lidsperrer

- 01-01-185-1** 8cm, solid blades




Lester-Burch Eye Speculum 8cm, solid blades
Lester-Burch Lidsperrer, 8cm, geschlossene Valven

- 01-01-190-1**



Williams Eye Speculum
Williams Lidsperrer

- 01-01-195-1** 11 mm - child
- 01-01-197-1** 14 mm - adult



Knapp-Clark Eye Speculum
Knapp-Clark Lidsperrer

01-01-200-1 16 mm - medium size



Castroviejo Eye Speculum
Castroviejo Lidsperrer

01-01-210-1 16 mm - medium size

01-01-212-1 18 mm - large size



Castroviejo Eye Speculum
Castroviejo Lidsperrer

01-01-214-1 new Version



Sauer Eye Speculum
Sauer Lidsperrer

01-01-220-1 3 mm - solid blades premature birth


01-01-222-1 8 mm - solid blades infant



Maumenee-Park Eye Speculum
Maumenee-Park Lidsperrer

01-01-240-1 fenestrated

01-01-242-1 solid blades




Pretori Eye Speculum 13,5mm
Pretori Lidsperrer 13,5mm

- 01-01-245-1** left
- 01-01-246-1** right




Murdoch Eye Speculum
Murdoch Lidsperrer

- 01-01-250-1** left - child size
- 01-01-251-1** left - child size - with locking screw
- 01-01-252-1** left - adult size
- 01-01-253-1** left - adult size - with locking screw
- 01-01-254-1** right - child size
- 01-01-256-1** right - adult size



Cook Eye Speculum solid blades
Cook Lidsperrer geschlossene Valven

- 01-01-260-1** 5 mm - infant
- 01-01-262-1** 8 mm - child - small
- 01-01-266-1** 11 mm - child - medium
- 01-01-268-1** 15 mm - adult- large



Mellinger Eye Speculum 15mm
Mellinger Lidsperrer 15mm

- 01-01-270-1** fenestrated adult
- 01-01-272-1** solid blades adult




Schott Eye Speculum
Schott Lidsperrer

- 01-01-280-1** 12 mm fenestrated child
- 01-01-282-1** 15 mm fenestrated adult



Bangerter Eye Speculum
Bangerter Lidsperrer

- 01-01-290-1** 11 mm - solid blades child
- 01-01-292-1** 15 mm - solid blades adult




Lieberman Brown Triple Eye Speculum
Lieberman Brown Triple Lidsperrer

- 01-01-294-1** 14 mm Valve
- 01-01-294-3** BLACK SERIES



Lieberman Eye Speculum O-Wire
Lieberman Lidsperrer O-Wire

- 01-01-298-1** 5 mm
- 01-01-298-3** BLACK SERIES
- 01-01-299-1** 15 mm
- 01-01-299-3** BLACK SERIES



Lieberman Eye Speculum K-Wire
Lieberman Lidsperrer K-Wire

- 01-01-300-1** 10 mm
- 01-01-300-3** BLACK SERIES
- 01-01-302-1** 15 mm
- 01-01-302-2** TITANIUM
- 01-01-302-3** BLACK SERIES




Lieberman Eye Speculum V-Wire
Lieberman Lidsperrer V-Wire

- 01-01-304-1** 10 mm
- 01-01-304-3** BLACK SERIES
- 01-01-306-1** 15 mm
- 01-01-306-2** TITANIUM
- 01-01-306-3** BLACK SERIES



Ginsberg Eye Speculum flared upper blades and fenestrated lower blades

01-01-308-1
01-01-308-3 BLACK SERIES




Lieberman Eye Speculum 15 mm Aspirating
Lieberman Lidsperrer, 15 mm Aspirating

01-01-310-1 K-Wire
01-01-311-1 V-Wire



Lieberman Eye Speculum
Lieberman Lidsperrer

01-01-312-1 10 mm - solid blades - child
01-01-312-3 BLACK SERIES
01-01-314-1 15 mm - solid blades - adult
01-01-314-3 BLACK SERIES



Kershner Eye Speculum nasal/temporal
Kershner Lidsperrer nasal/temporal

01-01-316-1
01-01-316-3 BLACK SERIES



Lieberman Eye Speculum nasal/temporal 15 mm
Lieberman Lidsperrer nasal/temporal 15 mm

01-01-318-1 solid blades
01-01-318-3 BLACK SERIES




Thorlakson Eye Speculum
Thorlakson Lidsperrer

01-01-330-1
01-01-330-2 TITANIUM
01-01-330-3 BLACK SERIES



Weiss Eye Speculum
Weiss Lidsperrer

01-01-340-1 14 mm - child
01-01-342-1 17 mm - adult



Weiss Eye Speculum
Weiss Lidsperrer

01-01-344-1 14 mm - solid blades - child
01-01-346-1 17 mm - solid blades - adult



Stevenson Lacrimal Sac Retractor
Stevenson Tränen­sackspreizer

01-01-350-1




Stevenson Lacrimal Sac Retractor
Stevenson Tränen­sackspreizer

01-01-352-1




Agricola Lacrimal Sac Retractor
Agricola Tränen­sackspreizer

01-01-355-1



Goldstein Lacrimal Sac Retractor
Goldstein Tränen­sackspreizer

01-01-357-1 3x3 teeth, opening wide 5mm
01-01-358-1 3x3 teeth, opening wide 17mm



Desmarres Lid Retractor
Desmarres Lidhalter

01-01-500-1 FH 8 mm
01-01-500-3 BLACK SERIES
01-01-502-1 FH 10 mm
01-01-502-3 BLACK SERIES
01-01-504-1 FH 12 mm
01-01-504-3 BLACK SERIES
01-01-506-1 FH 14 mm
01-01-506-3 BLACK SERIES
01-01-508-1 FH 16 mm
01-01-508-3 BLACK SERIES
01-01-510-1 FH 18 mm
01-01-510-3 BLACK SERIES



Jaffe Lid Retractor
Jaffe Lidhalter

- 01-01-520-1** 15 mm (pair)
01-01-522-1 10 mm (pair)

**Jaeger** Lid Plate 11 cm
Jaeger Lidplatte 11 cm

- 01-01-530-1**
01-01-530-2 TITANIUM
01-01-530-3 BLACK SERIES

**Helveston** Tissue Retractor
Helveston Spatel

- 01-01-550-1** Barbie- 7 mm - Steel
01-01-550-2 TITANIUM
01-01-550-3 BLACK SERIES
01-01-552-1 Big Barbie- 9 mm - Steel
01-01-552-2 TITANIUM
01-01-552-3 BLACK SERIES
01-01-554-1 Great Big Barbie- 11 mm - Steel
01-01-554-2 TITANIUM
01-01-554-3 BLACK SERIES

**Bangerer** Spatula
Bangerer Spatel

- 01-01-560-1** 6 mm
01-01-560-3 BLACK SERIES
01-01-562-1 10 mm
01-01-562-3 BLACK SERIES
01-01-564-1 12 mm
01-01-564-3 BLACK SERIES

**Bonn** Conjunctival and Orbital Retractor
Bonn Conjunctival and Orbital Retractor

- 01-01-570-1** used in retinal detachment
01-01-570-3 BLACK SERIES



Sautter Orbital Spatula
Sautter Orbitalspatel

- 01-01-573-1** small
- 01-01-573-3** BLACK SERIES
- 01-01-575-1** medium
- 01-01-575-3** BLACK SERIES
- 01-01-577-1** large
- 01-01-577-3** BLACK SERIES



Schepens Forked Orbital Retractor
Schepens Orbital Retraktor gegabelt

- 01-01-580-1**
- 01-01-580-3** BLACK SERIES




Gass Scleral Marker Ø 1,5 mm
Gass Sklera Marker Ø 1,5 mm

- 01-01-590-1**



Schepens Scleral Depressor
Schepens Sklera Depressor

- 01-01-592-1** #8 small
- 01-01-594-1** #10 medium
- 01-01-596-1** #12 large



Wilder Scleral Depressor pocket clip
Wilder Sklera Depressor mit Taschenclip

01-01-598-1



Schocket Scleral Depressor pocket clip
Schocket Sklera Depressor mit Taschenclip 2endig

01-01-600-1

01-01-600-2 TITANIUM



Lasik Flap Smoother
Lasik Glättstreicher für Lamellar Klappe

01-01-605-1




Caro/Avalos Lasik Flap Smoother
Caro/Avalos Glätter für Lamellar Klappe


01-01-610-1




Easy Clean Handle




Diamond Knife straight Blade 20° 0.5 mm
Diamantmesser gerade Klinge 20° 0.5mm
For Deep Sclerectomy
01-02-160-1




Diamond Knife straight Lance 40° Blade 1.0 mm
Diamantmesser gerade Lanze 40° Klinge 1.0mm
For Side Port Incisions
01-02-220-1




Diamond Knife straight Blade 30° 1.0 mm
Diamantmesser gerade Klinge 30° 1.0mm
For Side Port Incisions
01-02-170-1




Diamond Knife straight Lance 60°
Diamantmesser gerade Lanze 60°
For Side Port Incisions
01-02-230-1 Blade 1.0 mm
01-02-240-1 Blade 1.4 mm




Diamond Knife straight Blade 45° 1.0 mm
Diamantmesser gerade Klinge 45° 1.0mm
For Side Port Incisions
01-02-180-1



Diamond Knife straight Spear 90°
Diamantmesser gerade Speer 90°
For Side Port Incisions
01-02-250-1 Blade 1.0 mm
01-02-260-1 Blade 1.4 mm




Diamond Knife straight Tri-facet Blade 1.0 mm
Diamantmesser gerade Tri-facet Klinge 1.0mm
For Side Port Incisions
01-02-190-1




Diamond Knife straight TriSide Blade 0.55-1.0 mm
Diamantmesser gerade TriSide Klinge 0.55-1.0mm
For Side Port Incisions
01-02-280-1


Easy Clean Handle




Diamond Knife angled TriSide Blade 0.55-1.0 mm
Diamantmesser abgewinkelt TriSide Klinge 0.55-1.0mm
For Scleral Tunnel Incisions
For Opening of the Anterior Chamber
01-02-300-1




Diamond Knife angled Crescent Blade 1.4 mm
Diamantmesser abgewinkelt Tunnel Klinge 1.4mm
For Scleral Tunnel Incisions
01-02-310-1




Diamond Knife angled Crescent Blade 2.0 mm
Diamantmesser abgewinkelt Tunnel Klinge 2.0mm
For Scleral Tunnel Incisions
01-02-320-1




Diamond Knife angled Lance 40° Blade 1.0 mm
Diamantmesser abgewinkelt Lanze 40° Klinge 1.0mm
For Side Port Incisions
01-02-360-1




Diamond Knife angled Lance 60°
Diamantmesser abgewinkelt Lanze 60°
For Side Port Incisions
01-02-370-1 Blade 1.0 mm
01-02-380-1 Blade 1.4 mm



Diamond Knife angled Phaco Spear 90°
Diamantmesser abgewinkelt Speer 90°
For Clear Cornea Incisions
01-02-410-1 Blade 1.00 mm
01-02-420-1 Blade 1.40 mm
01-02-430-1 Blade 2.50 mm
01-02-450-1 Blade 2.80 mm
01-02-460-1 Blade 3.00 mm
01-02-470-1 Blade 3.20 mm




Diamond Knife angled Clear Cornea 90°
Diamantmesser abgewinkelt Clear Cornea 90°
For Clear Cornea Incisions
01-02-500-1 Blade 2.50 mm
01-02-510-1 Blade 2.65 mm
01-02-520-1 Blade 2.80 mm
01-02-530-1 Blade 3.00 mm
01-02-540-1 Blade 3.20 mm




Diamond Knife angled Trapezoid 90°
Diamantmesser abgewinkelt Trapezoid 90°
For Clear Cornea Incisions
01-02-580-1 Blade 2.5/3.0 mm
01-02-590-1 Blade 2.7/2.9 mm
01-02-595-1 Blade 2.8/3.2 mm
01-02-600-1 Blade 3.0/3.2 mm

Easy Clean Handle




Diamond Knife Pre-Set Lance 60°
Diamantmesser tiefenfixiert Lanze 60°
For Limbal Relaxation Incisions LRI
For Initial Groove Incisions

- 01-02-620-1** Cutting-depth 0.5 mm
- 01-02-630-1** Cutting-depth 0.55 mm
- 01-02-640-1** Cutting-depth 0.6 mm




Diamond Knife Pre-Set Truncated Lance 60°
Diamantmesser tiefenfixiert Truncated Lanze 60°
For Limbal Relaxation Incision LRI
For Initial Groove Incisions


- 01-02-660-1** Cutting-depth 0.4 mm
- 01-02-670-1** Cutting-depth 0.5 mm
- 01-02-680-1** Cutting-depth 0.6 mm




Diamond Knife Micrometer Arcuate Blade 15° 0.2 / 1.0 mm
Diamantmesser Mikrometer Arcuate Klinge 15° 0.2-1.0mm
For Astigmatic Surgery
01-02-690-1




Diamond Knife Micrometer TriSquare Blade 1.0 mm
Diamantmesser Mikrometer TriSquare Klinge 1.0mm
For Astigmatic Surgery
01-02-700-1




Diamond Knife Micrometer Blade 45° 1.0 mm
Diamantmesser Mikrometer Klinge 45° 1.0mm
For Refractive Surgery
01-02-710-1




Diamond Knife Micrometer Axe 45° Blade 1.6 mm
Diamantmesser Mikrometer Axe 45° Klinge 1.6mm
For Refractive Surgery
01-02-720-1



Diamond Knife Micrometer DoubleEdge 45° Blade 1.0 mm
Diamantmesser Mikrometer DoubleEdge 45° Klinge 1.0mm
For Refractive Surgery
01-02-730-1




Diamond Knife Micrometer Trak 45° Facet 0.25mm Blade 1.0 mm
Diamantmesser Mikrometer Trak 45° Facet 0.25mm Klinge 1.0mm
For Refractive Surgery
01-02-740-1




Diamond Knife Micrometer SingleEdgeVertical 45° Blade 1.0 mm
Diamantmesser Mikrometer SingleEdgeVertical 45° Klinge 1.0mm
For Refractive Surgery
01-02-750-1




7-Steps 0/0.3/0.4/0.5/0.6/0.8/1.0/2.5 mm




Blade 45° 1.0 mm, For Limbal Relaxation Incisions LRI
Diamantmesser 7-Stufen Klinge 45° 1.0mm
For Initial Groove Incisions
01-02-760-1




Tri-facet Blade 1.0 mm
Diamantmesser Tri-facet Klinge 1.0mm
For Limbal Relaxation Incisions LRI
For Initial Groove Incisions
01-02-770-1




DoubleEdge 45° Blade 1.0 mm
Diamantmesser DoubleEdge 45° Klinge 1.0mm
For Refractive Surgery
01-02-780-1




Lance 40° MicroFacet Blade 0.15-1.0 mm
Diamantmesser Lanze 40° MikroFacet Klinge 0.15-1.0mm
For Limbal Relaxation Incisions LRI
For Initial Groove Incisions
01-02-790-1




Phaco Spear 90°
Diamantmesser Phako Speer 90°
For Refractive Surgery
01-02-830-1 Blade 3.00 mm
01-02-840-1 Blade 3.20 mm



Clear Cornea 90°
Diamantmesser Clear Cornea 90°
01-02-850-1 Blade 2.50 mm
01-02-860-1 Blade 2.65 mm
01-02-870-1 Blade 2.80 mm
01-02-880-1 Blade 3.00 mm
01-02-890-1 Blade 3.20 mm




Diamond Knife Retina straight Blade 20° 0.5 mm
Diamantmesser Retina gerade Klinge 20° 0.5mm
For Cutting Of Subretinal Membranes
01-02-900-1



Diamond Knife Retina straight Blade 45° 0.9 mm
Diamantmesser Retina gerade Klinge 45° 0.9mm
For Cutting Of Subretinal Membranes
01-02-910-1

Scalpell Handle No.3
Skalpellgriff Nr.3

01-03-050-1



Scalpell Blades
Skalpellklingen
each 100pcs



01-03-051-1



01-03-052-1



01-03-053-1



01-03-054-1



01-03-055-1



01-03-056-1



01-03-057-1



01-03-058-1



01-03-059-1

Blade Handle 10 cm Ø 5 mm
Klingengriff 10 cm Ø 5 mm

01-03-100-1
01-03-100-2 TITANIUM



Scalpell Blades
Skalpellklingen
each 10pcs



01-03-101-1



01-03-102-1



01-03-103-1



01-03-104-1



01-03-105-1



01-03-106-1

Castroviejo Bladeholder 13 cm Ø 8 mm
Castroviejo Klingenthaler 13 cm Ø 8 mm

01-03-110-1

01-03-110-2 TITANIUM



Troutman-Chris Bladeholder 10 cm Ø 5.5 mm
Troutman-Chris Klingenthaler 10 cm Ø 5.5 mm

01-03-120-1

01-03-120-2 TITANIUM




Swiss Model Bladeholder 12 cm straight
Swiss Model Klingenthaler 12 cm gerade

01-03-150-2 TITANIUM




Barraquer Bladeholder 11 cm
Barraquer Klingenthaler 11 cm

01-03-160-1



RB Razor Blades breakable
Rasierklingen, brechbar

01-03-170-1 Package of 10x1 piece non-sterile




Graefe Cataract Knife
Gräfe Kataraktmesser

01-03-180-1 Blade 20 mm
01-03-190-1 Blade 25 mm




Ziegler Discission Knife
Ziegler Diszisionmesser

01-03-200-1 Blade 3 mm
01-03-210-1 Blade 5 mm




Kuhnt Corneal Scarifier
Kuhnt Hornhautschaber

01-03-230-1 Blade 3 x 5 mm



Bonn Scleral Knife
Bonn Skleramesser

01-03-240-1 Blade 5 mm
01-03-250-1 Blade 6 mm



Gill Corneal Knife curved
Gill Hornhautmesser, gebogen

01-03-260-1 Blade 3 mm



BROSS Micro Round Knife
for Visco Canalostomy

01-03-265-1 Blade 1 mm



Tooke Corneal Knife
Tooke Hornhautmesser

01-03-280-1 Blade 3 x 18 mm



Dardenne Sato Knife
Dardenne Satomesser

01-03-290-1 Blade 1.5 mm



Barkan Goniotomy Knife
Barkan Goniotomiemesser

01-03-310-1 Blade 2 mm



Wheeler Discission Knife
Wheeler Diszisionsmesser

01-03-320-1 Blade 1 x 18 mm



Mod.Hamburg Keratome
Mod.Hamburg Luftlanze

01-03-330-1



Golf-Club Foreign Body
Golf-Club Fremdkörpernadel

01-03-390-1 Spud, sharp edge




Davis Foreign Body spud
Davis Fremdkörpermeißel

01-03-410-1 sharp round





Ellis Foreign Body Spud
Ellis Fremdkörpermeißel

01-03-420-1 curved rounded blades





FNS Foreign Body Instrument
FNS Fremdkörperinstrument

01-03-440-1 Needle and Spud reversible




FNS Foreign Body Instrument
FNS Fremdkörperinstrument

01-03-450-1 Needle and Gouge reversible



FN Foreign Body Needle reversible
FN Fremdkörpernadel

01-03-460-1 without handle



FS Foreign Body Spud reversible
FS Fremdkörperflachmeißel

01-03-470-1 without handle




FS Foreign Body Gouge reversible
FS Fremdkörperhohlmeißel

01-03-480-1 without handle




Graefe Cystotome
Graefe Zystotom

01-03-490-1 malleable



LS Stiletto
LS Stilett

01-03-520-1 1.2 mm
01-03-520-2 1.2 mm **TITANIUM**



Wright Fascia Needle

01-03-600-1




Castroviejo Caliper
Castroviejo Zirkel


- 01-04-100-1** 0-15 mm straight
- 01-04-102-1** 0-20 mm straight
- 01-04-104-1** 0-20 mm Increments 0.5 mm - Steel
- 01-04-104-2** 0-20 mm Increments 0.5 mm - **TIT**
- 01-04-106-1** 0-20 mm curved

**Maloney** Intra-Operative Keratometer
Maloney Intra-Operativer Keratometer

- 01-04-110-2** **TITANIUM**

**Mendez** Degree Gauge C 10° from 0 - 180°
Mendez Grad Messinstrument K 10° von 0 - 180°

- 01-04-120-1**
- 01-04-122-1** Lu-Mendez, front part open


**Osher** Incision Caliper 1 - 10 mm
Osher Inzisions-Zirkel 1-10mm

- 01-04-130-1**




DK Caliper Block in 0.02 mm increments
DK Messblock Teilung 0.02 mm

01-04-140-1



Stahl Caliper for measuring Corneal Ø w/w
Stahl Kaliper z. Messung Cornea ø w/w

01-04-150-1



Helveston Scleral Marking Ruler
Helveston Scleral Marking Ruler

01-04-160-1 Steel
01-04-160-2 TITANIUM



Gonnin Amsler Marker 11 - 23 mm / 1 mm
Gonnin Amsler Marqueur 11-23 mm / 1 mm

01-04-170-2 TITANIUM




MH Ruler 150 mm inches and mm flexible
MH Meßstab 150 mm/6 inch in 1 mm Teilung flexibel

01-04-180-1



Jameson Caliper 0 - 80 mm / 0 - 3 inch
Jameson Meßschieber 0-80 mm/0-3 inch

01-04-190-1



Schioetz Tonometer with inclined scale
Schioetz Tonometer mit schräger Skala

01-04-200-1 weights 10gr and 7gr



Jameson Muscle Hook
Jameson Muskelhaken

- 01-05-110-1** large
- 01-05-110-3** BLACK SERIES
- 01-05-115-1** small
- 01-05-115-3** BLACK SERIES



Green Muscle Hook
Green Muskelhaken

- 01-05-120-1** Flat Handle
- 01-05-120-3** BLACK SERIES



Graefe Strabismus Hook
Graefe Schielhaken

- 01-05-130-1** #1 small
- 01-05-130-3** BLACK SERIES
- 01-05-135-1** #2 medium
- 01-05-135-3** BLACK SERIES
- 01-05-145-1** #3 large
- 01-05-145-3** BLACK SERIES



Gass Retinal Detachment Hook with hole
Gass Netzhautablösungshaken mit Loch

- 01-05-150-1**
- 01-05-150-3** BLACK SERIES



Stevens Tenotomy Hook
Stevens Sehnenhäkchen

01-05-160-1
01-05-160-3 BLACK SERIES



Helveston -Teaser- Hook 6 mm extra delicate
Helveston Schielhäkchen sehr fein

01-05-165-1
01-05-165-2 TITANIUM
01-05-165-3 BLACK SERIES



Helveston Finder Hook - extra delicate
Helveston Finder Haken - sehr fein

01-05-170-1 large - Steel
01-05-170-2 TITANIUM
01-05-170-3 BLACK SERIES
01-05-175-1 small - Steel
01-05-175-2 TITANIUM
01-05-175-3 BLACK SERIES



Helveston Muscle Hook - extra delicate
Helveston Muskelhaken - sehr fein

01-05-180-1 8 mm - Steel
01-05-180-2 TITANIUM
01-05-180-3 BLACK SERIES
01-05-185-1 10 mm - Steel
01-05-185-2 TITANIUM
01-05-185-3 BLACK SERIES
01-05-190-1 12 mm - Steel
01-05-190-2 TITANIUM
01-05-190-3 BLACK SERIES



Tyrell Iris Hook sharp
Tyrell Irishäkchen scharf

01-05-210-1
01-05-210-3 BLACK SERIES




Jaffe-Knolle Iris Hook blunt angled
Jaffe-Knolle Irishäkchen stumpf abgewinkelt

01-05-245-1
01-05-245-2 TITANIUM
01-05-245-3 BLACK SERIES



Crawford Retrieving Hook blunt 15 cm
Crawford Rückziehhaken stumpf 15cm

01-05-250-1
01-05-250-3 BLACK SERIES



Crawford Lacrimal Intubation Set
Crawford Tränenweg Set

01-05-251-1



Crawford Retrieving Hook for Lacrimal
Intubation probe

01-05-252-1
01-05-252-3 BLACK SERIES




Hirschmann Iris Hook H-tip angled 45°
Hirschmann Irishäkchen abgewinkelt 45°

01-05-255-1
01-05-255-3 BLACK SERIES




Guthrie Double-Fixation Hook
Guthrie Doppel-Fixierhaken

01-05-260-1 small
01-05-270-1 large




FC Sclera Twist Fixation Claw
FC Fixierkralle

01-05-280-1 right
01-05-285-1 left




Knapp Lacrimal Retractor sharp 4 prongs
Knapp Wundhaken scharf 4 Zähne

01-05-290-1




Knapp Lacrimal Retractor blunt 4 prongs
Knapp Wundhaken stumpf 4 Zähne

01-05-300-1




Bechert Lens Nucleus Rotator angled
Bechert Linsenkernelwender abgewinkelt

- 01-05-315-1 7 mm - Stainless Steel
- 01-05-315-2 TITANIUM
- 01-05-315-3 BLACK SERIES
- 01-05-320-1 9 mm - Stainless Steel
- 01-05-320-2 TITANIUM
- 01-05-320-3 BLACK SERIES




Clayman Lens Nucleus Rotator V-Tip
Clayman Linsenkernelwender V-Spitze

- 01-05-325-2 TITANIUM
- 01-05-325-3 BLACK SERIES



Drysdale Lens Nucleus Rotator angled
Drysdale Linsenkernelmanipulator abgewinkelt

- 01-05-330-1 Stainless Steel
- 01-05-330-2 TITANIUM
- 01-05-330-3 BLACK SERIES




Tennant Lens Nucleus Manipulator ball 0.5 mm
Tennant Linsenkernelmanipulator Kugel 0.5 mm

- 01-05-335-2 TITANIUM
- 01-05-335-3 BLACK SERIES



Osher Lens Nucleus Manipulator straight
Osher Linsenkernmanipulator gerade

01-05-340-2 TITANIUM
01-05-340-3 BLACK SERIES



Neuhann Nucleus Dividing Hook angled
Neuhann Kernspalter abgewinkelt

01-05-345-1



Nagahara Phaco Chopper blunt edge
Nagahara Phako Chopper Klinge stumpfes Eck

01-05-350-1 90° - Stainless Steel
01-05-350-2 TITANIUM
01-05-350-3 BLACK SERIES
01-05-355-1 60° - Stainless Steel
01-05-355-2 TITANIUM
01-05-355-3 BLACK SERIES



Rosen Phaco Splitter Blade blunt edge
Rosen Phaco Splitter Blade blunt edge

01-05-365-1 90°
01-05-365-2 TITANIUM
01-05-365-3 BLACK SERIES
01-05-370-1 60°
01-05-370-2 TITANIUM
01-05-370-3 BLACK SERIES




Micro Phaco Splitter
Micro Phako Spalter

01-05-370-1 short tip



Bunnkatu Nucleus Dividing and Separating Hook
Bunnkatu Linsenkern -Teiler und -Trenner

01-05-385-2 TITANIUM
01-05-385-3 BLACK SERIES



Mod.Mainz IOL Manipulator angled
Mod.Mainz IOL Positionshäkchen abgewinkelt

01-05-390-2 TITANIUM




Mod.Frankfurt IOL Manipulator straight Ø 0.2 mm
Mod.Frankfurt IOL Positionshäkchen gerade ø 0.2 mm

01-05-395-1



Sinskey IOL Manipulator straight \varnothing 0.2 mm
Sinskey IOL Positionshäkchen gerade \varnothing 0.2 mm

01-05-400-1
01-05-400-2 TITANIUM
01-05-400-3 BLACK SERIES




Sinskey IOL Manipulator angled \varnothing 0.2 mm
Sinskey IOL Positionshäkchen abgewinkelt \varnothing 0.2 mm

01-05-405-1
01-05-405-2 TITANIUM
01-05-405-3 BLACK SERIES




Sinskey IOL Manipulator straight \varnothing 0.25 mm
Sinskey IOL Positionshäkchen gerade \varnothing 0.25 mm

01-05-410-1
01-05-410-2 TITANIUM
01-05-410-3 BLACK SERIES




Sinskey IOL Manipulator angled \varnothing 0.25 mm
Sinskey IOL Positionshäkchen abgewinkelt \varnothing 0.25 mm

01-05-415-1
01-05-415-2 TITANIUM
01-05-415-3 BLACK SERIES




Lewicky IOL Manipulator angled vaulted Ø 0.15 mm
Lewicky IOL Positionshäkchen abgewinkelt ø 0.15 mm

01-05-420-1
01-05-420-2 TITANIUM
01-05-420-3 BLACK SERIES



Lester IOL Manipulator angled
Lester IOL Positionshäkchen abgewinkelt

01-05-430-1
01-05-430-2 TITANIUM
01-05-430-3 BLACK SERIES




Lester IOL Manipulator angled vaulted
Lester IOL Positionshäkchen abgewinkelt gewölbt

01-05-435-1
01-05-435-2 TITANIUM
01-05-435-3 BLACK SERIES




Kuglen Iris Hook and Lens Manipulator straight
Kuglen Iris und IOL Manipulationshäkchen gerade

01-05-445-1
01-05-445-2 TITANIUM
01-05-445-3 BLACK SERIES



Fenzl IOL Manipulator T-shaped tip straight
Fenzl IOL Positionshäkchen T-geformt gerade

01-05-455-1
01-05-455-2 TITANIUM
01-05-455-3 BLACK SERIES




Fenzl IOL Manipulator T-shaped tip angled
Fenzl IOL Positionshäkchen T-geformt abgewinkelt

01-05-460-1
01-05-460-2 TITANIUM
01-05-460-3 BLACK SERIES




Osher IOL Manipulator Y-Hook straight
Osher IOL Positionshäkchen Y-geformt gerade

01-05-465-2 TITANIUM
01-05-465-3 BLACK SERIES



Jaffe-Maltzmann IOL Manipulator Y-Hook angled
Jaffe-Maltzmann IOL Positionshäkchen Y-geformt
abgewinkelt

01-05-470-2 TITANIUM
01-05-470-3 BLACK SERIES



Rentsch Iris Retractor / IOL Manipulator Boat Hook
Rentsch Irisretraktor und IOL Manipulator
Boots-Häkchen

01-05-475-2 TITANIUM
01-05-475-3 BLACK SERIES



Kaminski Flap lifter
Flap lifter and retreatment spatula
Kaminski Flap lifter
Flap lifter und Spatel

01-05-480-2 TITANIUM
01-05-480-3 BLACK SERIES



BROSS Membrane Spatula, 1.5mm angled - 20G / 0.9mm

01-06-050-1

01-06-050-3 BLACK SERIES



Gimbel Lens Nucleus Spatula angled round Ø 0.25 mm
Gimbel Linsenkern und IOL Spatel abgewinkelt
rund ø 0.25 mm

01-06-100-1

01-06-100-3 BLACK SERIES



Jaffe Lens Spatula angled - flat 0.5 mm wide
Jaffe Positionsspatel flach 0.5 mm breit

01-06-140-1



Obstbaum Lens Spatula angled vaulted 0.25 mm wide
Obstbaum Positionsspatel flach 0.25 mm breit

01-06-150-1

01-06-150-3 BLACK SERIES



Koch Nucleus Spatula, gently curved duckbill shaped tip
Koch Linsenkernspatel, leicht gebogen

01-06-151-1 Stainless Steel

01-06-151-2 TITANIUM

01-06-151-3 BLACK SERIES



Kansas Nucleus Trisector
Kansas Linsenkerntrisektor

01-06-160-1 right

01-06-160-2 TITANIUM

01-06-170-1 left



Kansas Nucleus Bisector angled
Kansas Linsenkernbisektor abgewinkelt

01-06-180-1

01-06-180-2 TITANIUM



Kansas Nucleus Vectis right
Kansas Linsenkervectis rechts

01-06-190-1



Kansas Nucleus Vectis left
Kansas Linsenkervectis links

01-06-200-1



Culler Iris Spatula conical 1 mm
Culler Irisspatel konisch 1 mm

01-06-210-2 TITANIUM



Culler Iris Spatula conical 2 mm
Culler Irisspatel konisch 2 mm

01-06-220-2 TITANIUM



Wecker Iris Spatula silver malleable 2 mm
Wecker Irisspatel silber verformbar 2 mm

01-06-230-1



Barraquer-Fino Cyclodialysis Spatula
Barraquer-Fino Zyklodialysenspatel

01-06-260-1 0.25 x 12 mm

01-06-270-1 0.5 x 12 mm



Mod.Dortmund Phaco-Spatula angled forward
Mod.Dortmund Phakospatel aufgewinkelt

01-06-280-1 0.6 x 12 mm



Castroviejo-Barraquer Cyclodialysis Spatula
Castroviejo-Barraquer Zyklodialysenspatel

01-06-290-1 0.5 x 10 mm
01-06-300-1 1.0 x 10 mm
01-06-310-1 1.0 x 16 mm



Castroviejo Cyclodialysis Spatula 0.5 x 10/16 mm
Castroviejo Zyklodialysenspatel 0.5 x 10/16 mm

01-06-320-1



Castroviejo Synechia Spatula sharp tips
Castroviejo Doppel-Synechiaspatel scharfe Enden

01-06-330-1



PS Platinum Spatula
PS Platin Spatel

01-06-350-1



Kimura Platinum Spatula
Kimura Platin Spatel

01-06-360-1



Fukasaku Lasik Spatula angled fine
Fukasaku Lasikspatel abgewinkelt fein

01-06-370-1

01-06-370-3 BLACK SERIES



LS Lasikspatel abgewinkel medium

01-06-380-1

01-06-380-3 BLACK SERIES



Maddox Flap Elevator for Lasik flat / cylindrical
Maddox Lasikspatel abgewinkel 2endig flach und
zylindrisch

01-06-390-1



Lewis Lens Loop smooth
Lewis Starschlinge glatt

01-06-400-1 small
01-06-400-3 BLACK SERIES
01-06-410-1 large
01-06-410-3 BLACK SERIES



Snellen Lens Loop serrated
Snellen Starschlinge gerieft

01-06-420-1
01-06-420-3 BLACK SERIES



Arlt Lens Loop serrated
Arlt Starschlinge gerieft

01-06-430-1
01-06-430-3 BLACK SERIES




New Orleans Lens Loop smooth
New Orleans Starschlinge glatt

01-06-440-1
01-06-440-3 BLACK SERIES



Meyhoefer Chalazion Curette Flathandle
Meyhoefer Chalazionlöffel

01-06-490-1 Ø 1.0 mm
01-06-500-1 Ø 1.5 mm
01-06-510-1 Ø 1.8 mm
01-06-520-1 Ø 2.0 mm
01-06-530-1 Ø 2.5 mm
01-06-540-1 Ø 3.0 mm
01-06-550-1 Ø 3.5 mm



Hebra Chalazion Curette
01-06-610-1 Ø 2.0 mm
01-06-620-1 Ø 2.5 mm



Wells Enucleation Spoon large
01-06-640-1 small
01-06-650-1 large



Bunge Evisceration Spoon
Bunge Exenterationslöffel

01-06-660-1 small
01-06-670-1 medium
01-06-680-1 large



Barraquer Keratoplastic Spatula perforated
Barraquer Keratoplastiksiebspatel

01-06-690-1



Muldoon Lacrimal Dilator graduated tip
Muldoon Tränenkanaldilatator

01-07-100-1



Wilder Lacrimal Dilator #1 fine taper
Wilder Tränenkanaldilatator #1 fein

01-07-105-1

01-07-105-3 BLACK SERIES



Wilder Lacrimal Dilator #2 medium taper
Wilder Tränenkanaldilatator #2 mittel

01-07-110-1

01-07-110-3 BLACK SERIES



Wilder Lacrimal Dilator #3 heavy taper
Wilder Tränenkanaldilatator #3 stark

01-07-115-1

01-07-115-3 BLACK SERIES



Sisler Lacrimal Punctum Dilator
Sisler Dilatator für Tränenpünckchen

01-07-120-1



Ruedemann Lacrimal Dilator Infant very fine
Ruedemann Tränenkanaldilatator Säuglinge

01-07-125-1

01-07-125-3 BLACK SERIES



Castroviejo Lacrimal Dilator fine / medium taper
Castroviejo Tränenkanaldilatator 2endig

01-07-130-1



01-07-130-3 BLACK SERIES



LG Lacrimal Punctal Gauge 0.5 and 0.6 mm
LG Tränenpunktlehre 2endig 0.5 und 0.6 mm

01-07-135-2 TITANIUM







Bowman Lacrimal Probe cylindrical sterling silver

Size

01-07-140-1	#0000/000	Ø 0.4 mm / Ø 0.6 mm
01-07-145-1	#00/0	Ø 0.7 mm / Ø 0.8 mm
01-07-150-1	#0/1	Ø 0.8 mm / Ø 0.9 mm
01-07-155-1	#1/2	Ø 0.9 mm / Ø 1.1 mm
01-07-160-1	#2/3	Ø 1.1 mm / Ø 1.3 mm
01-07-165-1	#3/4	Ø 1.3 mm / Ø 1.4 mm
01-07-170-1	#4/5	Ø 1.4 mm / Ø 1.5 mm
01-07-175-1	#5/6	Ø 1.5 mm / Ø 1.6 mm
01-07-180-1	#6/7	Ø 1.6 mm / Ø 1.8 mm
01-07-185-1	#7/8	Ø 1.8 mm / Ø 1.9 mm



Bowman Lacrimal Probe bulbous sterling silver

Size

01-07-190-1	#0000/000	Ø 0.4 mm / Ø 0.6 mm
01-07-195-1	#00/0	Ø 0.7 mm / Ø 0.8 mm
01-07-200-1	#0/1	Ø 0.8 mm / Ø 0.9 mm
01-07-205-1	#1/2	Ø 0.9 mm / Ø 1.1 mm
01-07-210-1	#2/3	Ø 1.1 mm / Ø 1.3 mm
01-07-215-1	#3/4	Ø 1.3 mm / Ø 1.4 mm
01-07-220-1	#4/5	Ø 1.4 mm / Ø 1.5 mm
01-07-225-1	#5/6	Ø 1.5 mm / Ø 1.6 mm
01-07-230-1	#6/7	Ø 1.6 mm / Ø 1.8 mm
01-07-235-1	#7/8	Ø 1.8 mm / Ø 1.9 mm

Worst Pigtail Probe with suture holes
Worst Pigtail Sonde mit Nahtöhr

- 01-07-240-1** large
- 01-07-240-3** BLACK SERIES
- 01-07-241-1** small
- 01-07-241-3** BLACK SERIES



Harms Trabeculotomy Probe
Harms Trabekularsonde


- 01-07-245-1** right
- 01-07-245-2** right TITANIUM
- 01-07-250-1** left
- 01-07-250-2** left TITANIUM



Wadsworth-Todd Cautery with copper ball
Wadsworth-Todd Kauter mit Kupferkugel

- 01-07-265-1**





Cautery Sets (Disposable - 10/Box)

High Temperature Cauteries:

- 01-07-275-1 High-temp fine tip, 4-yr. shelf life
- 01-07-276-1 High-temp loop tip, 4-yr. shelf life
- 01-07-277-1 High-temp adjustable fine tip, 4-yr. shelf life

Low Temperature Cauteries:

- 01-07-280-1 Low-temp fine tip, 4-yr. shelf life
- 01-07-281-1 Low-temp elongated fine tip, 4-yr. shelf life
- 01-07-282-1 Low-temp adjustable fine tip, 4-yr. shelf life



Algerbrush Rust Ring Remover

For Corneal Rust Ring:

- 01-07-294-1 Algerbrush II .5 mm Complete Instrument
- 01-07-295-1 Algerbrush II 1 mm Complete Instrument

For Pterygium:

- 01-07-296-1 Complete Pterygium Instrument with 2.5 mm Round Fine Grit Diamond
- 01-07-297-1 Complete Pterygium Instrument with 3.5 mm Round Fine Grit Diamond
- 01-07-298-1 Complete Pterygium Instrument with 4.0 mm Round Fine Grit Diamond
- 01-07-299-1 Complete Pterygium Instrument with 5.0 mm Round Fine Grit Diamond

Spare parts on request !!!

Elliot Trephine Handle
with Scleral Trephine Ø 1.0 / 1.5 / 2 mm
01-07-315-1



Elliot Handle
01-07-320-1 with chuck for trephine



Elliot Scleral Trephine
01-07-325-1 Ø 1.0 mm



Elliot Scleral Trephine
01-07-330-1 Ø 1.5 mm
01-07-335-1 Ø 2.0 mm




**Castroviejo** Keratoplasty Trephine

01-07-340-1	Ø 5.0 mm
01-07-345-1	Ø 5.5 mm
01-07-350-1	Ø 6.0 mm
01-07-355-1	Ø 6.5 mm
01-07-360-1	Ø 7.0 mm
01-07-365-1	Ø 7.5 mm
01-07-370-1	Ø 8.0 mm
01-07-375-1	Ø 9.0 mm

**Franceschetti** Keratoplasty Trephine

01-07-380-1	Ø 7.0 mm
01-07-385-1	Ø 7.1 mm
01-07-390-1	Ø 8.0 mm
01-07-395-1	Ø 8.1 mm




Martin Keratoplasty Trephine

01-07-400-1	Ø 5.0 mm
01-07-405-1	Ø 5.5 mm
01-07-410-1	Ø 6.0 mm
01-07-415-1	Ø 6.5 mm
01-07-420-1	Ø 6.7 mm
01-07-425-1	Ø 7.0 mm
01-07-430-1	Ø 7.2 mm
01-07-435-1	Ø 7.5 mm
01-07-440-1	Ø 7.7 mm
01-07-445-1	Ø 8.0 mm
01-07-450-1	Ø 8.2 mm
01-07-455-1	Ø 8.5 mm
01-07-460-1	Ø 8.7 mm
01-07-465-1	Ø 9.0 mm
01-07-470-1	Ø 9.5 mm



Flieringa Scleral Fixation Ring

01-07-475-1	Ø 14 mm
01-07-480-1	Ø 15 mm
01-07-485-1	Ø 16 mm
01-07-490-1	Ø 17 mm
01-07-495-1	Ø 18 mm
01-07-500-1	Ø 19 mm
01-07-505-1	Ø 20 mm
01-07-510-1	Ø 21 mm
01-07-515-1	Ø 22 mm
01-07-520-1	Ø 23 mm
01-07-525-1	Ø 24 mm



Hoffer Optical Center Marker with cross hairs

01-08-100-1	3.00 mm
01-08-110-1	3.25 mm
01-08-120-1	3.50 mm
01-08-130-1	3.75 mm
01-08-140-1	4.00 mm
01-08-150-1	4.25 mm
01-08-160-1	4.50 mm
01-08-170-1	4.75 mm
01-08-180-1	5.00 mm
01-08-190-1	5.50 mm
01-08-200-1	6.00 mm
01-08-210-1	7.00 mm
01-08-220-1	9.00 mm
01-08-230-1	11.00 mm

LASEK Solution Reservoir
LASEK Lösungsmittelring

01-08-231-1	8.0 mm
01-08-232-1	8.5 mm
01-08-234-1	9.0 mm



Anis Radial Marker ID 4 mm, OD 12.5 mm
Anis Radialmarkeur ID 4 mm AD 12.5 mm

01-08-240-1	8 Blades
01-08-240-3	BLACK SERIES



Anis Radial Marker ID 4 mm, OD 12.5 mm
Anis Radialmarkeur ID 4 mm AD 12.5 mm

01-08-250-1 12 Blades
01-08-250-3 BLACK SERIES



Anis Suture Placement Marker ID 8 mm
Anis Nahtplatzierungsmarkeur ID 8 mm

01-08-260-1 8 Blades
01-08-260-3 BLACK SERIES



OSHER NEUMANN Corneal Marker
OSHER NEUMANN Hornhautmarker

01-08-261-1 12 Blades
01-08-261-3 BLACK SERIES
01-08-262-1 16 Blades
01-08-262-3 BLACK SERIES
01-08-263-1 8 Blades
01-08-263-3 BLACK SERIES



Bonn Globe Fixation Ring
Bonn Fixationsring PRK und Lasik

01-08-270-1 PRK and Lasik



Bonn Globe Fixation Ring 3/4 Clear Cornea and IOL
Bonn Fixationsring 3/4 Clear Cornea und IOL

01-08-280-1

01-08-280-3 BLACK SERIES



Bores Axis Marker correction of astigmatism
Bores Axialmarker für Astigmatismuskorrektur

01-08-290-1

01-08-290-3 BLACK SERIES



Thornton Swivel Fixation Ring Ø 13 mm 3/4
Thornton Fixationsring FG-Drehgelenk ø13 mm 3/4

01-08-300-1

01-08-300-2 TITANIUM



Thornton Swivel Fixation Ring Ø 16 mm 3/4
Thornton Fixationsring FG-Drehgelenk ø16 mm 3/4

01-08-310-1

01-08-310-2 TITANIUM



Gimbel Swivel Fixation Ring Ø 13 mm 3/4
Gimbel Fixationsring FG-Drehgelenk ø13 mm 3/4

01-08-311-2 TITANIUM



Grandon Transverse Incision Marker 3 mm
Grandon T-Inzision's Markeur 3 mm

01-08-320-2 TITANIUM



Thornton Limbal Incision Ruler 1-10 mm
Thornton Limbus Inzisionsmarkeur 1-10 mm

01-08-330-1
01-08-330-3 BLACK SERIES



Thornton Double Corneal Ruler
6-8 mm transverse incisions 0.5-5.0 mm
Thornton Doppel-Hornhautmarkeur
OZ 6-8 mm T-Inzision 0.5-5.0 mm

01-08-340-1




Lindstrom Accurate Incision Marker
Lindstrom Hornhautmarkeur

01-08-350-1 7.0 mm optical zone




Mendez Corneal Marker for Lasik
Mendez Hornhautmarkeur Lasik

01-08-360-1
01-08-360-3 BLACK SERIES




Chayet Corneal Marker for Lasik Ø 8.5 mm
Chayet Hornhautmarkeur Lasik Ø 8.5 mm

01-08-370-1
01-08-370-3 BLACK SERIES



Lu Corneal Marker for Lasik
Lu Hornhautmarkeur Lasik

01-08-380-1
01-08-380-3 BLACK SERIES




Ultra-Grip I/A Cannula-Handpiece
Ultra-Grip Spül- und Saughandgriff

01-09-100-1 Set with 5 Canullas
01-09-104-1 Handle only




I/A Cannula straight Aspiration Port Size 0.3 mm
I/A Kanüle gerade Saugöffnung 0.3 mm

01-09-110-1




I/A Cannula curved Aspiration Port Size 0.3 mm
I/A Kanüle gebogen Saugöffnung 0.3 mm

01-09-115-1




I/A Cannula 45° angled Aspiration Port Size 0.3 mm
I/A Kanüle 45° abgewinkelt Saugöffnung 0.3 mm

01-09-120-1




I/A Cannula 90° angled Aspiration Port Size 0.3 mm
I/A Kanüle 90° abgewinkelt Saugöffnung 0.3 mm

01-09-125-1



I/A Cannula Binkhorst Tip Aspiration Port front opening
I/A Kanüle Binkhorst Saugöffnung frontal

01-09-130-1




AI Air Injection Cannula
AI Luftkanüle abgewinkelt

01-09-135-1 30G - 5 mm angled

01-09-140-1 27G - 5 mm angled

01-09-145-1 30G - 7 mm angled


01-09-150-1 27G - 7 mm angled



McIntyre Anterior Chamber Cannula
McIntyre Vorderkammerkanüle abgewinkelt

01-09-155-1 30G angled


01-09-160-1 27G angled



Hydrodissection

Tübingen Air Injection Cannula curved
Tübingen Luftkanüle gebogen


01-09-165-1 30G



Hydrodissection

Sauter Air Injection Cannula curved
Sauter Luftkanüle spatelförmiges Ende

01-09-170-1 27G




Hydrodissection

Rainin Air Injection Cannula angled
Rainin Luftkanüle abgewinkelt

01-09-175-1 30G

01-09-180-1 27G



Castroviejo Cyclodialysis Cannula angled
Castroviejo Zyklodialysenkanüle abgewinkelt


01-09-185-1 10 mm

01-09-190-1 15 mm




Kelman Cyclodialysis Cannula
Kelman Zyklotodialysenkanüle

01-09-195-1 19G right
01-09-200-1 19G left




Troutman Alpha-Chymotrypsin Cannula
Troutman Alpha-Chymotrypsin-Kanüle

01-09-205-1 25G



Welsh Cortex Extractor angled flat olive
Welsh Kortex-Extraktor abgewinkelt flacher Olive

01-09-210-1 25G





Welsh Cortex Extractor
Welsh Kortex-Extraktor

01-09-215-1 25G right
01-09-220-1 25G left



Simcoe Cortex Extractor curved
Simcoe Kortex-Extraktor gebogen


01-09-225-1 23G



Atkinson Retrobulbar Cannula
Atkinson Retrobulbäre Kanüle

01-09-230-1 23G - 40 mm


01-09-235-1 23G - 60 mm



Gills/Welsh Cortex Extractor
Gills/Welsh Kortex-Extraktor abgewinkelt


01-09-240-1 25G - 45° angled - 5 mm

01-09-245-1 25G - 45° angled - 7 mm




Anel Lacrimal Cannula straight
Anel Tränenwegkanüle gerade

01-09-250-1 23G




Anel Lacrimal Cannula 90° curved
Anel Tränenwegkanüle 90° gebogen

01-09-255-1 23G




LC Lacrimal Cannula straight
LC Tränenwegkanüle gerade

01-09-260-1 23G




LC Lacrimal Cannula curved
LC Tränenwegkanüle gebogen

01-09-265-1 23G




Shahinian Lacrimal Cannula straight
Shahinian Tränenwegkanüle gerade

01-09-270-1 23G




Pioneer Lacrimal Cannula straight
Pioneer Tränenwegkanüle gerade

01-09-275-1 19 / 27G



McIntyre Lacrimal Cannula dual sideports 0.3 mm
McIntyre Tränenwegkanüle seitliche Öffnungen 0.3mm

01-09-280-1 23G



Steinkogler-Bangerter Lacrimal Cannula, lightly curved
with sideport opening
Steinkogler-Bangerter Tränenwegkanüle, leicht gebogen
mit seitlicher Öffnung


01-09-283-1 21G



Mod.Bonn Irrigating Cannula angled
Mod.Bonn Spülkanüle abgewinkelt spatelförmig



01-09-285-1 19G spatulated tip

01-09-290-1 20G spatulated tip




Simcoe Irrigating-Aspirating Cannula
Simcoe Spül-Saugkanüle

01-09-295-1 25/23G
01-09-300-1 25/23G reversed




Simcoe I/A Cannula for 12 o'clock aspiration
Simcoe I/A Kanüle für 12 Uhr Aspiration

01-09-305-1 23G right
01-09-310-1 23G left





Friedburg Irrigating-Aspirating Cannula
Friedburg Spül-Saugkanüle

01-09-315-1



Gills Irrigating-Aspirating Cannula angled
Gills Spül-Saugkanüle abgewinkelt

01-09-320-1 23G



Kelman Cystotome Cannula
Kelman Zystotomkanüle scharfe Schneide -Spitze

01-09-325-1 22G sharp edge & point
01-09-330-1 22G blunt edge & sharp point



Kratz Cystotome Cannula 22G sharp cutting edge and point
Kratz Zystotomkanüle scharfe Schneide -Spitze

01-09-335-1




Sato Cystotome Knife
Sato Zystotommesser Doppelschneide

01-09-340-1 23G




Kratz-Ziegler Knife Cannula 1.5 mm
Kratz-Ziegler Kanülemesser 1.5 mm

01-09-345-1 25G




HVL Injection Cannula for high viscous liquids
HVL Injektionskanüle für hochvisköse Flüssigkeiten

01-09-350-1 21G




Kratz Capsule Scraper
Kratz Kapselifeile abgewinkelt

01-09-355-1 23G 1 mm angled diamond dusted
01-09-360-1 23G 2 mm angled diamond dusted



Simcoe Posterior Capsule Polisher
Simcoe Hinterkapselpolierer abgewinkelt

01-09-365-1 23G angled diamond dusted




Jensen Capsule Polisher
Jensen Kapselpolierer

01-09-370-1 27G diamond dusted
01-09-375-1 23G diamond dusted




Knolle Posterior Capsule Polisher angled
Knolle Hinterkapselpolierer 2 mm abgewinkelt

01-09-380-1 30G - 2 mm diamond dusted
01-09-385-1 23G - 2 mm diamond dusted




Simcoe Posterior Capsule Polisher curved
Simcoe Hinterkapselpolierer gebogen aufgeraut

01-09-390-1 23G sandblasted




Graether Cannula
Graether Kragenknopfkanüle gerade

01-09-395-1 23G Micro Iris Retractor & posterior capsule polisher straight



Graether Cannula
Graether Kragenknopfkanüle abgewinkelt

01-09-400-1 23G Micro Iris Retractor & posterior capsule polisher angled




McIntyre Irrigating Cannula with Iris Hook angled
McIntyre Spülkanüle abgewinkelt mit Irishäkchen

01-09-405-1 23G




Knolle-Pearce Irrigating Vectis Loop
Knolle-Pearce Linsenschlinge mit Spülung

01-09-410-1 23G three front opening ports





Sheets Irrigating Vectis Loop
Sheets Linsenschlinge gerieft mit Spülung

01-09-415-1 20G serrated three front opening ports




Simcoe Irrigating Lens Nucleus Loop
Simcoe Linsenkernschlinge gerieft mit Spülung, rechts

01-09-420-1 25G fine serrated two front opening ports right




Simcoe Irrigating Lens Nucleus Loop
Simcoe Linsenkernschlinge gerieft mit Spülung, links

01-09-425-1 25G fine serrated two front opening ports left



Bell Erisiphake with silicone bulb
Bell Kapselsauger

01-09-430-1



SB Silicone Bulb for Bell Erisiphake
SB Silikonball für Bell Kapselsauger

01-09-435-1




Bishop-Harmon Anterior Chamber Irrigator
Bishop-Harmon Vorderkammerspülkanülensatz

01-09-440-1 Complete



SB Silicone Bulb for Anterior Chamber Irrigator
SB Silikonball für Vorderkammerspülkanülensatz

01-09-445-1




LC Adaptor Luer-lock
LC Verbindung Luer-lock für Silikonball

01-09-450-1 for Anterior Chamber Irrigator



Bishop-Harmon Irrigating Cannula angled
Bishop-Harmon Spülkanüle 19G abgewinkelt

01-09-455-1 19G



BAH Aspiration Handpiece male/female fitting
BAH Bimanueller Saugkanülenhandgriff male/male

01-09-600-2 TITANIUM - Yellow




BIH Irrigating Handpiece male/female fitting
BIH Bimanueller Spülkanülenhandgriff male/female

01-09-605-2 TITANIUM - Blue




Brauweiler Irrigating Cannula angled
Brauweiler Spülkanüle 21G abgewinkelt

01-09-610-1 21G



Koch Aspirating Cannula angled sandblasted
Koch Saugkanüle 21G abgewinkelt aufgerauht

01-09-615-1 21G







Heaflinger Aspiration Handpiece
Heaflinger Bimanueller Saugkanülenhandgriff




Heaflinger Irrigation Handpiece
Heaflinger Bimanueller Spülkanülenhandgriff




01-09-620-2 21G curved, 0.35 mm



01-09-635-2 21G curved dual 0.5 mm




01-09-625-2 21G curved, 0.30 mm with tungsten carbide tip




01-09-640-2 22G, overall length 105 mm



01-09-630-2 21G angled with tungsten carbide tip




01-09-645-2 21G, overall length 105 mm angled




Nagahara Phacochooper 0.9 mm working part
01-09-650-2 right hand
01-09-655-2 left hand



Rosen Phacochooper
01-09-660-2 wedge-shaped 0.9 mm working part
01-09-665-2 wide-shaped 1.7 mm working part



01-09-670-2 Hook Phacochooper 0.9 mm working part



01-09-675-2 Micro-Finger Phacochooper finger-shaped, 1.4 mm working part

Will be replaced by a new Version

Barraquer Cilia Forceps 11.5 cm
Barraquer Zilienpinzette 11.5 cm

01-10-100-1



Beer Cilia Forceps 9 cm
Beer Zilienpinzette 9 cm

01-10-106-1




Storz Cilia Forceps 8.5 cm
Storz Zilienpinzette 8.5 cm

01-10-108-1




Micro Adson Dressing Forceps 12 cm straight
Micro Adson Pinzette anatomisch 12 cm gerade

01-10-110-1 serrated - fine
01-10-110-2 TITANIUM



Micro Adson Tissue Forceps 12 cm straight
Micro Adson Pinzette chirurgisch 12 cm gerade


01-10-112-1 1x2 teeth - fine
01-10-112-2 **TITANIUM**



Graefe Dressing Forceps 10 cm, serrated
Graefe Pinzette anatomisch 10 cm gerade gerieft

01-10-114-1 1.0 mm straight
01-10-115-1 0.7 mm straight
01-10-116-1 0.5 mm straight


01-10-117-1 1.0 mm curved
01-10-118-1 0.7 mm curved
01-10-119-1 0.5 mm curved



Graefe Tissue Forceps 10 cm , 1x2 teeth
Graefe Pinzette chirurgisch 10 cm 1x2 Zähne


01-10-120-1 1.0 mm straight
01-10-121-1 0.7 mm straight
01-10-122-1 0.5 mm straight

01-10-123-1 1.0 mm curved
01-10-124-1 0.7 mm curved
01-10-125-1 0.5 mm curved




Bishop-Harmon Dressing Forceps 8.5 cm straight
Bishop-Harmon Pinzette anatomisch 8.5 cm gerade

01-10-130-1 serrated - standard
01-10-131-1 serrated - delicate
01-10-131-2 **TITANIUM**




Bishop-Harmon Tissue Forceps 8.5 cm straight
Bishop-Harmon Pinzette chirurgisch 8.5 cm gerade

- 01-10-132-1** 1x2 teeth - standard
- 01-10-133-1** 1x2 teeth - delicate
- 01-10-133-2** TITANIUM




Bonn Iris Forceps 7 cm serrated
Bonn Irispinzette 7cm gerieft

- 01-10-135-1** straight
- 01-10-135-3** BLACK SERIES
- 01-10-136-1** angled
- 01-10-136-3** BLACK SERIES




Bonn Iris Forceps 7 cm straight
Bonn Irispinzette 7cm gerade

- 01-10-143-1** B-teeth 0.12 mm
- 01-10-143-3** BLACK SERIES




Bonn Iris Forceps 7 cm angled 6 mm
Bonn Irispinzette 7 cm abgewinkelt 6 mm

- 01-10-146-1** B-teeth 0.12 mm
- 01-10-146-3** BLACK SERIES




Bonn Iris Forceps 10 cm straight
Bonn Irispinzette 10 cm gerade

01-10-150-1 B-teeth 0.12 mm
01-10-150-3 BLACK SERIES




Troutman-Barraquer Colibri Forceps 7.5 cm angled
Troutman-Barraquer Kolibripinzette 7.5 cm abgewinkelt

01-10-155-1 C-teeth 0.4 mm
01-10-155-3 BLACK SERIES




Gill Iris Forceps 7.5 cm curved
Gill Irispinzette 7.5 cm gebogen

01-10-165-1 1x2 B-teeth 0.3 mm
01-10-165-3 BLACK SERIES



Hess Iris Forceps 10 cm S-angled
Hess Irispinzette 10 cm S-abgewinkelt

01-10-175-1 B-teeth 0.3 mm
01-10-175-3 BLACK SERIES



Paufique Iris Forceps 8.5 cm straight with Tying Platform
Paufique Irispinzette 8.5 cm gerade mit Fadenplatte

- 01-10-180-1 B-teeth 0.12 mm
- 01-10-180-3 BLACK SERIES
- 01-10-182-1 B-teeth 0.30 mm
- 01-10-182-3 BLACK SERIES



Paufique Suturing Forceps 8.5 cm straight with Tying Platform
Paufique Nahtpinzette 8.5 cm gerade mit Fadenplatte

- 01-10-184-1 C-teeth 0.12 mm
- 01-10-184-3 BLACK SERIES
- 01-10-186-1 C-teeth 0.30 mm
- 01-10-186-3 BLACK SERIES




Bonn Iris Forceps 10 cm straight with Tying Platform
Bonn Irispinzette 10 cm gerade mit Fadenplatte

- 01-10-210-1 B-teeth 0.12 mm Stainless Steel
- 01-10-210-2 TITANIUM
- 01-10-210-3 BLACK SERIES




Bonn Iris Forceps 10 cm straight with Tying Platform
Bonn Irispinzette 10 cm gerade mit Fadenplatte

- 01-10-211-1 B-teeth 0.12 mm Stainless Steel
- 01-10-211-2 TITANIUM
- 01-10-211-3 BLACK SERIES




Bonn Iris Forceps 10 cm angled with
Tying Platform B-teeth 0.12 mm
Bonn Pinzette chirurgisch 10 cm abgewinkelt mit
Fadenplatte B-Zähne 0.12 mm

- 01-10-215-1 Stainless Steel
- 01-10-215-2 TITANIUM
- 01-10-215-3 BLACK SERIES




Kelman-McPherson Forceps 10 cm angled,
with Tying Platform, B-teeth 0.3 mm
Kelman-McPherson Pinzette 10 cm abgewinkelt
mit Fadenplatte B-Zähne 0.3 mm

- 01-10-220-1 Stainless Steel
- 01-10-220-2 TITANIUM
- 01-10-220-3 BLACK SERIES




Bonn Iris Forceps 7 cm straight with Tying Platform
Bonn Irispinzette 7 cm gerade mit Fadenplatte

- 01-10-225-1 B-teeth 0.12 mm
- 01-10-225-3 BLACK SERIES




Bonn Iris Forceps 7 cm angled 6 mm with Tying Platform
Bonn Irispinzette 7 cm abgewinkelt mit 6 mm Fadenplatte

- 01-10-226-1 B-teeth 0.12 mm
- 01-10-226-3 BLACK SERIES




Harms Colibri Forceps 7.5 cm with Tying Platform
Harms Kolibripinzette 7.5 cm mit Fadenplatte

- 01-10-230-1** B-teeth 0.12 mm
- 01-10-230-3** BLACK SERIES
- 01-10-234-1** B-teeth 0.40 mm
- 01-10-234-3** BLACK SERIES




Harms Colibri Forceps 7.5 cm with Tying Platform
Harms Kolibripinzette 7.5 cm mit Fadenplatte

- 01-10-230-2** B-teeth 0.12 mm **TITANIUM**
- 01-10-234-2** B-teeth 0.40 mm **TITANIUM**




Troutman-Barraquer Colibri Forceps 7.5 cm
Troutman-Barraquer Kolibripinzette 7.5 cm

- 01-10-235-1** B-teeth 0.12 mm
- 01-10-235-3** BLACK SERIES




Troutman-Barraquer Colibri Forceps 7.5 cm
Troutman-Barraquer Kolibripinzette 7.5 cm

- 01-10-235-2** B-teeth 0.12 mm **TITANIUM**




Troutman Colibri Forceps 7 cm Bayonet-shaped
Troutman Kolibripinzette 7 cm Bajonett-Form

- 01-10-240-1 B-teeth 0.12 mm
- 01-10-240-2 TITANIUM
- 01-10-240-3 BLACK SERIES




Barraquer Colibri Forceps 7.5 cm, Tying Platform
Barraquer Kolibripinzette 7.5 cm mit Fadenplatte

- 01-10-245-1 C-teeth 0.12 mm
- 01-10-245-2 TITANIUM
- 01-10-245-3 BLACK SERIES
- 01-10-247-1 C-teeth 0.25 mm
- 01-10-247-2 TITANIUM
- 01-10-247-3 BLACK SERIES




CF Colibri Forceps 7.5 cm with Tying Platform
CF Kolibripinzette 7.5 cm mit Fadenplatte

- 01-10-250-1 C-teeth 0.4 mm
- 01-10-250-2 TITANIUM
- 01-10-250-3 BLACK SERIES




Mod.Zürch Suturing Forceps 7.5 cm Tying Platform
Mod.Zürch Nahtpinzette 7.5 cm gerade mit Fadenplatte

- 01-10-255-1 C-teeth 0.12 mm
- 01-10-255-2 TITANIUM
- 01-10-255-3 BLACK SERIES




Mod.Zürich Suturing Forceps 7.5 cm angled Tying Platform
Mod.Zürich Nahtpinzette 7.5 cm abgewinkelt mit Fadenplatte

- 01-10-260-1 C-teeth 0.12 mm
- 01-10-260-2 TITANIUM
- 01-10-260-3 BLACK SERIES




Kirby-Artur Corneal Forceps 9 cm angled
Kirby-Artur Hornhautpinzette 9 cm abgewinkelt

- 01-10-275-1 C-teeth 0.2 mm
- 01-10-275-3 BLACK SERIES




Castroviejo Suturing Forceps 10 cm straight Tying Platform
Castroviejo Nahtpinzette 10 cm gerade mit Fadenplatte

- 01-10-280-1 C-teeth 0.12 mm
- 01-10-280-2 TITANIUM
- 01-10-280-3 BLACK SERIES
- 01-10-285-1 C-teeth 0.30 mm
- 01-10-285-2 TITANIUM
- 01-10-285-3 BLACK SERIES
- 01-10-290-1 C-teeth 0.50 mm
- 01-10-290-2 TITANIUM
- 01-10-290-3 BLACK SERIES
- 01-10-295-1 C-teeth 1.00 mm
- 01-10-295-2 TITANIUM
- 01-10-295-3 BLACK SERIES




Castroviejo Suturing Forceps 10 cm straight Tying Platform
Castroviejo Nahtpinzette 10 cm gerade mit Fadenplatte

- 01-10-282-1 C-teeth 0.12 mm
- 01-10-282-3 BLACK SERIES
- 01-10-287-1 C-teeth 0.30 mm
- 01-10-287-3 BLACK SERIES
- 01-10-292-1 C-teeth 0.50 mm
- 01-10-292-3 BLACK SERIES
- 01-10-297-1 C-teeth 1.00 mm
- 01-10-297-3 BLACK SERIES




Paufique Universal Suturing Forceps 8.5 cm
Paufique Allzweckpinzette 8.5 cm mit Fadenplatte

01-10-300-1 Tying Platform C-teeth 0.5 mm
01-10-300-3 BLACK SERIES





Pierse Corneal Forceps 8 cm, straight
Pierse Hornhautpinzette 8 cm gerade

01-10-310-1 1x1 Pierse teeth Ø 0.2 mm
01-10-310-2 1x1 Pierse teeth Ø 0.2 mm **TITANIUM**



Pierse Corneal Forceps 8 cm, curved
Pierse Hornhautpinzette 8 cm gebogen

01-10-315-1 1x1 Pierse teeth Ø 0.2 mm
01-10-315-2 1x1 Pierse teeth Ø 0.2 mm **TITANIUM**




Pierse Colibri Corneal Forceps 8 cm
Pierse Colibri Hornhautpinzette 8 cm

01-10-325-1 1x1 Pierse teeth Ø 0.2 mm
01-10-325-2 1x1 Pierse teeth Ø 0.2 mm **TIT**
01-10-330-1 1x1 Pierse teeth Ø 0.2 mm
01-10-330-2 1x1 Pierse teeth Ø 0.2 mm **TIT**



Buratto Lasik Flap Forceps 10 cm curved disc shaped
Buratto Lasik Hornhaut Pinzette 10 cm gebogen - gerieft

01-10-335-1 serrated jaws




Arita Meibomian Gland Compressor broad

01-10-337-1




Harms Tying Forceps 10 cm straight Tying Platform 7 mm
Harms Fadenpinzette 10 cm gerade mit Fadenplatte 7 mm

01-10-340-1 USP 7-0--9-0
01-10-340-3 BLACK SERIES




Harms Tying Forceps 10 cm curved Tying Platform 7 mm
Harms Fadenpinzette 10 cm gebogen mit Fadenplatte 7 mm

01-10-342-1 USP 7-0 -- 9-0
01-10-342-3 BLACK SERIES




Jaffe Tying Forceps 10 cm straight Tying Platform 5 mm
Jaffe Fadenpinzette 10 cm gerade mit Fadenplatte 5 mm

01-10-345-1 USP 8-0 -- 11-0 Stainless Steel
01-10-345-2 TITANIUM
01-10-345-3 BLACK SERIES




Jaffe Tying Forceps 10 cm curved Tying Platform 5 mm
Jaffe Fadenpinzette 10 cm gebogen mit Fadenplatte 5 mm

- 01-10-350-1 USP 8-0 -- 11-0 Stainless Steel
- 01-10-350-2 **TITANIUM**
- 01-10-350-3 BLACK SERIES




McPherson Tying Forceps 10 cm angled
Tying Platform 5 mm
McPherson Fadenpinzette 10 cm abgewinkelt
mit Fadenplatte 5 mm

- 01-10-355-1 USP 8-0 -- 11-0 Stainless Steel
- 01-10-355-2 **TITANIUM**
- 01-10-355-3 BLACK SERIES



K-McPherson Tying Forceps 10 cm angled
Tying Platform 8 mm
K-McPherson Fadenpinzette 10 cm abgewinkelt
mit Fadenplatte 8 mm

- 01-10-360-1 USP 8-0 -- 11-0 Stainless Steel
- 01-10-360-2 **TITANIUM**
- 01-10-360-3 BLACK SERIES




K-McPherson Tying Forceps 10 cm angled
Tying Platform 8 mm
K-McPherson Fadenpinzette 10 cm abgewinkelt
mit Fadenplatte 8 mm

- 01-10-361-1 Roundhandle Stainless Steel
- 01-10-361-2 **TITANIUM**
- 01-10-361-3 BLACKSERIES




Tennant Tying Forceps 10 cm straight
Tying Platform 6 mm
Tennant Fadenpinzette 10 cm gerade
mit Fadenplatte 6 mm

- 01-10-365-1 USP 8-0 -- 11-0 Stainless Steel
- 01-10-365-2 **TITANIUM**
- 01-10-365-3 BLACK SERIES




Tennant Tying Forceps 10 cm curved
Tying Platform 6 mm
Tennant Fadenpinzette 10 cm gebogen
mit Fadenplatte 6 mm

- 01-10-370-1 USP 8-0 -- 11-0 Stainless Steel
- 01-10-370-2 **TITANIUM**
- 01-10-370-3 BLACK SERIES



Paufique Tying Forceps 8.5 cm straight
Tying Platform 5.5 mm
Paufique Fadenpinzette 8.5 cm gerade
mit Fadenplatte 5.5 mm

- 01-10-375-1 USP 5-0 -- 8-0
- 01-10-375-3 BLACK SERIES
- 01-10-380-1 USP 7-0 -- 9-0
- 01-10-380-3 BLACK SERIES



Helveston Tying Forceps standard
Helveston Fadenpinzette standard

- 01-10-382-1
- 01-10-382-3 BLACK SERIES




Graefe Fixation Forceps 11 cm
Graefe Fixierpinzette 11 cm

01-10-385-1 without lock
01-10-387-1 with lock




Elschnig Superior Rectus Forceps 10 cm straight
Elschnig Zügelnahtpinzette 10 cm gerade

01-10-390-1 1x2 C-teeth delicate




Elschnig Superior Rectus Forceps 10 cm curved
Elschnig Zügelnahtpinzette 10 cm gebogen

01-10-392-1 1x2 C-teeth delicate




Lester Fixation Forceps 9 cm straight
Lester Fixierpinzette 9 cm gerade

01-10-395-1 1x2 B-teeth
01-10-397-1 2x3 B-teeth




Green Fixation Forceps 10 cm T-Shape
Green Fixierpinzette 10 cm T-Form

01-10-405-1 5 mm jaws fine teeth



Green Fixation Forceps 10 cm T-Shape
Green Fixierpinzette 10 cm T-Form 10mm

01-10-407-1 10 mm jaws fine teeth



Moody Fixation Forceps with lock
Moody Fixierpinzette mit Sperre


01-10-410-1 left

01-10-411-1 right



Kansas Nucleus Forceps 10 cm for extracting lens nucleus fragments delicate
Kansas Linsenkernpinzette 10 cm zur Entfernung der Fragmente abgewinkelt fein

01-10-440-1



Akahoshi Phaco Prechopper
4.5 mm long delicate jaws with conventional smooth sharp tips for cracking, dividing and rotating the nucleus without sculpting, 120 mm long

- 01-10-445-1 straight
- 01-10-445-2 straight **TITANIUM**
- 01-10-450-1 angled
- 01-10-450-2 angled **TITANIUM**



BROSS Inamura Capsulorhexis Forceps 12 cm
BROSS Inamura Kapsulorhexispinzette 12 cm

- 01-10-453-2 **TITANIUM**




Inamura Capsulorhexis Forceps 10 cm
Inamura Kapsulorhexispinzette 10 cm

- 01-10-454-2 **TITANIUM**





BROSS Capsulorhexis Forceps 10 cm
thin shanks 11 mm sharp tips
BROSS Kapsulorhexispinzette 10 cm
feine Schenkel 11 mm spitze Plattform

- 01-10-455-1 Stainless Steel
- 01-10-455-2 **TITANIUM**
- 01-10-455-3 BLACK SERIES




Masket Capsulorhexis Forceps 10 cm
vaulted shanks 11 mm sharp tips
Masket Kapsulorhexispinzette 10 cm
gewölbte Schenkel 11 mm spitze Plattform

- 01-10-460-1 Stainless Steel
- 01-10-460-2 **TITANIUM**
- 01-10-460-3 BLACK SERIES




Kraff Capsulorhexis Forceps 10 cm
 thin shanks 11 mm sharp tips
Kraff Kapsulorhexispinzette 10 cm
 feine Schenkel 11 mm spitze Plattform

- 01-10-465-1 Stainless Steel
- 01-10-465-2 TITANIUM
- 01-10-465-3 BLACK SERIES



Kraff Capsulorhexis Forceps 10 cm
 vaulted shanks 11 mm sharp tips
Kraff Kapsulorhexispinzette 10 cm
 gewölbte Schenkel 11 mm spitze Plattform

- 01-10-470-1 Stainless Steel
- 01-10-470-2 TITANIUM
- 01-10-470-3 BLACK SERIES



Giannetti Capsulorhexis MICS Forceps
Giannetti Kapsulorhexispinzette

- 01-10-472-1 1.5 mm incision
- 01-10-472-3 BLACK SERIES
- 01-10-473-1 2.2 mm incision
- 01-10-473-3 BLACK SERIES



Spaleck IOL inserting Forceps for silicone and acryl lens
Spaleck IOL einsetzpinzette für Silikon und Acryl Linsen

- 01-10-475-1 Stainless Steel
- 01-10-475-2 TITANIUM
- 01-10-475-3 BLACK SERIES



Livernois IOL Folding Forceps
Livernois IOL Faltpinzette

01-10-480-1
01-10-480-2 TITANIUM
01-10-480-3 BLACK SERIES



Nordan IOL Folding Forceps 12 cm, with lock
Nordan Linsenvorfaltpinzette 12 cm mit Sperre

01-10-485-1 Stainless Steel
01-10-485-2 TITANIUM



Doden Lens Implantation Forceps 12 cm
Doden Implantationspinzette 12 cm

01-10-487-1 Stainless Steel
01-10-487-3 BLACK SERIES



Faulkner IOL Holding Forceps 11 cm for foldable IOL's
Faulkner IOL Linsenhaltepinzette 11 cm für faltbare IOL's

01-10-490-1
01-10-490-3 BLACK SERIES



McDonald IOL Inserting Forceps 10 cm
cross-action for foldable IOL's
McDonald Linseneinsetzpinzette 10 cm
Kreuzaktion für faltbare IOL's

- 01-10-495-1 McDonald I
- 01-10-495-2 McDonald I TITANIUM
- 01-10-500-1 McDonald II
- 01-10-500-2 McDonald II TITANIUM



IOL Inserting Forceps II 12 cm for foldable IOL's
IOL Linseneinsetzpinzette II 12 cm für faltbare IOL's

- 01-10-510-2 with lock TITANIUM



IOL Inserting Forceps III 12 cm for foldable IOL's
IOL Linseneinsetzpinzette III 12 cm für faltbare IOL's

- 01-10-515-2 with lock TITANIUM



Bechert-McPherson Lens Inserting Forceps 11 cm
Bechert-McPherson Linsenhaltepinzette 11 cm

- 01-10-520-1 Plattform 11 mm
- 01-10-520-3 BLACK SERIES




Blaydes Lens Inserting Forceps 11 cm
Blaydes Linseneinsetzpinzette 11 cm

01-10-525-2 TITANIUM



Clayman Lens Implantation Forceps 12 cm tips
to grasp IOL haptics
Clayman Linseneinsetzpinzette 12 cm feine Enden
zum greifen der Haptiken

01-10-540-2 TITANIUM



Rappazzo IOL & Foreign Body Forceps 60°
Rappazzo IOL & Fremdkörperpinzette

01-10-550-1




IOL One-Hand Injector for IOL's
01-10-565-2




Troutman Superior Rectus Forceps 10 cm angled
Troutman Zügelnahtpinzette 10 cm abgewinkelt

01-10-575-1 1x2 C-teeth
01-10-575-3 BLACK SERIES




Jameson Muscle Forceps 10 cm
with lock 6 teeth - adult size
Jameson Muskelpinzette 10 cm mit Schloss
6 Zähne für Erwachsene

01-10-580-1 left
01-10-585-1 right



Jameson Muscle Forceps 10 cm
with lock 4 teeth - child size
Jameson Muskelpinzette 10 cm mit Schloss
4 Zähne für Kinder

01-10-590-1 left
01-10-595-1 right



Prince Muscle Forceps 10 cm
with lock 4 x 8 teeth - adult size
Prince Muskelpinzette 10 cm mit Schloss
4 x 8 Zähne für Erwachsene

01-10-600-1 left
01-10-605-1 right



Bangerter Muscle Clamp 10 cm with push button
10 teeth very delicate
Bangerter Muskelklemme 10 cm mit Drückknopf
10 Zähne sehr fein

01-10-615-1



Berke Ptosis Forceps 11 cm with lock
Berke Ptosispinzette 11 cm mit Schloss

01-10-620-1 20 mm

01-10-625-1 27 mm



Desmarres Chalazion Forceps 9 cm oval
Desmarres Chalazionspinzette 9 cm oval

- 01-10-630-1** small, 20 mm
- 01-10-632-1** medium, 25 mm
- 01-10-635-1** large, 30 mm



Lambert Chalazion Forceps 9 cm round
Lambert Chalazionspinzette 9 cm rund

- 01-10-640-1** small
- 01-10-640-2** TITANIUM
- 01-10-642-1** standard
- 01-10-642-2** TITANIUM



Heath Chalazion Forceps 9 cm oval
Heath Chalazionspinzette 9 cm oval

- 01-10-645-1** standard
- 01-10-645-2** TITANIUM



Lordan Chalazion Forceps 9 cm
Lordan Chalazionspinzette 9 cm

- 01-10-650-1**




Snellen Entropion Forceps 10 cm
Snellen Entropiumpinzette 10 cm

01-10-655-1 left
01-10-660-1 right



Barraquer Colibri Foreign Body Forceps 7 cm
Barraquer Kolibri Fremdkörperpinzette 7 cm

01-10-670-1 with Platform delicate
01-10-670-3 BLACK SERIES



JF Jewelers Forceps
JF Juwelierpinzette

01-10-685-1 #1 standard
01-10-690-1 #3c fine
01-10-695-1 #3 standard



JF Jewelers Forceps
JF Juwelierpinzette

01-10-700-1 #4 delicate
01-10-705-1 #5 extra delicate
01-10-705-2 #5 extra delicate **TITANIUM**



JF Suture Removal Forceps #5 angled extra delicate
JF Fadenlöspinzette #5 abgewinkelt sehr Fein

01-10-710-1
01-10-710-2 **TITANIUM**



JF Suture Removal Forceps #7 curved
JF Fadenlöspinzette #7 Sichelform

01-10-715-1
01-10-715-2 TITANIUM




JF Suture Removal Forceps #5a
JF Fadenlöspinzette #5a

01-10-720-1



Barraquer Colibri Suture Removal Forceps 7 cm
Barraquer Colibri Fadenlöspinzette 7 cm

01-10-730-1 Tying Platform



Bonaccolto Conjunctiva Forceps 10 cm
longitudinally and cross serrated jaws
Bonaccolto Bindehautpinzette 10 cm
längs / kreuzgerieft

01-10-735-1 fine
01-10-740-1 standard



Nugent Utility Forceps 10 cm angled
Nugent Mehrzweckpinzette 10 cm abgewinkelt

01-10-745-1 smooth jaws
01-10-750-1 criss-cross serrated jaws



Watzke Sleeve Spreading Forceps 12 cm cross-action
Watzke Ärmelspreizpinzette 12 cm Kreuzaktion

01-10-755-1



Carter Sphere Introducer and Holder 12.5 cm
Carter E nukleationsplomben Einführer & Halter 12.5 cm

01-10-770-1



Micro Scissors Sideways angled
Mikro Scheren seitlich abgewinkelt

- 01-11-010-1** Length of Blade 8 mm, sharp
- 01-11-010-3** BLACK SERIES
- 01-11-011-1** Length of Blade 8 mm, sharp with Ball
- 01-11-011-3** BLACK SERIES
- 01-11-012-1** Length of Blade 8 mm, blunt
- 01-11-012-3** BLACK SERIES
- 01-11-013-1** Length of Blade 8 mm, blunt with Ball
- 01-11-013-3** BLACK SERIES
- 01-11-014-1** Length of Blade 12 mm, sharp
- 01-11-014-3** BLACK SERIES
- 01-11-015-1** Length of Blade 12 mm, sharp with Ball
- 01-11-015-3** BLACK SERIES
- 01-11-016-1** Length of Blade 12 mm, blunt
- 01-11-016-3** BLACK SERIES
- 01-11-017-1** Length of Blade 12 mm, blunt with Ball
- 01-11-017-3** BLACK SERIES
- 01-11-018-1** Length of Blade 14 mm, sharp
- 01-11-018-3** BLACK SERIES
- 01-11-019-1** Length of Blade 14 mm, sharp with Ball
- 01-11-019-3** BLACK SERIES
- 01-11-020-1** Length of Blade 14 mm, blunt
- 01-11-020-3** BLACK SERIES
- 01-11-021-1** Length of Blade 14 mm, blunt with Ball
- 01-11-021-3** BLACK SERIES



Micro Scissors Round Handle

Mikro Scheren Rundgriff

01-11-025-1	9 cm, Blade 10 mm, straight, sharp
01-11-025-3	BLACK SERIES
01-11-026-1	9 cm, Blade 10 mm, straight, Vannas - sharp
01-11-026-3	BLACK SERIES
01-11-027-1	9 cm, Blade 10 mm, curved, sharp
01-11-027-3	BLACK SERIES
01-11-028-1	9 cm, Blade 10 mm, curved, Vannas - sharp
01-11-028-3	BLACK SERIES
01-11-029-1	9 cm, Blade 10 mm, straight, blunt
01-11-029-3	BLACK SERIES
01-11-030-1	9 cm, Blade 10 mm, curved, blunt
01-11-030-3	BLACK SERIES
01-11-031-1	12 cm, Blade 10 mm, straight, sharp
01-11-031-3	BLACK SERIES
01-11-032-1	12 cm, Blade 10 mm, straight, Vannas - sharp
01-11-032-3	BLACK SERIES
01-11-033-1	12 cm, Blade 10 mm, curved, sharp
01-11-033-3	BLACK SERIES
01-11-034-1	12 cm, Blade 10 mm, curved, Vannas - sharp
01-11-034-3	BLACK SERIES
01-11-035-1	12 cm, Blade 10 mm, straight, blunt
01-11-035-3	BLACK SERIES
01-11-036-1	12 cm, Blade 10 mm, curved, blunt
01-11-036-3	BLACK SERIES



Barraquer Iris Scissors 5.5 cm
Barraquer Iridektomieschere 5.5 cm

- 01-11-100-1** Length of Blade 7 mm blunt/blunt
- 01-11-110-1** Length of Blade 5 mm blunt/blunt





Wecker Iris Scissors 11 cm
Wecker Iridektomieschere 11 cm

- 01-11-120-1** Length of Blade 11 mm sharp/sharp
- 01-11-121-1** Length of Blade 7 mm sharp/sharp
- 01-11-130-1** Length of Blade 11 mm sharp/blunt
- 01-11-131-1** Length of Blade 7 mm sharp/blunt
- 01-11-140-1** Length of Blade 11 mm blunt/blunt
- 01-11-141-1** Length of Blade 7 mm blunt/blunt



Castroviejo Anterior Synechia Scissors 10.5 cm
Castroviejo Augenschere 10.5 cm
für die vorderen Synechien

- 01-11-145-1** delicate
- 01-11-145-3** BLACK SERIES



Vannas Capsulotomy Scissors 8 cm

Vannas Kapsulotomieschere 8 cm


- 01-11-148-1** Length of Blade 5 mm straight, sharp
- 01-11-148-3** BLACK SERIES
- 01-11-149-1** Length of Blade 5 mm straight, blunt
- 01-11-149-3** BLACK SERIES
- 01-11-150-1** Length of Blade 6 mm straight, sharp
- 01-11-150-3** BLACK SERIES
- 01-11-151-1** Length of Blade 6 mm straight, blunt
- 01-11-151-3** BLACK SERIES
- 01-11-158-1** Length of Blade 5 mm curved, sharp
- 01-11-158-3** BLACK SERIES
- 01-11-159-1** Length of Blade 5 mm curved, blunt
- 01-11-159-3** BLACK SERIES
- 01-11-160-1** Length of Blade 6 mm curved, sharp
- 01-11-160-3** BLACK SERIES
- 01-11-161-1** Length of Blade 6 mm curved, blunt
- 01-11-161-3** BLACK SERIES
- 01-11-162-1** Length of Blade 5 mm angled to side
- 01-11-162-3** BLACK SERIES
- 01-11-163-1** Length of Blade 6 mm angled to side
- 01-11-163-3** BLACK SERIES
- 01-11-164-1** Length of Blade 10 mm angled to side
- 01-11-164-3** BLACK SERIES



Vannas Capsulotomy Scissors 8 cm


Vannas Kapsulotomieschere 8 cm

- 01-11-170-1** Length of Blade 6 mm angled forward
- 01-11-170-3** BLACK SERIES





Mod.Tübingen Capsulotomy Scissors 8.5 cm, 10mm Blade
Mod.Tübingen Kapsulotomieschere 8.5 cm, 10mm Schneide

- 01-11-180-1 straight / sharp
- 01-11-180-3 BLACK SERIES
- 01-11-181-1 straight / blunt
- 01-11-181-3 BLACK SERIES
- 01-11-190-1 curved / sharp
- 01-11-190-3 BLACK SERIES
- 01-11-191-1 curved / blunt
- 01-11-191-3 BLACK SERIES




Mod.Tübingen Capsulotomy Scissors 8.5 cm
Mod.Tübingen Kapsulotomieschere 8.5 cm

- 01-11-200-1 Length of Blade 10 mm angled forward
- 01-11-200-3 BLACK SERIES




Gills-Welsh Capsulotomy Scissors 8.5 cm
Gills-Welsh Kapsulotomieschere 8.5 cm

- 01-11-210-1 Length of Blade 12 mm straight
- 01-11-210-3 BLACK SERIES
- 01-11-220-1 Length of Blade 12 mm curved
- 01-11-220-3 BLACK SERIES



Gills-Welsh Capsulotomy Scissors 8.5 cm
Gills-Welsh Kapsulotomieschere 8.5 cm

- 01-11-230-1 Length of Blade 12 mm angled forward
- 01-11-230-3 BLACK SERIES
- 01-11-240-1 Length of Blade 12 mm angled to side
- 01-11-240-3 BLACK SERIES




Vannas-Cohan Capsulotomy Scissors 6.5 cm
Vannas-Cohan Kapsulotomieschere 6.5 cm

01-11-242-1 straight
01-11-242-3 BLACK SERIES
01-11-243-1 curved
01-11-243-3 BLACK SERIES



Cornic Capsulotomy Scissors 10 cm
Cornic Kapsulotomieschere 10 cm

01-11-250-1 Extra Thin Blade 10 mm curved
01-11-250-3 BLACK SERIES




Micro Capsulotomy Scissors 12 cm
Micro Kapsulotomieschere 12 cm

01-11-252-1 7 mm Blade, straight
01-11-252-3 BLACK SERIES
01-11-253-1 7 mm Blade, curved
01-11-253-3 BLACK SERIES



Vannas Scissors 15 cm, 11 mm Blade
Vannas Schere 15 cm, 11 mm Blatt

01-11-254-1 straight
01-11-254-3 BLACK SERIES
01-11-256-1 curved
01-11-256-3 BLACK SERIES




Castroviejo Scissors strong curved 10 cm
Castroviejo Schere stark gebogen 10 cm

- 01-11-257-1 Flathandle
- 01-11-257-3 BLACK SERIES
- 01-11-258-1 Roundhandle
- 01-11-258-3 BLACK SERIES




Castroviejo Corneoscleral Scissors 10 cm
Castroviejo Hornhautschere 10 cm

- 01-11-260-1 right 10 mm blades
- 01-11-260-3 BLACK SERIES
- 01-11-261-1 right 10 mm blades with lock
- 01-11-262-1 right 7 mm blades
- 01-11-262-3 BLACK SERIES
- 01-11-263-1 right 10 mm blades, width Handle
- 01-11-263-3 BLACK SERIES
- 01-11-264-1 right 13 mm blades, width Handle
- 01-11-264-3 BLACK SERIES
- 01-11-270-1 left 10 mm blades
- 01-11-270-3 BLACK SERIES
- 01-11-271-1 left 10 mm blades with lock
- 01-11-272-1 left 7 mm blades
- 01-11-272-3 BLACK SERIES
- 01-11-273-1 left 10 mm blades, width Handle
- 01-11-273-3 BLACK SERIES
- 01-11-274-1 left 13 mm blades, width Handle
- 01-11-274-3 BLACK SERIES




Castroviejo Universal Corneal Scissors 11 cm
Castroviejo Universal Hornhautschere 11 cm

- 01-11-280-1 curved - blunt tips
- 01-11-280-3 BLACK SERIES
- 01-11-281-1 straight - blunt tips
- 01-11-281-3 BLACK SERIES
- 01-11-282-1 curved - sharp tips
- 01-11-282-3 BLACK SERIES
- 01-11-283-1 straight - sharp tips
- 01-11-283-3 BLACK SERIES
- 01-11-284-1 curved - sharp / blunt tips
- 01-11-284-3 BLACK SERIES
- 01-11-285-1 straight - sharp / blunt tips
- 01-11-285-3 BLACK SERIES



Katzin Corneal Transplant Scissors 10 cm
Katzin Keratoplastikschiere 10 cm

- 01-11-290-1 right 7 mm blades
- 01-11-290-3 BLACK SERIES
- 01-11-291-1 right 7 mm blades with lock
- 01-11-292-1 right 10mm blades
- 01-11-292-3 BLACK SERIES
- 01-11-293-1 right 13 mm blades
- 01-11-293-3 BLACK SERIES
- 01-11-300-1 left 7 mm blades
- 01-11-300-3 BLACK SERIES
- 01-11-301-1 left 7 mm blades with lock
- 01-11-302-1 left 10 mm blades
- 01-11-302-3 BLACK SERIES
- 01-11-303-1 left 13 mm blades
- 01-11-303-3 BLACK SERIES



Troutman Corneal Transplant Scissors 10 cm
Troutman Keratoplastikschiere 10 cm

- 01-11-305-1 right 7 mm blades
- 01-11-305-3 BLACK SERIES
- 01-11-306-1 right 7 mm blades with lock
- 01-11-307-1 left 7 mm blades
- 01-11-307-3 BLACK SERIES
- 01-11-308-1 left 7 mm blades with lock



Castroviejo Scissors 11 cm, strongly curved
Castroviejo Federschiere 11 cm starkgebogen

- 01-11-310-1 sharp tips
- 01-11-310-3 BLACK SERIES



Castroviejo Universal Corneal Scissors 11 cm

Castroviejo Universal Keratomieschere 11 cm

01-11-320-1 strongly curved - blunt tips

01-11-320-3 BLACK SERIES



Mackensen Keratoplasty Scissors 8 cm

Mackensen Keratoplastikscherer 8 cm

01-11-321-1 curved

01-11-321-3 BLACK SERIES



Castroviejo Keratoplasty Scissors 9 cm, angled to side

Castroviejo Keratoplastikscherer 9 cm kniegewinkelt

01-11-330-1 standard blades

01-11-330-3 BLACK SERIES

01-11-340-1 small blades

01-11-340-3 BLACK SERIES



Westcott-Shepard Tenotomy Scissors 11 cm, curved, blunt serrated

Westcott-Shepard Sehnenscherer 11 cm, gebogen, stumpf, gezahnt

01-11-345-1 right

01-11-345-3 BLACK SERIES

01-11-346-1 left

01-11-346-3 BLACK SERIES



Standard Handle



Width Handle



Westcott Tenotomy Scissors
Westcott Sehnenscherer

- 01-11-350-1 11 cm, curved, blunt, medium blades
- 01-11-350-3 BLACK SERIES
- 01-11-360-1 11 cm, curved, blunt, standard blades right Hand
- 01-11-360-2 11 cm, curved, blunt, standard blades right Hand **TITANIUM**
- 01-11-360-3 BLACK SERIES
- 01-11-370-1 11 cm, curved, blunt standard blades left Hand
- 01-11-370-3 BLACK SERIES
- 01-11-371-1 13 cm width Handle, curved, blunt standard blades right Hand
- 01-11-371-3 BLACK SERIES
- 01-11-372-1 13 cm width Handle, curved, blunt standard blades left Hand
- 01-11-372-3 BLACK SERIES
- 01-11-373-1 11 cm, straight, blunt standard blades right Hand
- 01-11-373-3 BLACK SERIES
- 01-11-374-1 11 cm, straight, blunt standard blades left Hand
- 01-11-374-3 BLACK SERIES
- 01-11-375-1 13 cm width Handle, straight, blunt standard blades right Hand
- 01-11-375-3 BLACK SERIES
- 01-11-376-1 13 cm width Handle, straight, blunt standard blades left Hand
- 01-11-376-3 BLACK SERIES

Fine Stitch Scissors 11 cm
Fine Fadenscherer 11 cm

- 01-11-380-1 straight, sharp
- 01-11-380-3 BLACK SERIES
- 01-11-382-1 straight, sharp / blunt
- 01-11-382-3 BLACK SERIES
- 01-11-384-1 straight, blunt
- 01-11-384-3 BLACK SERIES
- 01-11-390-1 curved, sharp
- 01-11-390-3 BLACK SERIES
- 01-11-392-1 curved, sharp / blunt
- 01-11-392-3 BLACK SERIES
- 01-11-394-1 curved, blunt
- 01-11-394-3 BLACK SERIES



MSS Micro Iris Scissors 10 cm, 9 mm Blade
MSS Mikro Fadenschere 10 cm, 9 mm Blatt

- 01-11-400-1** straight, sharp
- 01-11-400-3** BLACK SERIES
- 01-11-402-1** straight, sharp / blunt
- 01-11-402-3** BLACK SERIES
- 01-11-404-1** straight, blunt
- 01-11-404-3** BLACK SERIES
- 01-11-410-1** curved, sharp
- 01-11-410-3** BLACK SERIES
- 01-11-412-1** curved, sharp / blunt
- 01-11-412-3** BLACK SERIES
- 01-11-414-1** curved, blunt
- 01-11-414-3** BLACK SERIES
- 01-11-415-1** strong curved, sharp
- 01-11-415-3** BLACK SERIES
- 01-11-416-1** strong curved, blunt
- 01-11-416-3** BLACK SERIES



Westcott Stitch Scissors
Westcott Fadenschere


- 01-11-420-1** 11 cm curved, sharp tips
- 01-11-420-3** BLACK SERIES
- 01-11-421-1** 11 cm straight, sharp tips
- 01-11-421-3** BLACK SERIES
- 01-11-422-1** 13 cm width Handle, curved, sharp tips
- 01-11-422-3** BLACK SERIES
- 01-11-423-1** 13 cm width Handle, straight, sharp tips
- 01-11-423-3** BLACK SERIES



McPherson-Westcott Scissors

McPherson-Westcott Schere

- 01-11-428-1** 10 cm, sharp tips - straight
- 01-11-428-3** BLACK SERIES
- 01-11-429-1** 10 cm, blunt tips - straight
- 01-11-429-3** BLACK SERIES
- 01-11-430-1** 10 cm, sharp tips - curved
- 01-11-430-3** BLACK SERIES
- 01-11-431-1** 10 cm, blunt tips - curved
- 01-11-431-3** BLACK SERIES




- 01-11-432-1** 13 cm, width Handle, sharp tips - straight
- 01-11-432-3** BLACK SERIES
- 01-11-433-1** 13 cm, width Handle, blunt tips - straight
- 01-11-433-3** BLACK SERIES
- 01-11-434-1** 13 cm, width Handle, sharp tips - curved
- 01-11-434-3** BLACK SERIES
- 01-11-435-1** 13 cm, width Handle, blunt tips - curved
- 01-11-435-3** BLACK SERIES



BROSS IOL Cutting Scissors, blunt tips, serrated lower blade

- 01-11-437-1** straight
- 01-11-437-3** BLACK SERIES



- 01-11-438-1** angled to side
- 01-11-438-3** BLACK SERIES



Mod.Heidelberg Strabismus Scissors 9 cm
Mod.Heidelberg Strabismusschere 9 cm

- 01-11-440-1 short blades - straight
- 01-11-440-3 BLACK SERIES
- 01-11-450-1 short blades - curved
- 01-11-450-3 BLACK SERIES



Knapp Strabismus Scissors 10.5 cm
Knapp Strabismusschere 10.5 cm

- 01-11-460-1 straight
- 01-11-460-3 BLACK SERIES
- 01-11-460-6 Tungsten Carbide Line
- 01-11-460-7 SuperCut
- 01-11-470-1 curved
- 01-11-470-3 BLACK SERIES
- 01-11-470-6 Tungsten Carbide Line
- 01-11-470-7 SuperCut



Mod.Bonn Eye Scissors 9 cm
Mod.Bonn Augenschere 9 cm

- 01-11-490-1 straight - sharp tips
- 01-11-490-2 TITANIUM
- 01-11-490-3 BLACK SERIES
- 01-11-490-7 SuperCut
- 01-11-500-1 straight - blunt tips
- 01-11-500-3 BLACK SERIES
- 01-11-510-1 curved - sharp tips
- 01-11-510-2 TITANIUM
- 01-11-510-3 BLACK SERIES
- 01-11-510-7 SuperCut
- 01-11-520-1 curved - blunt tips
- 01-11-520-3 BLACK SERIES



ESL Eye Scissors 10 cm

with large rings - sharp tips

ESL Augenschere mit großen Ringen 10 cm ,
spitze Spitzen

- 01-11-530-1** straight
- 01-11-530-3** BLACK SERIES
- 01-11-540-1** curved
- 01-11-540-3** BLACK SERIES



SL Eye Scissors 11 cm for lids, sharp tips

SL Augenschere 11 cm für Lider, spitze Spitzen

- 01-11-550-1** straight
- 01-11-550-2** TITANIUM
- 01-11-550-3** BLACK SERIES
- 01-11-550-6** Tungsten Carbide
- 01-11-550-7** SuperCut
- 01-11-560-1** curved
- 01-11-560-2** TITANIUM
- 01-11-560-3** BLACK SERIES
- 01-11-560-6** Tungsten Carbide
- 01-11-560-7** SuperCut



SF Eye Scissors extra flat 10.5 cm, sharp tips

SF Augenschere extra flach 10.5 cm, spitze Spitzen

- 01-11-570-1** straight
- 01-11-570-3** BLACK SERIES
- 01-11-580-1** curved
- 01-11-580-3** BLACK SERIES



Suture Stitch Scissors 10.5 cm, Needle points

Suture Fadenschere 10.5 cm, Nadelspitz

- 01-11-590-1 curved
- 01-11-590-2 TITANIUM
- 01-11-590-3 BLACK SERIES
- 01-11-600-1 straight
- 01-11-600-3 BLACK SERIES



Stevens Tenotomy Scissors 10.5 cm blunt tips

Stevens Sehnenschere 10.5 cm, stumpfe Spitzen

- 01-11-610-1 straight
- 01-11-610-3 BLACK SERIES
- 01-11-610-7 SuperCut
- 01-11-620-1 curved
- 01-11-620-3 BLACK SERIES
- 01-11-620-7 SuperCut



Stevens Tenotomy Scissors 11.5 cm blunt tips

Stevens Sehnenschere 11.5 cm, stumpfe Spitzen

- 01-11-630-1 straight
- 01-11-630-3 BLACK SERIES
- 01-11-630-6 Tungsten Carbide
- 01-11-630-7 SuperCut
- 01-11-640-1 curved
- 01-11-640-3 BLACK SERIES
- 01-11-640-6 Tungsten Carbide
- 01-11-640-7 SuperCut



Spencer Stitch Scissors, sharp / blunt

Spencer Ligaturschere, spitz / stumpf

01-11-650-1	9 cm, straight
01-11-650-3	BLACK SERIES
01-11-651-1	11 cm, straight
01-11-651-3	BLACK SERIES
01-11-652-1	12 cm, straight
01-11-652-3	BLACK SERIES
01-11-653-1	9 cm, curved
01-11-653-3	BLACK SERIES
01-11-654-1	12 cm, curved
01-11-654-3	BLACK SERIES
01-11-655-1	11 cm, angled to side
01-11-655-3	BLACK SERIES



ES Eye Scissors 11 cm, angled to side, sharp tips

ES Augenschere 11 cm, kniegewinkelt, spitze Spitzen

01-11-660-1	
01-11-660-3	BLACK SERIES
01-11-660-7	SuperCut
01-11-661-1	with ball tip
01-11-661-3	BLACK SERIES
01-11-662-1	curved to side
01-11-662-3	BLACK SERIES



ES Enucleation Scissors 13 cm curved
ES Eukliationsschere 13 cm gebogen

01-11-670-1
01-11-670-3 BLACK SERIES
01-11-670-7 SuperCut
01-11-675-1 medium curved
01-11-675-3 BLACK SERIES
01-11-675-7 SuperCut



ES Enucleation Scissors 13 cm strong curved
ES Eukliationsschere 13 cm stark gebogen

01-11-680-1
01-11-680-3 BLACK SERIES



Kalt Needle Holder 14 cm with lock

Kalt Nadelhalter 14 cm mit Sperre

- 01-12-100-1** straight, serrated jaws
- 01-12-105-1** curved, serrated jaws
- 01-12-110-1** straight, smooth jaws
- 01-12-115-1** curved, smooth jaws



Barraquer-Troutman Needle Holder 10 cm

Barraquer-Troutman Nadelhalter 10 cm

- 01-12-120-1** wo/L, straight, very delicate
- 01-12-120-3** BLACK SERIES
- 01-12-125-1** wo/L, curved, very delicate
- 01-12-125-3** BLACK SERIES
- 01-12-130-1** w/L, straight, very delicate
- 01-12-130-3** BLACK SERIES
- 01-12-135-1** w/L, curved, very delicate
- 01-12-135-3** BLACK SERIES



Barraquer Needle Holder, straight, without lock

Barraquer Nadelhalter, gerade, ohne Sperre

- 01-12-140-1** 11 cm delicate
- 01-12-140-2** TITANIUM
- 01-12-140-3** BLACK SERIES
- 01-12-145-1** 11 cm, very delicate
- 01-12-145-3** BLACK SERIES
- 01-12-150-1** 13 cm, delicate
- 01-12-150-3** BLACK SERIES
- 01-12-155-1** 13 cm, very delicate
- 01-12-155-3** BLACK SERIES




Barraquer Needle Holder curved, without lock
Barraquer Nadelhalter, gebogen, ohne Sperre

- 01-12-160-1 11 cm, delicate
- 01-12-160-2 TITANIUM
- 01-12-160-3 BLACK SERIES
- 01-12-165-1 11 cm, very delicate
- 01-12-165-3 BLACK SERIES
- 01-12-169-1 12 cm, very delicate
- 01-12-169-3 BLACK SERIES
- 01-12-170-1 13 cm, delicate
- 01-12-170-2 TITANIUM
- 01-12-170-3 BLACK SERIES
- 01-12-175-1 13 cm, very delicate
- 01-12-175-3 BLACK SERIES



Barraquer Needle Holder, straight, with lock
Barraquer Nadelhalter, gerade, mit Sperre

- 01-12-180-1 11 cm, delicate
- 01-12-180-3 BLACK SERIES
- 01-12-185-1 11 cm, very delicate
- 01-12-185-3 BLACK SERIES
- 01-12-190-1 13 cm delicate
- 01-12-190-3 BLACK SERIES
- 01-12-195-1 13 cm, very delicate
- 01-12-195-3 BLACK SERIES



Barraquer Needle Holder, curved, with lock
Barraquer Nadelhalter, gebogen, mit Sperre

- 01-12-200-1 11 cm delicate
- 01-12-200-3 BLACK SERIES
- 01-12-205-1 11 cm very delicate
- 01-12-205-3 BLACK SERIES
- 01-12-210-1 13 cm delicate
- 01-12-210-2 TITANIUM
- 01-12-210-3 BLACK SERIES
- 01-12-215-1 13 cm, very delicate
- 01-12-215-3 BLACK SERIES



Castroviejo Needle Holder 10 cm straight, very delicate
Castroviejo Nadelhalter 10 cm gerade, sehr fein

- 01-12-220-1 without lock
- 01-12-220-3 BLACK SERIES
- 01-12-225-1 with lock
- 01-12-225-3 BLACK SERIES



Castroviejo Needle Holder 10 cm curved, very delicate
Castroviejo Nadelhalter 10 cm gebogen, sehr fein

- 01-12-230-1 without lock
- 01-12-230-3 BLACK SERIES
- 01-12-235-1 with lock
- 01-12-235-3 BLACK SERIES



Castroviejo Standard Needle Holder 14 cm, straight
Castroviejo Standard Nadelhalter 14 cm, gerade

- 01-12-240-1 without lock
- 01-12-240-2 TITANIUM
- 01-12-240-3 BLACK SERIES
- 01-12-245-1 with lock
- 01-12-245-2 TITANIUM
- 01-12-245-3 BLACK SERIES



Castroviejo Standard Needle Holder 14 cm, curved
Castroviejo Standard Nadelhalter 14 cm, gebogen


- 01-12-250-1 without lock
- 01-12-250-2 TITANIUM
- 01-12-250-3 BLACK SERIES
- 01-12-255-1 with lock
- 01-12-255-2 TITANIUM
- 01-12-255-3 BLACK SERIES



Castroviejo Needle Holder 14 cm straight, delicate

Castroviejo Nadelhalter 14 cm gerade, fein

- 01-12-260-1 without lock
- 01-12-260-2 TITANIUM
- 01-12-260-3 BLACK SERIES
- 01-12-265-1 with lock
- 01-12-265-2 TITANIUM
- 01-12-265-3 BLACK SERIES



Castroviejo Needle Holder 14 cm curved, delicate

Castroviejo Nadelhalter 14 cm gebogen, fein

- 01-12-270-1 without lock
- 01-12-270-2 TITANIUM
- 01-12-270-3 BLACK SERIES
- 01-12-275-1 with lock
- 01-12-275-2 TITANIUM
- 01-12-275-3 BLACK SERIES



Castroviejo-Jaffe Needle Holder 13 cm straight

Castroviejo-Jaffe Nadelhalter 13 cm gerade


- 01-12-280-1 without lock, delicate
- 01-12-280-3 BLACK SERIES
- 01-12-285-1 without lock, very delicate
- 01-12-285-3 BLACK SERIES
- 01-12-290-1 with lock, delicate
- 01-12-290-3 BLACK SERIES
- 01-12-295-1 with lock, very delicate
- 01-12-295-3 BLACK SERIES




Castroviejo-Jaffe Needle Holder 13 cm curved

Castroviejo-Jaffe Nadelhalter 13 cm gebogen


- 01-12-300-1 without lock, delicate
- 01-12-300-3 BLACK SERIES
- 01-12-305-1 without lock, very delicate
- 01-12-305-3 BLACK SERIES
- 01-12-310-1 with lock, delicate
- 01-12-310-3 BLACK SERIES
- 01-12-315-1 with lock, very delicate
- 01-12-315-3 BLACK SERIES



TC = Tungsten-Carbide



smooth profile



micro profile

Castroviejo TC Needle Holder 14 cm, straight
Castroviejo TC Nadelhalter 14 cm, gerade

- 01-12-500-6** with lock, smooth profile
- 01-12-505-6** without lock, smooth profile
- 01-12-510-6** with lock, micro profile
- 01-12-515-6** without lock, micro profile



Castroviejo TC Needle Holder 14 cm, curved
Castroviejo TC Nadelhalter 14 cm, gebogen

- 01-12-520-6** with lock, smooth profile
- 01-12-525-6** without lock, smooth profile
- 01-12-530-6** with lock, micro profile
- 01-12-535-6** without lock, micro profile



Barraquer TC Needle Holder 13 cm
Barraquer TC Nadelhalter 13 cm

- 01-12-540-6** straight, without lock, smooth profile
- 01-12-545-6** straight, with lock, smooth profile
- 01-12-550-6** straight, without lock, micro profile
- 01-12-555-6** straight, with lock, micro profile
- 01-12-560-6** curved, without lock, smooth profile
- 01-12-565-6** curved, with lock, smooth profile
- 01-12-570-6** curved, without lock, micro profile
- 01-12-575-6** curved, with lock, micro profile



Bross-Castroviejo TC Needle Holder, 12 cm, with lock
Bross-Castroviejo TC Nadelhalter, 12 cm, mit Sperre

- 01-12-580-6** curved, smooth profile
- 01-12-585-6** straight, smooth profile
- 01-12-590-6** curved, micro profile
- 01-12-595-6** straight, micro profile





Vitrectomy Instrument



Vitrectomy Instrument with irrigation




Vitrectomy Membrane Peeler Forceps
20G/0.9 mm 45° angled
01-13-105-1 without irrigation
01-13-110-1 with irrigation



Vitreous Body Scissors
20G/0.9 mm straight
01-13-145-1 without irrigation
01-13-150-1 with irrigation



Vitrectomy Membrane Peeler Forceps
20G/0.9 mm 90° angled
01-13-115-1 without irrigation
01-13-120-1 with irrigation



Vitreous Body Scissors
20G/0.3 mm , 30° angled
01-13-155-1 without irrigation
01-13-160-1 with irrigation




Vitrectomy Forceps
20G/0.9 mm straight serrated
01-13-125-1 without irrigation
01-13-130-1 with irrigation



Vitreous Body Scissors
20G/0.9 mm , 60° angled
01-13-165-1 without irrigation



Vitrectomy Forceps
20G/0.9 mm straight with grip tips
01-13-135-1 without irrigation
01-13-140-1 with irrigation



Squeeze handle instruments
for Vitreo-Retinal surgery 20G and 23G



Foreign Body Forceps
with cap jaws - straight
01-13-170-2 20G
01-13-172-2 23G



Membrane Peeler
smooth - 45°
01-13-202-2 20G
01-13-204-2 23G




Capsulorhexis Forceps - 2.5 mm jaws
01-13-174-2 20G, curved
01-13-176-2 23G, curved
01-13-178-2 20G, straight
01-13-180-2 23G, straight



Membrane Peeler
grooved - 45°
01-13-206-2 20G
01-13-208-2 23G




Gripping Forceps
2.5 mm - straight
01-13-182-2 20G
01-13-184-2 23G



Scissor - straight pointed tips
01-13-210-2 20G
01-13-212-2 23G



Gripping Forceps
2.0 mm smooth rounded tips - straight
01-13-186-2 20G
01-13-188-2 23G



Scissor - straight blunt tips
01-13-214-2 20G
01-13-216-2 23G



Forceps serrated - straight
01-13-190-2 20G
01-13-192-2 23G



Scissor - curved (R12)
to follow the contour of the retina
01-13-218-2 20G
01-13-220-2 23G



Forceps with 1x2 teeth - straight
01-13-194-2 20G
01-13-196-2 23G



Scissor - angled 45°
01-13-222-2 20G
01-13-224-2 23G



Forceps with gripping jaws - straight
01-13-198-2 20G
01-13-200-2 23G

RETINAL INSTRUMENTS 20-23-25 GAUGE



SAVE MONEY

The new system of detachable tips allows the surgeon to keep a reusable handpiece with the option to use a reusable or disposable tip, in different models.

- REUSABLE HANDPIECE
- REUSABLE TIPS
- DISPOSABLE TIPS

The Handpieces are reusable and have interchangeable disposable and reusable tips. Different tip colors help identify the instrument with ease. This system allows the surgeon to change only the tips after a 20-25ga surgery, saving money and time.



Reusable Tip



Disposable Tip



Squeeze Handle:
- High Constructive Technology
- Color Coded
- High Performance



Finger Lever Handle:
- High Performance
- Color Coded
- Low cost


Spear Swabs with handle sterile
Spear Saugtupfer mit Griff steril

01-14-100-1 Box 20 x 10 pieces



Universal Syringe 5 ml Luer Lock - heat resistant
Universal Spritze 5ml mit Luer Lock

01-14-120-1



Barraquer Sable Brush
Barraquer Operationspinsel

01-14-180-1



Skin Marking pens
Hautmarkierstifte

01-14-182-1



West Lacrimal Sac Chisel 16 cm straight 5 mm
West Tränensackflachmeißel 16cm gerade 5mm

01-14-190-1



West Lacrimal Sac Gouge 16 cm straight 5 mm
West Tränensackhohlmeißel 16cm gerade 5mm

01-14-200-1



Freer Lacrimal Sac Chisel 16 cm curved 4 mm
Freer Tränensackflachmeißel 16 cm gebogen 4 mm

01-14-210-1




Mead Mallet with interchangeable disks
Mead Hammer mit auswechselbaren Backen

01-14-220-1



Mead Spare Facings (pair) for Mallet
Mead Ersatzbacken (Paar) für Hammer

01-14-230-1



Freer Septum Elevator 18 cm sharp/blunt
Freer Septumelevator 18 cm scharf/stump

01-14-240-1



Joseph Periosteal Elevator 16 cm
Joseph Periostealelevator 16 cm

01-14-250-1



Carter Periosteal Elevator 21 cm
sharp spoon / blunt spatula
Carter Periostealelevator 21 cm
scharfer Löffel / stumpfer Spatel

01-14-260-1




Hartmann-Halle Nasal Speculum 15 cm
Hartmann-Halle Nasenspekulum 15 cm

01-14-270-1 Infant
01-14-280-1 Child
01-14-290-1 Adult



Descemet Membrane Punch **TITANIUM**
Descemet Membranstanze **TITANIUM**

01-14-300-2 Ø 1 mm, 0.75 x 0.5 mm deep
01-14-310-2 Ø 0.6 mm, 0.3 x 0.5 mm deep



Kelly Descemet's Membrane Punch Ø 1 mm
Kelly Descemet's Membranstanze Ø 1 mm

01-14-320-1




Luntz-Dodick Trabeculectomy Punch
Ø 1 mm, 0.5 mm deep
Luntz-Dodick Trabeculectomiestanze
Ø 1 mm, 0.5 mm

01-14-330-1 Steel
01-14-330-2 **TITANIUM**



Haitz Canaliculus Punch
Haitz Kanalikulus-Stanze

01-14-341-1 angled, 1.0 mm bite





Citelli Bone Nibbling Rongeur
Citelli Knochenstanze

01-14-380-1 3.0 mm



Kerrison Bone Nibbling Rongeur
Kerrison Knochenstanze

01-14-390-1 3.0 mm
01-14-395-1 3.5 mm
01-14-400-1 4.0 mm
01-14-405-1 4.5 mm
01-14-410-1 5.0 mm

Mini-Friedmann Bone Rongeur Forceps 12 cm
Mini-Friedmann Hohlmeißelzange 12 cm

01-14-420-1



Dieffenbach Serrefine straight
Dieffenbach Bulldogklemme gerade

01-14-430-1 6.0 cm

01-14-430-2 6.0 cm **TITANIUM**

01-14-440-1 3.5 cm

01-14-442-1 2.5 cm



Dieffenbach Serrefine curved
Dieffenbach Bulldogklemme gebogen

01-14-450-1 6.0 cm

01-14-450-2 6.0 cm **TITANIUM**


01-14-460-1 3.5 cm

01-14-462-1 2.5 cm



Hartmann Mosquito Forceps 10 cm
Hartmann Arterienklemme 10 cm

01-14-470-1 straight
01-14-470-2 straight TITANIUM
01-14-480-1 curved
01-14-480-2 curved TITANIUM




Halsted-Mosquito Artery Forceps straight
Halsted-Mosquito Arterienklemme, gerade

01-14-490-1 12,5 cm
01-14-490-2 12,5 cm TITANIUM
01-14-495-1 14.0 cm



Halsted-Mosquito Artery Forceps curved
Halsted-Mosquito Arterienklemme, gebogen

01-14-500-1 12 cm
01-14-500-2 12 cm TITANIUM
01-14-505-1 14 cm



Backhouse Towel Clamp
Backhouse Tuchklemme

01-14-510-1 9 cm
01-14-520-1 13 cm



Jones Towel Clamp
Jones Tuchklemme

01-14-530-1 5.5 cm
01-14-530-2 5.5 cm
01-14-540-1 8.5 cm

TITANIUM





Honan Balloon Complete

01-14-800-1



Adson Forceps, 12 cm
 Adson Pinzette, 12 cm

01-15-050-1




EU Connector

01-15-050-2



US Connector



Jeweler Forceps, straight, pointed tips
 Juwelierpinzette, gerade, spitz

01-15-100-1



EU Connector

01-15-100-2



US Connector

Jeweler Forceps, straight, 0.5 mm tips
 Juwelierpinzette, gerade, 0.5 mm spitze

01-15-102-1



EU Connector

01-15-102-2



US Connector



Jeweler Forceps, angled, pointed tips
 Juwelierpinzette, gewinkelt, spitz

01-15-104-1



EU Connector

01-15-104-2



US Connector

Jeweler Forceps, angled, 0.5 mm tips
 Juwelierpinzette, gewinkelt, 0.5 mm spitze

01-15-106-1



EU Connector

01-15-106-2



US Connector



Flachstecker
Flat Plug



2-Pin-Stecker
2-pin plug

	2- Banana	3 m	01-15-500-1	3 m	01-15-500-2
	* Non-compliant with EU standards	5 m	01-15-501-1	5 m	01-15-501-2
	ERBE	3 m	01-15-502-1	3 m	01-15-502-2
		5 m	01-15-503-1	5 m	01-15-503-2
	Martin, Berchtold, Aesculap GK 55, GK 60	3 m	01-15-504-1	3 m	01-15-504-2
		5 m	01-15-505-1	5 m	01-15-505-2
	Valleylab, Lamidey, EMC	3 m	01-15-506-1	3 m	01-15-506-2
		5 m	01-15-507-1	5 m	01-15-507-2
	Erbe, EMC / Dolley	3 m	01-15-508-1	3 m	01-15-508-2
		5 m	01-15-509-1		
	Codman CMC II	3 m	01-15-510-1	3 m	01-15-510-2
		5 m	01-15-511-1	5 m	01-15-511-2
	Valleylab autobipolar	3 m	01-15-512-1		

Ophthalmic Containers for Steam Sterilization

Sterilisation cases with silicone cushioning stripes, cover and bottom perforated with long-life textile filters and some of them also with single use paper filters. All cases are in Stainless Steel. With interchangeable identification tags (Colour and text at option).



Size 1
220 x 165 x 55 mm (75 mm)



Size 6
400 x 225 x 55 mm (75 mm)



Size 2
420 x 125 x 55 mm (75 mm)




Size 7
400 x 290 x 55 mm (75 mm)



Size 3
420 x 175 x 55 mm (75 mm)



Size 8
450 x 260 x 55 mm (75 mm)



Size 4
300 x 220 x 55 mm (75 mm)



Size 9
400 x 150 x 55 mm (75 mm)



Size 5
285 x 280 x 55 mm (75 mm)



Size 10
495 x 270 x 55 mm (75 mm)



Size

One Rack

Two Racks

Size		One Rack	Two Racks
1	AM	95-01-100-1 220 x 165 x 55 mm	95-01-102-1 220 x 165 x 75 mm
	IM	200 x 140 x 45 mm	200 x 140 x 65 mm
2	AM	95-01-104-1 420 x 125 x 55 mm	95-01-106-1 420 x 125 x 75 mm
	IM	400 x 100 x 45 mm	400 x 100 x 65 mm
3	AM	95-01-108-1 420 x 175 x 55 mm	95-01-110-1 420 x 175 x 75 mm
	IM	400 x 150 x 45 mm	400 x 150 x 65 mm
4	AM	95-01-112-1 300 x 220 x 55 mm	95-01-114-1 300 x 220 x 75 mm
	IM	280 x 195 x 45 mm	280 x 195 x 65 mm
5	AM	95-01-116-1 285 x 280 x 55 mm	95-01-118-1 285 x 280 x 75 mm
	IM	265 x 255 x 45 mm	265 x 255 x 65 mm
6	AM	95-01-120-1 400 x 225 x 55 mm	95-01-122-1 400 x 225 x 75 mm
	IM	380 x 190 x 45 mm	380 x 190 x 65 mm
7	AM	95-01-124-1 400 x 290 x 55 mm	95-01-126-1 400 x 290 x 75 mm
	IM	380 x 260 x 45 mm	380 x 260 x 65 mm
8	AM	95-01-128-1 450 x 260 x 55 mm	95-01-130-1 450 x 260 x 75 mm
	IM	430 x 230 x 45 mm	430 x 230 x 65 mm
9	AM	95-01-132-1 400 x 150 x 55 mm	95-01-134-1 400 x 150 x 75 mm
	IM	380 x 125 x 45 mm	380 x 125 x 65 mm
10	AM	95-01-136-1 495 x 270 x 55 mm	95-01-138-1 495 x 270 x 75 mm
	IM	470 x 250 x 45 mm	470 x 250 x 65 mm

Identification Labels



V2A Rust free
 without marking 68 x 22 mm 95-01-200-1
 with Laser marking 68 x 22 mm 95-01-202-1

Aluminium, without marking, 68 x 22 mm



Red 95-01-204-1



Blue 95-01-206-1



Green 95-01-208-1



Black 95-01-210-1



Brown 95-01-212-1




Yellow 95-01-214-1



with single-line marking
 max. 17 signs 95-01-216-1
 aluminium all available colors



with 2-line marking
 max. 17 signs each line 95-01-218-1
 aluminium all available colors



Paper label with indicator 68 x 22 mm 95-01-220-1


Accessories



Plastic Security Seal

95-01-222-1

Silicon Mat




Size1	195 x 135 x 15 mm	95-01-224-1
Size2	390 x 91 x 15 mm	95-01-226-1
Size3	392 x 144 x 15 mm	95-01-228-1
Size4	275 x 185 x 15 mm	95-01-230-1
Size5	255 x 245 x 15 mm	95-01-232-1
Size6	370 x 185 x 15 mm	95-01-234-1
Size7	365 x 252 x 15 mm	95-01-236-1
Size8	420 x 220 x 15 mm	95-01-238-1
Size9	374 x 120 x 15 mm	95-01-240-1
Size10	460 x 231 x 15 mm	95-01-242-1
Universal	600 x 300 x 15 mm	95-01-244-1

Single use paper filter (1pack=100pcs)



Size1	95-01-246-1
Size2	95-01-248-1
Size3	95-01-250-1
Size4	95-01-252-1
Size5	95-01-254-1
Size6	95-01-256-1
Size7	95-01-258-1
Size8	95-01-260-1
Size9	95-01-262-1
Size10	95-01-264-1

Long Life textile filters



Size1	95-01-266-1
Size2	95-01-268-1
Size3	95-01-270-1
Size4	95-01-272-1
Size5	95-01-274-1
Size6	95-01-276-1
Size7	95-01-278-1
Size8	95-01-280-1
Size9	95-01-282-1
Size10	95-01-284-1


Kleinsetcontainer Mini Container

Bitte beachten Sie, nur der Deckel des Containers ist in Farbe.
 Please note that only the lid of the container is coloured.



außen / outside
 300 x 140 x 40 mm
 innen / inside
 280 x 130 x 27 mm

Boden ungelocht
 Bottom non perforated



95-02-100-1
 95-02-104-1
 95-02-108-1
 95-02-112-1
 95-02-116-1
 95-02-120-1



Boden gelocht
 Bottom perforated



95-02-102-1
 95-02-106-1
 95-02-110-1
 95-02-114-1
 95-02-118-1
 95-02-122-1



außen / outside
 300 x 140 x 70 mm
 innen / inside
 280 x 130 x 55 mm



95-02-130-1
 95-02-134-1
 95-02-138-1
 95-02-142-1
 95-02-146-1
 95-02-150-1



95-02-132-1
 95-02-136-1
 95-02-140-1
 95-02-144-1
 95-02-148-1
 95-02-152-1



außen / outside
 300 x 140 x 100 mm
 innen / inside
 280 x 130 x 84 mm



95-02-160-1
 95-02-164-1
 95-02-168-1
 95-02-172-1
 95-02-176-1
 95-02-180-1



95-02-162-1
 95-02-166-1
 95-02-170-1
 95-02-174-1
 95-02-178-1
 95-02-182-1

Zubehör Accessories


Siebkörbe aus Edelstahlblech sind auch mit Füßen erhältlich
Stainless steel screen basket also with feet available



Drahtsiebkorb
für Minicontainer
Wire basket
for Mini container

Artikelnr.
Article-No.
95-02-200-1

Maße
Dimensions
265 x 120 x 45 mm



Siebkorb
aus Edelstahlblech
Stainless steel
screen basket

Artikelnr.
Article-No.
95-02-202-1

Maße
Dimensions
268 x 124 x 22 mm



Siebkorb
aus Edelstahlblech
Stainless steel
screen basket

Artikelnr.
Article-No.
95-02-204-1

Maße
Dimensions
271 x 117 x 45 mm



Spannrahmen
Drape Retainer

Artikelnr.
Article-No.
95-02-206-1

Maße
Dimensions
280 x 125 mm



Einschlagtuch
Drape

Artikelnr.
Article-No.
95-02-208-1

Maße
Dimensions
600 x 400 mm

Container mit Deckel Containers with lid

Bitte beachten Sie, nur der Deckel des Containers ist in Farbe.
 Please note that only the lid of the container is coloured.



außen / outside
 310 x 190 x 40 mm
 innen / inside
 280 x 182 x 35 mm

Boden ungelocht
 Bottom non perforated



95-03-100-1
 95-03-104-1
 95-03-108-1
 95-03-112-1
 95-03-116-1
 95-03-120-1



Boden gelocht
 Bottom perforated



95-03-102-1
 95-03-106-1
 95-03-110-1
 95-03-114-1
 95-03-118-1
 95-03-122-1



außen / outside
 310 x 190 x 65 mm
 innen / inside
 280 x 182 x 55 mm



95-03-130-1
 95-03-134-1
 95-03-138-1
 95-03-142-1
 95-03-146-1
 95-03-150-1



95-03-132-1
 95-03-136-1
 95-03-140-1
 95-03-144-1
 95-03-148-1
 95-03-152-1



außen / outside
 310 x 190 x 80 mm
 innen / inside
 280 x 182 x 70 mm



95-03-160-1
 95-03-164-1
 95-03-168-1
 95-03-172-1
 95-03-176-1
 95-03-180-1



95-03-162-1
 95-03-166-1
 95-03-170-1
 95-03-174-1
 95-03-178-1
 95-03-182-1

Zubehör Accessories

Siebkörbe aus Edelstahlblech sind auch mit Füßen erhältlich
Stainless steel screen basket also with feet available



Drahtsiebkorb
für Minicontainer
Wire basket
for Mini container

Artikelnr.
Article-No.
95-03-200-1


Maße
Dimensions
265 x 170 x 45 mm



Siebkorb
aus Edelstahlblech
Stainless steel
screen basket

Artikelnr.
Article-No.
95-03-202-1


Maße
Dimensions
265 x 175 x 43 mm



Tablett aus
Edelstahlblech
Stainless steel
tray

Artikelnr.
Article-No.
95-03-204-1

Maße
Dimensions
265 x 175 mm



Spannrahmen
Drape Retainer

Artikelnr.
Article-No.
95-03-206-1

Maße
Dimensions
310 x 190 mm



Einschlagtuch
Drape

Artikelnr.
Article-No.
95-03-208-1

Maße
Dimensions
850 x 850 mm


Barriere Filter Einsatz / Barrier Filter Element



Das Bio-Barriere System
 The Bio-Barriere system



Sicherheits Vorrichtung zum Verplomben des Containers
 Security mechanism for the sealing of the container



Schutz vor mechanischen Verdeckter Lufteinlass zum Schäden
 Hidden air-inlet for protection against mechanical damages




Kleinsetcontainer Mini Container

außen / outside 300 x 135 x 47 mm	außen / outside 300 x 135 x 62 mm	außen / outside 300 x 135 x 71 mm
95-04-100-1	95-04-102-1	95-04-104-1




Dental Container

außen / outside 310 x 190 x 50 mm	außen / outside 310 x 190 x 60 mm
95-04-200-1	95-04-202-1



1/2 Dental Container

außen / outside 205 x 145 x 50 mm	außen / outside 205 x 145 x 60 mm
95-04-300-1	95-04-302-1





Bross Steribox is made of a transparent, highly stable and temperature resistant special plastic material tested up to 150 °C. It is therefore especially suitable for autoclave vacuum sterilisation at 121 °C or 134°C as the method most recommended by the relevant authorities.

It is also suitable for plasma sterilisation and for gas sterilisation at a maximum of 65 °C and twice the

ventilation time as that of a metal container. Lid and bottom part close hermetically. A built-in bacteria barrier prevents the penetration of germs and bacteria without additional gasket.

The Bross Steribox is CE-Market and meets the requirements of standards EN 868-8, DIN 58952 and 58953 applicable to sterilisation containers.

***for export only.
Not available for Germany**

Description	Exterior dimensions	Interior dimensions	Article No.
Steribox A8 Long-term filter	90 x 90 x 55 mm	65 x 60 x 47 mm	95-10-100-1
Steribox A6 Long-term filter	190 x 110 x 60 mm	165 x 92 x 49 mm	95-10-102-1
Steribox A4 Long-term filter	290 x 190 x 60 mm	265 x 165 x 49 mm	95-10-104-1

* All Steriboxes are completely assembled. Accessories made of silicone or metal are not included.
We recommend to change the sterile filter after 150 autoclaving processes, at latest after one year.

Accessoires

Description	Pack Unit	Article No.
Incription Seal STEAM (change to brown)	2000	95-10-200-1
Incription Seal PLASMA (change to pink)	2000	95-10-202-1
safety seals	100	95-10-204-1
Date Stamp	1	95-10-206-1
Silicone Cushion Mat A8	1	95-10-208-1
Silicone Cushion Mat A6	1	95-10-210-1
Silicone Cushion Mat A4	1	95-10-212-1
Long-Term-Sterile Filter A8	1	95-10-214-1
Long-Term-Sterile Filter A6	1	95-10-216-1
Long-Term-Sterile Filter A4	1	95-10-218-1

Silicone Grid Trays


Grid Tray A6	1	95-10-220-1
Grid Tray A4	1	95-10-222-1

*Toolsafe Grid Trays made of silicone will be used as underlayers for instruments.
In Combination with Cushion Mats or Metal Grid Inserts, the Grid Tray will additional be used to fix the instruments and can also be used as second storage level.

Metal Grid Trays

Grid Insert A6	1	95-10-224-1
Grid Insert A4	1	95-10-226-1

*Metal Grid Inserts for automatic processing the metal grid inserts offer an optimal solution for all instruments due to the fact that water gains an uniform access to all items being cleaned.



Mini Trays

Mini Tray - Base, Lid, & Mat	3.8 x 6.8 x 1.9	95-20-100-1
Mini Tray - Base (deep), Lid, & Mat	3.8 x 6.8 x 2.5	95-20-110-1
Mini Tray - Base (deep), Lid, & Mat	3.8 x 6.8 x 3.2	95-20-120-1

Large Microsurgical Instrument Trays

Large Tray - Base, Lid, & Mat	25.4 x 38.1 x 1.9	95-20-130-1
Large Tray - Base, Insert Tray, Lid, & 2 mats	25.4 x 38.1 x 3.8	95-20-140-1
Large Tray - Base (deep), Lid, & Mat	25.4 x 38.1 x 3.8	95-20-150-1

Phaco Trays

Phaco Tray - Base, Lid, & 2 Bars / 1 Slot	10.1 x 19 x 3.8	95-20-160-1
Phaco Tray - Base, Lid, & 2 Bars / 3 Slot	10.1 x 19 x 3.8	95-20-170-1
Phaco Tray - Base, Lid, & 2 Bars / 1 Slot	15.2 x 25.4 x 3.8	95-20-180-1
Phaco Tray - Base, Lid, & 2 Bars / 3 Slot	15.2 x 25.4 x 3.8	95-20-190-1

Diamond Knife Trays

Knife Tray, 1 Slot - Base, Lid, & 2 Bars	6.4 x 15.2 x 1.9	95-20-200-1
Knife Tray, 2 Slot - Base, Lid, & 2 Bars	6.4 x 15.2 x 1.9	95-20-210-1
Knife Tray, 3 Slot - Base, Lid, & 2 Bars	6.4 x 15.2 x 1.9	95-20-220-1
Knife Tray, 5 Slot - Base, Lid, & 2 Bars	10.2 x 16.2 x 1.9	95-20-230-1
Knife Tray, 1 Slot - Base, Lid, & 2 Bars	6.4 x 19 x 1.9	95-20-240-1
Knife Tray, 2 Slot - Base, Lid, & 2 Bars	6.4 x 19 x 1.9	95-20-250-1
Knife Tray, 3 Slot - Base, Lid, & 2 Bars	6.4 x 19 x 1.9	95-20-260-1

Microsurgical Instrument Trays

Micro Tray - Base, Lid, & Mat	6.4 x 15.2 x 1.9	95-20-300-1
Micro Tray - Base, Lid, & Mat	6.4 x 15.2 x 3.2	95-20-310-1
Micro Tray - Base, Lid, & Mat	6.4 x 19 x 1.9	95-20-320-1
Micro Tray - Base, Lid, & Mat	10.2 x 16.2 x 1.9	95-20-330-1
Micro Tray - Base, Lid, & Mat	10.1 x 19 x 1.9	95-20-340-1
Micro Tray - Base, Lid, & Mat	10.1 x 19 x 3.8	95-20-350-1
Micro Tray - Base, Lid, & Mat	15.2 x 25.4 x 1.9	95-20-360-1
Micro Tray - Base, Insert Tray, Lid, & 2 mats	15.2 x 25.4 x 3.8	95-20-370-1
Micro Tray - Base (deep), Lid, & Mat	15.2 x 25.4 x 3.8	95-20-380-1


SUGGESTED INSTRUMENT SETS



CONTENT


English	Deutsch	Page
Retinal Detachment	Ablatio	166.
Astigmatic Keratotomy AK	Astigmatische Keratotomie AK	167.
Chalazion	Chalazion	168.
Dacryocystorhinostomy D.C.R. -Lacrimonal-	Dacryocystorhinostomie D.C.R. -Tränenweg -	169.
Enucleation	Enukleation	170.
Stitch Removal	Fadenentfernung	171.
Foreign Body	Fremdkörper	172.
Glaucoma	Glaukom	173.
Cataract E.C.	Katarakt E.C.	174.
Cataract Phacoemulsification	Katarakt Phakoemulsifikation	175.
Keratoplasty	Keratoplastik	176.
LASIK	LASIK	177.
Lid - Conjunctiva Trichiasis - Ptosis	Lid-Bindehaut Trichiasis	178.
Emergency	Notfall	179.
Pterygium	Pterygium	180.
Radiale Keratotomy RK	Radial Keratotomie RK	181.
Strabismus / Muscle	Strabismus	182.
Vitreotomy -Pars Plana-	Vitrektomie -Pars Plana-	183.

Ablatio - Retinal Detachment



- | | |
|---|---|
| 1. 01-01-104-1 Barraquer-Colibri Eye Speculum | 19. 01-05-130-1 Graefe Strabismus Hook |
| 2. 01-14-480-1 Hartmann Mosquito Forceps | 20. 01-05-135-1 Graefe Strabismus Hook |
| 3. 01-14-510-1 Backhaus Towel Clamp 2x | 21. 01-05-145-1 Graefe Strabismus Hook |
| 4. 01-14-440-1 Dieffenbach Serrefine 2x | 22. 01-10-250-1 Colibri Forceps |
| 5. 01-03-120-1 Troutman-Chris Bladeholder | 23. 01-10-114-1 Graefe Dressing Forceps |
| 6. 01-03-104-1 Micro Blades, sterile, 10 per Box | 24. 01-10-120-1 Graefe Tissue Forceps |
| 7. 01-10-390-1 Elschnig Superior Rectus Forceps | 25. 01-10-285-1 Castroviejo Suturing Forceps |
| 8. 01-11-360-1 Westcott Tenotomy Scissors | 26. 01-10-585-1 Jameson Muscle Forceps |
| 9. 01-11-410-1 Micro Iris Scissors | 27. 01-10-580-1 Jameson Muscle Forceps |
| 10. 01-11-530-1 Eye Scissors | 28. 01-10-350-1 Jaffe Tying Forceps |
| 11. 01-11-620-1 Stevens Tenotomy Scissors | 29. 01-10-375-1 Paufigue Tying Forceps |
| 12. 01-11-100-1 Barraquer Iris Scissors | 30. 01-10-210-1 Bonn Iris Forceps |
| 13. 01-04-200-1 Schioetz Tonometer | 31. 01-12-160-1 Barraquer Needle Holder |
| 14. 01-01-575-1 Sautter Orbital Spatula | 32. 01-12-200-1 Barraquer Needle Holder |
| 15. 01-06-360-1 Kimura Platinum Spatula | 33. 01-12-245-1 Castroviejo Standard Needle Holder |
| 16. 01-05-150-1 Gass Retinal Detachment Hook | 34. 01-04-180-1 Ruler 150mm |
| 17. 01-05-110-1 Jameson Muscle Hook | 35. 01-04-104-1 Castroviejo Caliper |
| 18. 01-05-120-1 Green Muscle Hook | |

Astigmatische Keratotomie AK




1. **01-01-132-1** Barraquer Phaco Eye Speculum
2. **01-09-170-1** Sauter Air Injection Cannula 27G curved **2x**
3. **01-08-310-1** Thornton Swivel Fixation Ring FH \varnothing 16mm 3/4
4. **01-08-290-1** Bore's Axis Marker for correction of astigmatism
5. **01-04-120-1** Mendez Degree Gauge C 10° from 0 - 180°
6. **01-08-350-1** Lindstrom Accurate Incision Marker 7.0mm optical zone
7. **01-08-210-1** Hoffer Optical Center Marker with cross hairs 7.0mm
8. **01-08-220-1** Hoffer Optical Center Marker with cross hairs 9.0mm
9. **01-02-690-1** Premium Diamond Knife Micrometer Arcuate Balde 15° 0.2-1.0mm
10. **01-10-230-1** Harms Colibri Forceps 7.5cm with Tying Platform B-teeth 0.12mm

Chalazion




1. **01-01-530-1** Jaeger Lid Plate 11cm
2. **01-11-540-1** Eye Scissors with large rings 10cm curved - sharp tips
3. **01-10-300-1** Paufigue Universal Suturing Forceps 8.5cm with Tying Platform C-Teeth 0.5mm
4. **01-10-630-1** Desmarres Chalazion Forceps 9cm oval - small
5. **01-10-632-1** Desmarres Chalazion Forceps 9cm oval - medium
6. **01-06-490-1** Meyhoefer Chalazion Curette FH \varnothing 1.0mm
7. **01-06-520-1** Meyhoefer Chalazion Curette FH \varnothing 2.0mm
8. **01-06-530-1** Meyhoefer Chalazion Curette FH \varnothing 2.5mm
9. **01-12-255-1** Castroviejo Standard Needle Holder 14cm with lock curved

Dacryocystorhinostomie D.C.R.-Tränenweg / Dacryocystorhinostomy D.C.R.-Lacrimal




- | | |
|--|--|
| 1. 01-09-440-1 Bishop-Harmon Anterior Chamber Irrigator | 16. 01-07-140-1 Bowman Lachrymal Probe #0000/000 |
| 2. 01-01-350-1 Stevenson Lacrimal Sac Retractor | 17. 01-07-145-1 Bowman Lachrymal Probe #00/0 |
| 3. 01-14-510-1 Backhaus Towel Clamp 2x | 18. 01-07-155-1 Bowman Lachrymal Probe #1/2 |
| 4. 01-14-470-1 Hartmann Mosquito Forceps 2x | 19. 01-07-165-1 Bowman Lachrymal Probe #3/4 |
| 5. 01-14-480-1 Hartmann Mosquito Forceps 2x | 20. 01-07-175-1 Bowman Lachrymal Probe #5/6 |
| 6. 01-14-400-1 Kerrison Bone Nibbling Rongeur | 21. 01-10-115-1 Graefe Dressing Forceps |
| 7. 01-14-290-1 Hartmann-Halle Nasal Speculum | 22. 01-10-300-1 Paufique Universal Suturing Forceps |
| 8. 01-14-420-1 Mini-Friedmann Bone Rongeur Forceps | 23. 01-10-735-1 Bonaccolto Conjunctiva Forceps |
| 9. 01-14-220-1 Mallet with interchangeable disks | 24. 01-09-250-1 Anel Lacrimal Cannula 23G straight |
| 10. 01-07-130-1 Castroviejo Lacrimal Dilator | 25. 01-09-255-1 Anel Lacrimal Cannula 23G 90° curved |
| 11. 01-14-250-1 Joseph Periosteal Elevator | 26. 01-10-285-1 Castroviejo Suturing Forceps |
| 12. 01-14-190-1 West Lacrimal Sac Chisel | 27. 01-05-300-1 Knapp Lacrimal Retractor |
| 13. 01-14-200-1 West Lacrimal Sac Gouge | 28. 01-11-640-1 Stevens Tenotomy Scissors |
| 14. 01-14-210-1 Freer Lacrimal Sac Chisel | 29. 01-11-530-1 Eye Scissors, large rings 10cm straight |
| 15. 01-07-240-1 Worst Pigtail Probe | 30. 01-12-255-1 Castroviejo Standard Needle Holder 14cm |
| | 31. 01-11-410-1 Micro Iris Scissors 10cm curved |

Enukleation / Enucleation



- | | | | |
|-----|--------------------|--|-----------|
| 1. | 01-14-510-1 | Backhaus Towel Clamp 9cm | 2x |
| 2. | 01-14-495-1 | Halsted-Mosquito Artery Forceps 14cm straight | 2x |
| 3. | 01-14-505-1 | Halsted-Mosquito Artery Forceps 14cm curved | 2x |
| 4. | 01-01-268-1 | Cook Eye Speculum 15mm solid blades adult | |
| 5. | 01-10-385-1 | Graefe Fixation Forceps 11cm | |
| 6. | 01-10-300-1 | Paufique Universal Suturing Forceps 8.5cm | |
| 7. | 01-11-540-1 | Eye Scissors with large rings 10cm curved | |
| 8. | 01-11-640-1 | Stevens Tenotomy Scissors 11.5cm curved | |
| 9. | 01-06-650-1 | Wells Enucleation Spoon large | |
| 10. | 01-06-660-1 | Bunge Evisceration Spoon RH small | |
| 11. | 01-05-145-1 | Graefe Strabismus Hook FH #3 large | |
| 12. | 01-10-615-1 | Bangerter Muscle Clamp 10cm | |
| 13. | 01-11-670-1 | Enucleation Scissors 13cm curve | |
| 14. | 01-12-255-1 | Castroviejo Standard Needle Holder 14cm with lock curved | |
| 15. | 01-10-770-1 | Carter Sphere Introducer and Holder 12.5cm | |

Fadenentfernung / Stitch Removal



1. **01-01-104-1** Barraquer-Colibri Eye Speculum
2. **01-10-730-1** Barraquer Colibri Suture Removal Forceps
3. **01-11-160-1** Vannas Capsulotomy Scissors

Fremdkörper / Foreign Body



1.




2.



3.


1. **01-01-104-1** Barraquer-Colibri Eye Speculum 15mm adult
2. **01-06-490-2** Meyhoefer Chalazion Curette RH ø 1.0mm
3. **01-03-440-1** Foreign Body Instrument Needle and Spud reversible

Glaukom / Glaucoma




- | | | | |
|-----------------|--|----|--|
| 1. 01-01-180-1 | Corcelle Eye Speculum 8mm | | |
| 2. 01-14-470-1 | Hartmann Mosquito Forceps 10cm | | |
| 3. 01-14-510-1 | Backhaus Towel Clamp 9cm | | |
| 4. 01-10-390-1 | Elschnig Superior Rectus Forceps 10cm | | |
| 5. 01-12-245-1 | Castroviejo Standard Needle Holder | | |
| 6. 01-12-165-1 | Barraquer Needle Holder RH 11cm | | |
| 7. 01-11-510-1 | Mod. Bonn Eye Scissors 9cm curved | | |
| 8. 01-11-360-1 | Westcott Tenotomy Scissors 11cm | | |
| 9. 01-11-190-1 | Mod. Tübingen Capsulotomy Scissors | | |
| 10. 01-11-100-1 | Barraquer Iris Scissors 5.5cm | | |
| 11. 01-04-104-1 | Castroviejo Caliper 0-20mm | | |
| 12. 01-09-440-1 | Bishop-Harmon Anterior Chamber Irrigator | | |
| 13. 01-09-160-1 | McIntyre Anterior Chamber Cannula | 2x | |
| 14. 01-09-170-1 | Sauter Air Injection Cannula | 2x | |
| 15. 01-02-900-1 | Premium Diamond Knife | | |
| 16. 01-02-770-1 | Premium Diamond Knife | | |
| 17. 01-07-320-1 | Elliot Handle with chuck for trephine | | |
| 18. 01-07-325-1 | Elliot Scleral Trephine ø 1.0mm | | |
| 19. 01-07-245-1 | Harms Trabeculotomy Probe right | | |
| 20. 01-07-250-1 | Harms Trabeculotomy Probe left | | |
| 21. 01-14-320-1 | Kelly Descemet's Membrane | | |
| 22. 01-10-121-1 | Graefe Tissue Forceps 10cm straight | | |
| 23. 01-10-120-1 | Graefe Tissue Forceps 10cm straight | | |
| 24. 01-10-245-1 | Barraquer Colibri Forceps 7.5cm | | |
| 25. 01-10-345-1 | Jaffe Tying Forceps 10cm straight | 2x | |
| 26. 01-10-143-1 | Bonn Iris Forceps 7cm straight | | |
| 27. 01-06-280-1 | Mod. Dortmund Phaco-Spatula | | |
| 28. 01-04-200-1 | Schioetz Tonometer with inclined scale | | |

Katarakt E.C. / Cataract E.C.




- | | |
|---|---|
| 1. 01-01-132-1 Barraquer-Phaco Eye Speculum | 23. 01-11-200-1 Mod.Tübingen Capsulotomy Scissors |
| 2. 01-14-470-1 Hartmann Mosquito Forceps | 24. 01-11-160-1 Vannas Capsulotomy Scissors |
| 3. 01-14-510-1 Backhaus Towel Clamp 2x | 25. 01-10-460-1 Masket Capsulorhexis Forceps |
| 4. 01-10-390-1 Elschmig Superior Rectus Forceps | 26. 01-10-245-1 Barraquer Colibri Forceps |
| 5. 01-12-100-1 Kalt Needle Holder | 27. 01-10-345-1 Jaffe Tying Forceps |
| 6. 01-12-245-1 Castroviejo Standard Needle Holder | 28. 01-10-360-1 K-McPherson Tying Forceps |
| 7. 01-12-165-1 Barraquer Needle Holder | 29. 01-06-280-1 Mod.Dortmund Phaco-Spatula |
| 8. 01-11-510-1 Mod.Bonn Eye Scissors | 30. 01-06-420-1 Snellen Lens Loop |
| 9. 01-11-360-1 Westcott Tenotomy Scissors | 31. 01-10-146-1 Bonn Iris Forceps |
| 10. 01-10-120-1 Graefe Tissue Forceps | 32. 01-11-100-1 Barraquer Iris Scissors |
| 11. 01-14-440-1 Dieffenbach Serrefine 2x | 33. 01-04-104-1 Castroviejo Caliper |
| 12. 01-03-120-1 Troutman-Chris Bladeholder | 34. 01-05-255-1 Hirschmann Iris Hook |
| 13. 01-03-104-1 Micro Blades, sterile, 10 per Box | 35. 01-10-525-1 Blaydes Lens Inserting Forceps - Titanium |
| 14. 01-09-160-1 McIntyre Cannula | 36. 01-05-405-1 Sinskey IOL Manipulator - Titanium |
| 15. 01-09-370-1 Jensen Capsule Polisher | 37. 01-04-110-1 Maloney Intra-Operative Keratometer - Titanium |
| 16. 01-09-170-1 Sauter Air Injection 2x | |
| 17. 01-09-410-1 Knolle-Pearce Irrigating Vectis | OPTION: |
| 18. 01-09-295-1 Simcoe Irrigating-Aspirating Cannula | 38. 01-02-180-1 Diamond Knife straight Blade 45° 1.0mm |
| 19. 01-09-440-1 Bishop-Harmon Anterior Chamber Irrigator | 39. 01-02-760-1 Diamond Knife 7 Steps Blade 45° 1.0mm |
| 20. 01-11-260-1 Castroviejo Corneoscleral Scissors | 40. 01-06-160-1 Kansas Nucleus Trisector |
| 21. 01-11-270-1 Castroviejo Corneoscleral Scissors | 41. 01-06-200-1 Kansas Nucleus Vectis |
| 22. 01-11-410-1 Micro Iris Scissors | 42. 01-06-170-1 Kansas Nucleus Trisector |
| | 43. 01-06-190-1 Kansas Nucleus Vectis |
| | 44. 01-10-440-1 Kansas Nucleus Forceps |

Katarakt Phakoemulsifikation / Cataract Phacoemulsification




- | | |
|---|---|
| 1. 01-01-132-1 Barraquer Phaco Eye Speculum | 24. 01-10-345-1 Jaffe Tying Forceps |
| 2. 01-01-302-1 Lieberman Eye Speculum | 25. 01-10-360-1 K-McPherson Tying Forceps |
| 3. 01-14-470-1 Hartmann Mosquito Forceps | 26. 01-05-355-2 Nagahara Phaco Chopper Blade 60° |
| 4. 01-14-510-1 Backhaus Towel Clamp 9cm 2x | 27. 01-05-320-2 Jaffe-Bechert Lens Nucleus Rotator |
| 5. 01-10-390-1 Elschmig Superior Rectus Forceps | 28. 01-05-255-1 Hirschmann Iris Hook |
| 6. 01-12-100-1 Kalt Needle Holder 14cm | 29. 01-04-104-1 Castroviejo Caliper 0-20mm |
| 7. 01-12-245-1 Castroviejo Standard Needle Holder | 30. 01-10-525-2 Blaydes Lens Inserting Forceps |
| 8. 01-12-165-1 Barraquer Needle Holder RH | 31. 01-05-405-2 Sinskey IOL Manipulator angled |
| 9. 01-11-510-1 Mod.Bonn Eye Scissors | 32. 01-04-110-2 Maloney Intra-Operative Keratometer |
| 10. 01-11-360-1 Westcott Tenotomy Scissors | 33. 01-10-490-1 Faulkner Lens Holding Forceps |
| 11. 01-10-120-1 Graefe Tissue Forceps | 34. 01-10-485-2 Nordan Lens Folding Forceps |
| 12. 01-14-440-1 Dieffenbach Serrefine 2x | 35. 01-10-495-1 McDonald Lens Inserting Forceps |
| 13. 01-03-330-1 Mod.Hamburg Keratome | |
| 14. 01-03-120-1 Troutman-Chris Bladeholder | Option: |
| 15. 01-03-104-1 Micro Blades, sterile, 10 per Box | 36. 01-08-310-1 Thornton Swivel Fixation Ring |
| 16. 01-09-160-1 McIntyre Anterior Chamber Cannula | 37. 01-08-300-1 Thornton Swivel Fixation Ring |
| 17. 01-09-370-1 Jensen Capsule Polisher | 38. 01-09-635-2 Heaflinger Irrigation Handpiece |
| 18. 01-09-170-1 Sauter Air Injection Cannula 2x | 39. 01-09-620-2 Heaflinger Aspiration Handpiece |
| 19. 01-09-440-1 Bishop-Harmon Anterior Chamber Irrigator | 40. 01-13-100-1 SI Capsulorhexis Forceps IO |
| 20. 01-11-200-1 Mod.Tübingen Capsulotomy Scissors | 41. 01-02-320-1 Diamond Knife for Scleral Tunnel |
| 21. 01-11-160-1 Vannas Capsulotomy Scissors | 42. 01-02-640-1 Diamond Knife for LRI |
| 22. 01-10-460-1 Masket Capsulorhexis Forceps | 43. 01-02-770-1 Diamond Knife 7 Step for LRI |
| 23. 01-10-245-1 Barraquer Colibri Forceps | 44. 01-02-530-1 Diamond Knife for Clear Cornea Incisions |
| | 45. 01-02-300-1 Diamond Knife for Wound Widening |

Keratoplastik / Keratoplasty




- | | |
|---|--|
| 1. 01-01-104-1 Barraquer-Colibri Eye Speculum | 18. 01-11-300-1 Katzin Corneal Transplant Scissors |
| 2. 01-14-470-1 Hartmann Mosquito Forceps | 19. 01-11-280-1 Corneal Scissors 10cm |
| 3. 01-14-510-1 Backhaus Towel Clamp 9cm 2x | 20. 01-11-160-1 Vannas Capsulotomy Scissors 8cm |
| 4. 01-10-120-1 Graefe Tissue Forceps 10cm | 21. 01-07-380-1 Franceschetti Keratoplasty Trephine ø7.0mm |
| 5. 01-10-390-1 Elschmig Superior Rectus Forceps | 22. 01-07-390-1 Franceschetti Keratoplasty Trephine ø 8.0mm |
| 6. 01-12-100-1 Kalt Needle Holder 14cm | 23. 01-06-690-1 Barraquer Keratoplastic Spatula |
| 7. 01-12-245-1 Castroviejo Standard Needle Holder | 24. 01-08-260-1 Anis Suture Placement Marker |
| 8. 01-12-165-1 Barraquer Needle Holder | 25. 01-10-245-1 Barraquer Colibri Forceps 7.5cm |
| 9. 01-07-475-1 Flieringa Scleral Fixation Ring ø 14mm | 26. 01-10-290-1 Castroviejo Suturing Forceps 10cm |
| 10. 01-07-485-1 Flieringa Scleral Fixation Ring ø 16mm | 27. 01-10-345-1 Jaffe Tying Forceps 10cm 2x |
| 11. 01-07-495-1 Flieringa Scleral Fixation Ring ø 18mm | 28. 01-09-170-1 Sauter Air Injection Cannula 2x |
| 12. 01-07-505-1 Flieringa Scleral Fixation Ring ø 20mm | 29. 01-04-110-2 Maloney Intra-Operative Keratometer |
| 13. 01-03-240-1 Bonn Scleral Knife | |
| 14. 01-09-440-1 Bishop-Harmon Anterior Chamber Irrigator | Option: |
| 15. 01-11-510-1 Mod.Bonn Eye Scissors | 30. 01-02-180-1 Diamond Knife straight Blade 45° |
| 16. 01-11-360-1 Westcott Tenotomy Scissors | |
| 17. 01-11-290-1 Katzin Corneal Transplant Scissors | |

Lasik / Lasik




1. **01-01-310-1** Lieberman Aspirating Eye Speculum 15mm
2. **01-01-302-1** Lieberman Eye Speculum K-Wire 15mm
3. **01-08-300-1** Thornton Swivel Fixation Ring FH \varnothing 13mm 3/4
4. **01-08-370-1** Chayet Corneal Marker for Lasik \varnothing 8.5mm
5. **01-06-370-1** Fukasaku Lasik Spatula angled fine
6. **01-06-390-1** Flap Elevator for Lasik double ended flat and cylindrical
7. **01-09-170-1** Sauter Air Injection Cannula 27G curved
8. **01-09-610-1** Brauweiler Irrigating Cannula 21G angled
9. **01-10-335-1** Buratto Lasik Flap Forceps 10cm curved disc shaped serrated jaws
10. **01-01-605-1** Lasik Flap Smoother

Lid-Bindehaut Trichiasie-Ptosis / Lid-Conjunctiva Trichiasis-Ptosis




- | | |
|---|---|
| 1. 01-09-440-1 Bishop-Harmon Anterior Chamber Irrigator | 19. 01-05-135-1 Graefe Strabismus Hook medium |
| 2. 01-01-530-1 Jaeger Lid Plate 11cm | 20. 01-05-260-1 Guthrie Double-Fixation Hook small |
| 3. 01-14-540-1 Jones Towel Clamp 8.5cm | 2x 21. 01-05-300-1 Knapp Lacrimal Retractor RH |
| 4. 01-14-470-1 Hartmann Mosquito Forceps 10cm straight | 2x 22. 01-06-490-1 Meyhoefer Chalazion Curette FH \varnothing 1.0mm |
| 5. 01-14-480-1 Hartmann Mosquito Forceps 10cm curved | 23. 01-06-520-1 Meyhoefer Chalazion Curette FH \varnothing 2.0mm |
| 6. 01-11-620-1 Stevens Tenotomy Scissors 10.5cm curved | 24. 01-06-530-1 Meyhoefer Chalazion Curette FH \varnothing 2.5mm |
| 7. 01-11-550-1 Eye Scissors 11cm for lids straight | 25. 01-06-550-1 Meyhoefer Chalazion Curette FH \varnothing 3.5mm |
| 8. 01-03-240-1 Bonn Scleral Knife Blade 5mm | 26. 01-10-385-1 Graefe Fixation Forceps 11cm |
| 9. 01-11-360-1 Westcott Tenotomy Scissors 11cm curved | 27. 01-10-735-1 Bonaccolto Conjunctiva Forceps 10cm |
| 10. 01-11-410-1 Micro Iris Scissors 10cm curved | 28. 01-10-108-1 Storz Cilia Forceps 8.5cm |
| 11. 01-12-255-1 Castroviejo Standard Needle Holder 14cm | 29. 01-10-285-1 Castroviejo Suturing Forceps 10cm |
| 12. 01-10-620-1 Berke Ptosis Forceps 11cm with lock | 30. 01-10-290-1 Castroviejo Suturing Forceps 10cm |
| 13. 01-12-100-1 Kalt Needle Holder 14cm with lock straight | 31. 01-10-342-1 Harms Tying Forceps 10cm |
| 14. 01-04-104-1 Castroviejo Caliper 0-20mm | 32. 01-10-632-1 Desmarres Chalazion Forceps 9cm |
| 15. 01-01-504-1 Desmarres Lid Retractor FH 12mm | 33. 01-10-640-1 Lambert Chalazion Forceps 9cm |
| 16. 01-01-506-1 Desmarres Lid Retractor FH 14mm | |
| 17. 01-01-508-1 Desmarres Lid Retractor FH 16mm | |
| 18. 01-05-130-1 Graefe Strabismus Hook #1 small | |

Notfall / Emergency




- | | |
|--|---|
| 1. 01-01-318-1 Lieberman Eye Speculum | 22. 01-05-110-1 Jameson Muscle Hook |
| 2. 01-01-100-1 Barraquer-Colibri Eye Speculum, birth | 23. 01-05-135-1 Graefe Strabismus Hook |
| 3. 01-01-102-1 Barraquer-Colibri Eye Speculum, child | 24. 01-05-160-1 Stevens Tenotomy Hook |
| 4. 01-01-132-1 Barraquer Phaco Eye Speculum | 25. 01-01-502-1 Desmarres Lid Retractor FH 10mm |
| 5. 01-14-540-1 Jones Towel Clamp 8.5cm 2x | 26. 01-01-504-1 Desmarres Lid Retractor FH 12mm |
| 6. 01-14-470-1 Hartmann Mosquito Forceps 10cm 2x | 27. 01-01-506-1 Desmarres Lid Retractor FH 14mm |
| 7. 01-14-480-1 Hartmann Mosquito Forceps 10cm 2x | 28. 01-01-600-1 Schocket Scleral Depressor |
| 8. 01-14-440-1 Dieffenbach Serrefine 3.5cm 2x | 29. 01-10-114-1 Graefe Dressing Forceps 10cm |
| 9. 01-09-440-1 Bishop-Harmon Anterior Chamber Irrigator | 30. 01-10-285-1 Castroviejo Suturing Forceps 10cm |
| 10. 01-09-170-1 Sauter Air Injection Cannula | 31. 01-10-290-1 Castroviejo Suturing Forceps 10cm |
| 11. 01-03-120-1 Troutman-Chris Bladeholder | 32. 01-10-210-1 Bonn Iris Forceps 10cm |
| 12. 01-03-170-1 Razor Blades breakable | 33. 01-10-360-1 K-McPherson Tying Forceps 10cm |
| 13. 01-05-280-1 Sclera Twist Fixation Claw | 34. 01-10-375-1 Paufigue Tying Forceps 8.5cm |
| 14. 01-11-360-1 Westcott Tenotomy Scissors 11cm | 35. 01-10-132-1 Bishop-Harmon Tissue Forceps 8.5cm |
| 15. 01-11-410-1 Micro Iris Scissors 10cm | 36. 01-10-395-1 Lester Fixation Forceps 9cm |
| 16. 01-11-510-1 Mod.Bonn Eye Scissors 9cm | 37. 01-10-250-1 Colibri Forceps 7.5cm |
| 17. 01-11-530-1 Eye Scissors with large rings | 38. 01-12-165-1 Barraquer Needle Holder |
| 18. 01-11-640-1 Stevens Tenotomy Scissors 11.5cm | 39. 01-12-210-1 Barraquer Needle Holder |
| 19. 01-11-100-1 Barraquer Iris Scissors 5.5cm | 40. 01-12-245-1 Castroviejo Standard Needle Holder |
| 20. 01-04-200-1 Schioetz Tonometer | 41. 01-04-180-1 Ruler 150mm inches |
| 21. 01-05-150-1 Gass Retinal Detachment Hook | 42. 01-04-104-1 Castroviejo Caliper 0-20mm |

Pterygium / Pterygium




1. **01-01-185-1** Lancaster Eye Speculum 8cm solid blades
2. **01-01-502-1** Desmarres Lid Retractor FH 10mm
3. **01-11-530-1** Eye Scissors with large rings 10cm straight - sharp tips
4. **01-03-240-1** Bonn Scleral Knife Blade 5mm
5. **01-12-250-1** Castroviejo Standard Needle Holder 14cm curved
6. **01-10-300-1** Paufigue Universal Suturing Forceps 8.5cm with Tying Platform C-Teeth 0.5mm
7. **01-10-114-1** Graefe Dressing Forceps 10cm straight serrated - standard
8. **01-10-120-1** Graefe Tissue Forceps 10cm straight 1x2 teeth - standard

Radiale Keratotomie RK / Radial Keratotomy RK




- 1. **01-01-132-1** Barraquer Phaco Eye Speculum
- 2. **01-09-170-1** Sauter Air Injection Cannula 27G curved **2x**
- 3. **01-08-100-1** Hoffer Optical Center Marker with cross hairs 3.0mm
- 4. **01-08-110-1** Hoffer Optical Center Marker with cross hairs 3.25mm
- 5. **01-08-140-1** Hoffer Optical Center Marker with cross hairs 4.0mm
- 6. **01-08-180-1** Hoffer Optical Center Marker with cross hairs 5.0mm
- 7. **01-08-210-1** Hoffer Optical Center Marker with cross hairs 7.0mm
- 8. **01-08-240-1** Anis Radial Marker 8 Blades ID 4mm OD 12.5mm
- 9. **01-08-310-1** Thornton Swivel Fixation Ring FH ø16mm 3/4
- 10. **01-10-230-1** Harms Colibri Forceps 7.5mm with Tying Platform B-teeth 0.12mm
- 11. **01-02-710-1** Premium Diamond Knife Micrometer Balde 45° 1.0mm

Strabismus / Strabismus - Muscle




- | | |
|---|---|
| 1. 01-01-210-1 Castroviejo Eye Speculum | 21. 01-10-210-1 Bonn Iris Forceps 10cm |
| 2. 01-14-470-1 Hartmann Mosquito Forceps | 22. 01-10-290-1 Castroviejo Suturing Forceps |
| 3. 01-14-510-1 Backhaus Towel Clamp | 23. 01-10-285-1 Castroviejo Suturing Forceps |
| 4. 01-14-440-1 Dieffenbach Serrefine | 24. 01-10-375-1 Paufigue Tying Forceps |
| 5. 01-09-440-1 Bishop-Harmon Anterior Chamber Irrigator | 25. 01-10-615-1 Bangerter Muscle Clamp |
| 6. 01-11-530-1 Eye Scissors | 26. 01-12-255-1 Castroviejo Standard Needle Holder |
| 7. 01-11-620-1 Stevens Tenotomy Scissors | 27. 01-12-245-1 Castroviejo Standard Needle Holder |
| 8. 01-11-410-1 Micro Iris Scissors | |
| 9. 01-11-360-1 Westcott Tenotomy Scissors | |
| 10. 01-01-550-1 Helveston Tissue Retractor -Barbie- 7mm | |
| 11. 01-01-552-1 Helveston Tissue Retractor -Big Barbie- 9mm | |
| 12. 01-01-554-1 Helveston Tissue Retractor -Great Big Barbie- 11mm | |
| 13. 01-05-165-1 Helveston -Teaser- Hook extra delicate | |
| 14. 01-05-170-1 Helveston Finder Hook large extra delicate | |
| 15. 01-05-175-1 Helveston Finder Hook small extra delicate | |
| 16. 01-05-180-1 Helveston Muscle Hook 8mm extra delicate | |
| 17. 01-05-185-1 Helveston Muscle Hook 10mm extra delicate | |
| 18. 01-05-190-1 Helveston Muscle Hook 12mm extra delicate | |
| 19. 01-04-170-2 Gonnin Amsler Marker 11-23mm / 1mm | |
| 20. 01-04-104-1 Castroviejo Caliper 0-20mm | |

Vitrektomie-Pars Plana / Vitrectomy-Pars Plana



- | | | | |
|---|-----------|---|-----------|
| 1. 01-01-310-1 Lieberman Aspirating Eye Speculum | | 19. 01-10-735-1 Bonaccolto Conjunctiva Forceps 10cm | |
| 2. 01-14-540-1 Jones Towel Clamp 8.5cm | 2x | 20. 01-10-290-1 Castroviejo Suturing Forceps 10cm | 2x |
| 3. 01-14-440-1 Dieffenbach Serrefine 3.5cm straight | 2x | 21. 01-10-345-1 Jaffe Tying Forceps 10cm | 2x |
| 4. 01-14-510-1 Backhaus Towel Clamp 9cm | 2x | 22. 01-10-245-1 Barraquer Colibri Forceps 7.5cm | |
| 5. 01-14-495-1 Halsted-Mosquito Artery Forceps | 2x | 23. 01-12-245-1 Castroviejo Standard Needle Holder | |
| 6. 01-14-470-1 Hartmann Mosquito Forceps | 2x | 24. 01-12-250-1 Castroviejo Standard Needle Holder, curved | |
| 7. 01-11-530-1 Eye Scissors with large rings | | 25. 01-13-145-1 Vitrectomy Forceps straight serrated | |
| 8. 01-11-540-1 Eye Scissors with large rings | | 26. 01-13-125-1 Vitrectomy Scissors straight | |
| 9. 01-11-630-1 Stevens Tenotomy Scissors | | 27. 01-09-440-1 Bishop-Harmon Anterior Chamber Irrigator | |
| 10. 01-11-360-1 Westcott Tenotomy Scissors 11cm curved | | 28. 01-09-350-1 Injection Cannula 21G | |
| 11. 01-11-150-1 Vannas Capsulotomy Scissors 8cm | | 29. 01-10-765-1 Insertion Forceps 10cm | |
| 12. 01-01-575-1 Sautter Orbital Spatula medium | | | |
| 13. 01-01-550-1 Helveston Tissue Retractor Barbie 7mm | | | |
| 14. 01-03-520-1 LS Stiletto 1.2mm | | | |
| 15. 01-05-110-1 Jameson Muscle Hook | 2x | | |
| 16. 01-05-150-1 Gass Retinal Detachment Hook | 2x | | |
| 17. 01-04-104-1 Castroviejo Caliper 0-20mm | | | |
| 18. 01-10-114-1 Graefe Dressing Forceps | | | |
| | | Option: | |
| | | 30. 01-02-910-1 Premium Diamond Knife | |



Bross Medical GmbH
-Headquarter-
Albstrasse 12
78532 Tuttlingen
GERMANY

-Administration-
Dreifaltigkeitsbergstrasse 12
78588 Denkingen
GERMANY

Phone: +49 (0) 7424 7037330
Fax: +49 (0) 7424 7037339
Mail: Ben@bross-medical.com
simone@bross-medical.com