

a1 medical



O P H T H A L M O L O G Y

COMPANY

a1 medical is a manufacturer of highest quality micro-surgical instruments used in Ophthalmology, offering customers the freedom of choice. These hallmarks have proven to address a gap in the market, since incorporation of the company in 2012. Consequently more and more surgeons and nurses worldwide experience and value the consistent high quality level of a1 medical products and services.

The a1 medical team combines decades of experience in medical technology, with passion for instruments, skills, team spirit and mastery of latest manufacturing technologies. We deliver on time, are responsive and committed to medical progress. Driving improved approaches for medical routines by realizing novel instrument designs together with our partners and customers is therefore our continued effort.

We are independent, rooted in the Tuttlingen region, manufacture and distribute highest quality level instruments Made in Germany, enabling medical professionals to help patients as best as humanly possible.

a1 medical provides you with the reliability patients deserve.

NOTE

This catalogue has been made for you with great care. However, we would like to point out that we assume no liability for the correctness of every detail in the catalogue.

The manufacturer reserves the right to change technical details. No part of this catalogue may be reproduced (printed, photocopied or any other means) or used without the written permission of a1 medical GmbH.

Photos may vary slightly from actual size.

We would like to thank you in advance for your understanding.

TABLE OF CONTENTS

| | Numeric Prefix | CHAPTER | Page |
|-----|-------------------|--|---------|
| 1. | | Ceramic Coating Information | 7-9 |
| 2. | ES- | Chisels, Eye Speculums, Fixation Hooks, Fixation Rings, Lid Plates, Muscle Hooks, Needles, Retractors, Scleral Depressors, Strabismus Hooks | 10-33 |
| 3. | C- | Blade Handles, Burrs, Dissectors, Diamond Knives, Foreign Body Instruments, IOL Cutter, Knives, Punches, Spuds, Trephines, Needles, Corneal Scarifiers | 34-49 |
| 4. | HS- | Fixation Rings, Hooks, Keratometer, Loops, Marker, Manipulators, Marking Pen, Phaco Choppers, Preassure Plate, Spatulas & Curettes, Spoon, Titanium Hooks, a1 medical Single- and Double Hooks | 50-87 |
| 5. | F- | Forceps, MIC Forceps, Titanium Forceps, Bipolar Forceps, ErgoGrip Forceps, TC Forceps | 88-151 |
| 6. | N- | Blade Breakers, Needle Holders & Titanium Needle Holders, ErgoGrip Needle Holders, TC Needle Holders | 152-167 |
| 7. | S- | Micro Scissors, Super Cut Scissors, Ring Scissors | 168-187 |
| 8. | VR- | Vitrectomy Instruments | 188-193 |
| 9. | L- | Bimanual I/A Instruments, Cannulas, Lacrimal Instruments, Detachable I/A Instruments | 194-233 |
| 10. | M- | Marking Instruments | 234-237 |
| 11. | MIS- | Boles, Eye Shields, Miscellaneous, Trays & Cassettes | 238-252 |
| 12. | | Content | 253-260 |
| 13. | | Alphabetical Index | 261-266 |

MADE IN GERMANY

a1 medical

CERAMIC COATING

- ... LONGER PRODUCT LIFE
- ... EASIER CLEANING
- ... REDUCED EYE-FATIGUE
- ... 100 % BIOCOMPATIBLE
- ... ANTI ALLERGIC



CERAMIC COATING

ORDER INFORMATION

Please add „C“ to the end of the article number. e.g. F-0010C

... improves instruments:

- cost reduction and longer product life through
 - higher surface hardness** – up to 10 times harder than stainless steel
 - lower friction** – improves wear and tear of scissor blades
 - more resistance** against corrosion
- reduced eye-fatigue through less glare due to light absorbing black surface
- easier cleaning and reprocessing due to lower surface energy and non-stick effect
- anti-allergic due to absence of Chromium and Nickel in surface layer

... established in cardiology and neurology

... consists of TiAlN, applied in a magnetron sputtering process

... is 100% quality controlled and comes with lifetime warranty

EYE SPECULUMS

RETRACTORS

FIXATION RINGS

FIXATION HOOKS

MUSCLE HOOKS


CHISELS

LID PLATES

SCLERAL DEPRESSORS

FASCIA NEEDLE

STRABISMUS HOOKS



Barraquer Wire Speculum (standard),
10.0mm blades

ES-0010

Barraquer Wire Speculum (standard),
15.0mm blades

ES-0020

Barraquer Wire Speculum (heavy wire),
15.0mm blades

ES-0030



Barrett Eye Speculum adjustable,
14.5mm open blades,
for XEN implantation, used with ES-0725 Vera Hook

ES-0311

For Ceramic Coating please add "C"




Temporal Approach, Wire Speculum,
with spring on nasal side,
15.0mm open blades

ES-0050

Temporal Approach, Wire Speculum,
with spring on nasal side,
15.0mm closed blades

ES-0060




Barraquer Wire Speculum,
10.0mm solid blades

ES-0070

Barraquer Wire Speculum,
15.0mm solid blades

ES-0080

For Ceramic Coating please add "C"




Kratz-Barraquer Wire Speculum,
15.0mm, standard wire

ES-0090

Kratz-Barraquer Wire Speculum,
11.0mm, solid blades

ES-0100



Barraquer Wire Eye Speculum,
14.0mm open blades,
length 4.0cm

ES-0122



Lieberman Aspirating Speculum complete,
15.0mm blades

ES-0130

Replacement silicone tubing

ES-0140



Lieberman Eye Speculum, with Aspiration and Irrigation ports,
15.0mm blades

ES-0150

Replacement silicone tubing

ES-0160

For Ceramic Coating please add "C"



Sauer Lid Speculum, for Infants,
11.0mm blades

ES-0170

Sauer Lid Speculum, for Infants,
6.0mm blades


ES-0180

Sauer Lid Speculum, for premature Infants,
4.0mm blades

ES-0190

Sauer Lid Speculum, for premature Infants,
2.0mm blades


ES-0200



Infant Lid Speculum,
4.0mm blades


ES-0210

For Ceramic Coating please add "C"



Barraquer Infant Wire Speculum,
9.0mm blades

ES-0220



Cook Eye Speculum #1,
8.0mm blades

ES-0230


Cook Eye Speculum #2,
11.0mm blades

ES-0240

Cook Eye Speculum #3,
15.0mm blades

ES-0250

For Ceramic Coating please add "C"




Murdoch Eye Speculum #2,
12.0mm blades

ES-0260

Murdoch Eye Speculum #3,
15.0mm blades

ES-0270




Mellinger Eye Speculum,
open wire, 16.0mm blades

ES-0280

Mellinger Eye Speculum, for children
open wire, 11.0mm blades

ES-0290

For Ceramic Coating please add "C"



Williams Eye Speculum with lock screw,
11.0mm blades

ES-0300

Williams Eye Speculum with lock screw,
15.0mm blades

ES-0310



Lancaster Eye Speculum with locking nut,
16.0mm solid blades

ES-0320



Lieberman Eye Speculum, adjustable mechanism,
10.0mm open wire blades

ES-0330

Lieberman Eye Speculum, adjustable mechanism,
15.0mm open wire blades

ES-0340

Lieberman Eye Speculum, adjustable mechanism,
7.0mm open wire blades

ES-0341

For Ceramic Coating please add "C"



Lieberman Eye Speculum, adjustable mechanism,
15.0mm V-shaped wire blades

ES-0350



Lieberman Eye Speculum, adjustable mechanism,
10.0mm solid blades

ES-0360

Lieberman Eye Speculum, adjustable mechanism,
15.0mm solid blades

ES-0370

Kraff Eye Seculum temporal, adjustable,
14.0mm solid blades

ES-0371



Ginsberg Eye Speculum,
15.0mm blades

ES-0380

For Ceramic Coating please add "C"



Lieberman Reversible Eye Speculum, adjustable mechanism,
for nasal or temporal placement,
15.0mm solid blades

ES-0390



Lieberman Eye Speculum, Brown Triple,
12.0mm blades

ES-0400



Lieberman Eye Speculum,
15.0mm round blades

ES-0410

Lieberman Eye Speculum,
8.0mm round blades,
length 6.0cm

ES-0420

For Ceramic Coating please add "C"



Castroviejo Eye Speculum, adjustable mechanism,
14.0mm blades

ES-0430

Castroviejo Eye Speculum, adjustable mechanism,
16.0mm blades

ES-0440



Fishkind Eye Speculum, adjustable mechanism,
15.0mm blade


ES-0450



Stevenson Lacrimal Sac Retractor,
with 3x3 blunt prongs,
15.0mm blades

ES-0460

For Ceramic Coating please add "C"



Flieringa Scleral Fixation Ring, set of 9 sizes
Ø 14.0mm - 22.0mm (1.0mm increments)

ES-0500

To receive individual rings please add the requested 'mm' at the end of the article number, e.g. ES-050020



McNeil-Goldman Scleral Fixation Ring & Blepharostat,
small 15.0mm

ES-0510

McNeil-Goldman Scleral Fixation Ring & Blepharostat,
medium 16.0mm

ES-0520

McNeil-Goldman Scleral Fixation Ring & Blepharostat,
large 17.0mm

ES-0530

For Ceramic Coating please add "C"



Schepens Scleral Depressor 7.0mm wide tip,
small thimble, #8

ES-0540

Schepens Scleral Depressor 7.0mm wide tip,
medium thimble, #10

ES-0550

Schepens Scleral Depressor 7.0mm wide tip,
large thimble, #12

ES-0560



Gass Scleral Marker,
Ø 1.5mm marking tip

ES-0570



Schocket Scleral Depressor,
with 6.0mm wide bar & 2.5mm teardrop shaped ends,
length 13.5cm

ES-0580

For Ceramic Coating please add "C"



Jaffe Wire Lid Retractor,
(pair) 15.0mm blades

ES-0590

Jaffe Wire Lid Retractor,
(pair) 13.0mm blades

ES-0600



Desmarres Lid Retractor, Fig. 0,
11.0mm

ES-0610

Desmarres Lid Retractor, Fig.1,
14.0mm

ES-0620

Desmarres Lid Retractor, Fig. 2,
15.0mm

ES-0630

For Ceramic Coating please add "C"



Jaeger Lid Plate,
length 11.0cm

ES-0640



Wright Fascia Needle,
1.0 mm x 6.0 mm oval hole,
length 13.0cm

ES-0650

For Ceramic Coating please add "C"



Needle & Gouge,
foreign body spud

ES-0660



Helvoston Tissue Retractor,
7.0mm

ES-0670

Helvoston Tissue Retractor,
9.0mm

ES-0680

Helvoston Tissue Retractor,
11.0mm


ES-0690

For Ceramic Coating please add "C"



Agricola Lacrimal Sac Retractor,
3x3 semi-sharp prongs,
length 4.5cm

ES-0700



Goldstein Lacrimal Sack Retractor

ES-0715

For Ceramic Coating please add "C"



Guthrie Fixation Hook, double, sharp,
small, length 13.0cm

ES-0710

Guthrie Fixation Hook, double, sharp,
large, length 13.0cm

ES-0720



Vera Fixation Hook for XEN Gel Stent implantation,
used with Barret speculum ES-0311

ES-0725

For Ceramic Coating please add "C"



Rollet Lacrimal Sac Retractor, 4 prongs,
sharp, length 13.5cm

ES-0730

Rollet Lacrimal Sac Retractor, 4 prongs
blunt, length 13.5cm

ES-0740



Axenfeld Hook, 2 prongs,
sharp, length 14.0cm

ES-0750

Axenfeld Hook, 2 prongs,
blunt, length 14.0cm

ES-0760

Axenfeld Hook, 3 prongs,
sharp, length 14.0cm

ES-0770

Axenfeld Hook, 3 prongs,
blunt, length 14.0cm

ES-0780

For Ceramic Coating please add "C"



Graefe Muscle Hook, small,
Size 1

ES-0790



Graefe Muscle Hook, medium,
Size 2

ES-0800



Graefe Muscle Hook, large,
Size 3

ES-0810



Helveston Muscle Hook extra delicate, small,
with blunt tip, 8.0mm

ES-0820



Helveston Muscle Hook, extra delicate, medium,
with blunt tip, 10.0mm

ES-0830



Helveston Muscle Hook, extra delicate, large,
with blunt tip, 12.0mm

ES-0840



Jameson Muscle Hook,
small

ES-0850



Jameson Muscle Hook,
medium

ES-0860



Jameson Muscle Hook,
large

ES-0870



Green Muscle Hook

ES-0880

Stevens Tenotomy Hook,
blunt small

ES-0890



Gass Retinal Detachment Hook,
with oval opening

ES-0900

Helveston „Teaser“ Hook, extra delicate,
blunt

ES-0910

For Ceramic Coating please add "C"



Jaeger Strabismus Hook,
Fig. 1

ES-0920

Jaeger Strabismus Hook,
Fig. 2

ES-0930



Manson Double Ended Strabismus Hook

ES-0940



Holder for Absorbent Strips,
serrated jaws

ES-0950

For Ceramic Coating please add "C"



Urretts-Zavalía Depressor, for fixation & pick

ES-0960



Iris Retractor 10.0mm blue micro hooks,
6 retractors in autoclavable case, reusable

ES-1050

Iris Retractor 10.0mm blue micro hooks,
1 set of 5 retractors, disposable, sterile

ES-1060

For Ceramic Coating please add "C"

TREPHINES
SPUDS & GOUGES
NEEDLES
CORNEAL SCARIFIERS
DISSECTORS
BURRS
PUNCHES
IOL CUTTERS
KNIVES
BLADE HANDLES
DIAMOND KNIVES
FOREIGN BODY
INSTRUMENTS



Tooke Corneal Knife,
straight blade & curved cutting edge
3.0x18.0mm, length 11.5cm

C-0010



Gill Corneal Knife,
straight & curved cutting edge

C-0020

Gill Corneal Knife,
slightly curved, curved cutting edge

C-0030



Fukasaku „Hokey“ Knife,
length 11.5cm

C-0040

For Ceramic Coating please add "C"



Bonn Scleral Knife,
4.0mm blade
length 12.5cm

C-0050



Paton Corneal Dissector,
length 11.5cm

C-0060



Kuhnt Corneal Scarifier,
Fig. 1, length 11.5cm

C-0070



Kuhnt Corneal Scarifier,
Fig. 2, length 11.5cm

C-0080



Kuhnt Corneal Scarifier,
Fig. 3, length 11.5cm

C-0090



Kuhnt Corneal Scarifier,
Fig. 4, length 11.5cm

C-0100

Kuhnt Corneal Scarifier,
Fig. 5, length 11.5cm

C-0110

For Ceramic Coating please add "C"



Graefe Cataract Knife,
Fig. 1, 1.5 x 25.0mm,
length 12.5cm

C-0120

Graefe Cataract Knife,
Fig. 2, 2.0 x 30.0mm,
length 12.5cm

C-0130

Graefe Cataract Knife,
Fig. 3, 2.3 x 33.0mm,
length 12.5cm

C-0140



Wheeler Discission Knife,
1.0 x 18.0mm,
length 11.5cm

C-0150

For Ceramic Coating please add "C"



Ziegler Iris Knife,
Fig. 1, 5.0mm,
length 12.0cm

C-0160

Ziegler Iris Knife,
Fig. 2, 6.0mm,
length 12.0cm

C-0170

Ziegler Iris Knife,
Fig. 3, 7.0mm,
length 12.0cm

C-0180

Ziegler Capsulotomy Knife,
length 12.0cm

C-0190



Sato Corneal Knife, straight,
length 12.5cm

C-0200



Foreign Body Needle, straight,
length 12.5cm

C-0210

Foreign Body Needle, curved,
length 12.5cm

C-0220

For Ceramic Coating please add "C"



Bowman Discission Needle, straight,
length 12.5cm

C-0230



Plange Foreign Body Needle,
length 12.0cm

C-0240



Francis Foreign Body Spud,
length 12.0cm

C-0250



La Forces Foreign Body Spud,
length 12.5cm

C-0260

For Ceramic Coating please add "C"



Algerbrush II Corneal Rust Ring Remover, complete C-0280

Burr 0.5mm C-0290

Burr 1.0mm C-0300

Pterygium Drill, round, 3.5mm burr, battery powered, complete C-0310

Replacement burr with chuck C-0320



Diamond Step Knife,
single, 45°, 1.0mm cutting edge

C-0330



Diamond Step Knife,
tri-facet, 1.0mm cutting edge

C-0340



Diamond Step Knife,
10-facet, 1.0mm cutting edge

C-0350



Crescent Diamond Knife, angled,
crescent shaped blade with 2.5mm cutting edge

C-0360



MICS Diamond Knife,
for 1.8mm incision

C-0361

MICS Diamond Knife,
for 2.2mm incision

C-0362



Diamond Phaco Knife, straight, keratome shaped blade, 2.5mm C-0370

Diamond Phaco Knife, angled, keratome shaped blade, 2.8mm C-0380

Diamond Phaco Knife, straight, keratome shaped blade, 3.0mm C-0390

Diamond Phaco Knife, angled, keratome shaped blade, 3.0mm C-0400

Diamond Phaco Knife, straight, keratome shaped blade, 3.2mm C-0410

Diamond Phaco Knife, angled, keratome shaped blade, 3.2mm C-0420



Phaco Diamond Knife, angled, keratome shaped blade & 4 cutting edges, 2.5mm C-0430

Phaco Diamond Knife, angled, keratome shaped blade & 4 cutting edges, 2.8mm C-0440

Phaco Diamond Knife, angled, keratome shaped blade & 4 cutting edges, 3.0mm C-0450

Phaco Diamond Knife, angled, keratome shaped blade & 4 cutting edges, 3.2mm C-0460

Phaco Diamond Knife, angled, keratome shaped blade & 4 cutting edges, trapezoid from 2.7mm to 3.2mm C-0470



Micrometer Diamond Knife, double-vertical & 45°-cutting edge with 200 microns base cutting edge

C-0480

Micrometer Diamond Knife, Thornton triple-vertical & 15°-cutting edge with 100 microns base cutting edge

C-0490



Diamond Knife with 30 ° cutting edge, 1.0mm wide

C-0500

Diamond Knife with 45 ° cutting edge, 1.0mm wide

C-0510

Diamond Knife with trie-faced, 1.0mm wide

C-0520

Diamond Knife Spear, 1.0mm wide

C-0530



Universal Trephine Handle

C-0550



Trephine Blade „short model”,
6.0mm - 9.5mm Ø (0.25mm increments)

C-0560 + size

Trephine Blade „long model”,
10.0mm - 15.0mm Ø (0.5mm increments)

C-0570 + size



Bard-Parker Scalpel Handle #3

C-0670



Bard-Parker Scalpel Handle #9

C-0680



Chuck Scalpel Handle round, serrated,
10.0cm

C-0690

Bard-Parker Scalpel blade, sterile,
disposable PK/100 #11

C-0700

Bard-Parker Scalpel blade, sterile,
disposable PK/100 #15

C-0710

Bard-Parker Scalpel blade, sterile,
disposable PK/100 #65

C-0711

Bard-Parker Scalpel blade, sterile,
disposable PK/100 #64

C-0712


Bard-Parker Scalpel blade, sterile,
disposable PK/100 #69

C-0714



Luntz-Dodick Trabeculectomy Punch,
1.0 x 0.5mm

C-0720



Kelly Descemet Membrane Punch,
1.0 x 0.75mm

C-0730



Gass Corneoscleral Punch,
1.5 x 0.75mm

C-0740



Soft IOL Cutter,
6.0mm long serrated blades,
for cutting soft IOL's

C-0750



IOL Acrylic Lens Cutter,
titanium

C-0760

a1 medical HOOKS,
HOOKS, SPATULAS,
MANIPULATORS,
FIXATION RINGS,
KERATOMES, LOOPS, MARKERS,
MARKING PEN,
PHACO CHOPPERS,
PREASSURE PLATES, SPOONS,










SINGLE HOOK



Choose a favorite color:

Get in touch with our newly designed, patented handle, which has been developed to enhance the performance of your clients! Feel the difference, with our new ergonomic handle. The titanium base material is the best there is. It is half the weight of stainless steel, causes less trauma and is antimagnetic. It helps to reduce fatigue during extended operations thanks to its ultralight weight. Our patented handle is ergonomically round and has perfect balance during surgery. Its projecting grooves not only provide an excellent grip, but the handle is also easy to maintain and sterilize. And finally, to help you save on costs, we have designed the handle so that the instrument tips can be replaced.

Total length 11.0cm

| | | |
|---|--|-------------|
| | Titanium Handle only, for Single Hook System, used with tip variations | HS-0001TSX* |
|  | Kuglen Iris Hook and Lens Manipulator, push-pull model, angled | HS-0000TSX* |
| | Replacement tip HS-0000TS-TIP | A-1** |
|  | Rosen Phaco Splitter, blunt tips, 10.0mm angled | HS-0010TSX* |
| | Replacement tip HS-0010TS-TIP | A-2** |
|  | Rosen Phaco Splitter, blunt tips, 10.0mm angled 60° | HS-0020TSX* |
| | Replacement tip HS-0020TS-TIP | A-3** |
|  | Haefliger Phaco cleaver, 45° angled, right | HS-0030TSX* |
| | Replacement tip HS-0030TS-TIP | A-4** |
|  | Nagahara Phaco-Chopper, 1.25mm tip, offset blade, angled 60°, right | HS-0040TSX* |
| | Replacement tip HS-0040TS-TIP | A-5** |

* replace "X" with either B (for blue), M (for magenta), Y (for yellow), G (for green), N (for neutral)

** Will be factory assembled on a returned hook handle. Please contact your a1-medical dealer for details.



Nagahara Phaco-Chopper, 1.5mm tip, offset blade, angled 60°, right

HS-0050TSX*

Replacement tip HS-0050TS-TIP

A-6**



Pisacano Nucleus Rotator, vertical flat disc

HS-0060TSX*

Replacement tip HS-0060TS-TIP

A-7**



Pisacano Nucleus Rotator, horizontal flat disc

HS-0070TSX*

Replacement tip HS-0070TS-TIP

A-8**



Bechert Nucleus Rotator, blunt forked tip, Y-horizontal, 7.0mm angled

HS-0080TSX*

Replacement tip HS-0080TS-TIP

A-9**



Drysdale Nucleus Manipulator

HS-0090TSX*

Replacement tip HS-0090TS-TIP

A-10**

Drysdale Nucleus Manipulator, extra delicate

HS-0091TSX

Replacement tip HS-0091TS-TIP

A-39**



Lester IOL-Manipulator, straight

HS-0110TSX*

Replacement tip HS-0110TS-TIP

A-11**

* replace "X" with either B (for blue), M (for magenta), Y (for yellow), G (for green), N (for neutral)

** Will be factory assembled on a returned hook handle. Please contact your a1-medical dealer for details.



| | | |
|---|--|-------------|
|  | Lester IOL-Manipulator, angled, 10.0mm | HS-0120TSX* |
| | Replacement tip HS-0120TS-TIP | A-12** |
|  | Rentsch IOL Manipulator, straight | HS-0130TSX* |
| | Replacement tip HS-0130TS-TIP | A-13** |
|  | Sinskey II Lens Manipulating Hook, Ø 0.2mm, blunt tip, straight | HS-0140TSX* |
| | Replacement tip HS-0140TS-TIP | A-14** |
|  | Sinskey II Lens Manipulating Hook, Ø 0.2mm, blunt tip, angled | HS-0150TSX* |
| | Replacement tip HS-0150TS-TIP | A-15** |
|  | Sinskey Lens Manipulating Hook, Ø 0.25mm, blunt tip, straight | HS-0160TSX* |
| | Replacement tip HS-0160TS-TIP | A-16** |
|  | Sinskey Lens Manipulating Hook, Ø 0.25mm, blunt tip, angled | HS-0170TSX* |
| | Replacement tip HS-0170TS-TIP | A-17** |
|  | Fenzl Lens Manipulating Hook, Ø 0.15mm tip, straight | HS-0180TSX* |
| | Replacement tip HS-0180TS-TIP | A-18** |
|  | Fenzl Lens Manipulating Hook, Ø 0.15mm tip, angled | HS-0190TSX* |
| | Replacement tip HS-0190TS-TIP | A-19** |

* replace "X" with either B (for blue), M (for magenta), Y (for yellow), G (for green), N (for neutral)

** Will be factory assembled on a returned hook handle. Please contact your a1-medical dealer for details.)



| | | |
|--|--|-------------|
| | Maltzman-Fenzl Lens Manipulator Hook, straight | HS-0200TSX* |
| | Replacement HS-0200TS-TIP | A-20** |
| | Maltzman-Fenzl Lens Manipulator Hook, angled | HS-0210TSX* |
| | Replacement tip HS-0210TS-TIP | A-21** |
| | Barraquer Iris Spatula, smooth Ø 0.5mm, angled, 12.0mm | HS-0220TSX* |
| | Replacement tip HS-0220TS-TIP | A-22** |
| | Barraquer Iris Spatula, smooth Ø 2.5mm, angled, 12.0mm | HS-0230TSX* |
| | Replacement tip HS-0230TS-TIP | A-23** |
| | Osher Y-Hook Lens Manipulator, straight | HS-0240TSX* |
| | Replacement tip HS-0240TS-TIP | A-24** |
| | Osher Y-Hook Lens Manipulator, angled | HS-0250TSX* |
| | Replacement tip HS-0250TS-TIP | A-25** |
| | Ligabue Nucleus Manipulator & Splitter, 12.0mm curved arm flat sandplasted oive shaped 90° rotated | HS-0256TSX* |
| | Replacement tip HS-0256TS-TIP | A-45** |

* replace "X" with either B (for blue), M (for magenta), Y (for yellow), G (for green), N (for neutral)










** Will be factory assembled on a returned hook handle. Please contact your a1-medical dealer for details.



DOUBLE HOOK

You can either choose a tip from our selection below or choose a tip according to your needs..

Total length 12.0cm

| | | | |
|---|---|--|-------------|
| | | Titanium Handle only, for Double Hook System, used with tip variations | HS-0002TSX* |
| 1. | 2. | | |
|  |  | Rowen Phaco Spatula & Beichert Nucleus Rotator | HS-0260TSX* |
| | | 1. Replacement tip HS-0260TS-Tip Rowen | A-26** |
| | | 2. Replacement tip HS-0260TS-Tip Beichert | A-27** |
|  |  | Nucleus Rotator & Dividing Hook | HS-0270TSX* |
| | | 1. Replacement tip HS-0270TS-Tip Rotator | A-28** |
| | | 2. Replacement tip HS-0270TS-Tip Dividing | A-29** |
|  |  | Nucleus Rotator (push-pull) & Dividing Hook | HS-0280TSX* |
| | | 1. Replacement tip HS-0280TS-Tip Push-Pull | A-30** |
| | | 2. Replacement tip HS-0280TS-Tip Dividing | A-31** |
|  | | Chang Phaco Chopper, for left hand | HS-0310TSX* |
| | | 1. Replacement tip HS-0310TS-Tip Microfinger | A-32** |
| | | 2. Replacement tip HS-0310TS-Tip Quick-Chop | A-33** |
|  | | Chang Phaco Chopper, for right hand | HS-0320TSX* |
| | | 1. Replacement tip HS-0320TS-Tip Microfinger | A-34** |
| | | 2. Replacement tip HS-0320TS-Tip Quick-Chop | A-35** |
|  | | Masket Phaco Spatula | HS-0330TSX* |
| | | 1. Replacement tip HS-0330TS-Tip Splitter | A-36** |
| | | 2. Replacement tip HS-0330TS-Tip Spatula | A-3** |

* replace "X" with either B (for blue), M (for magenta), Y (for yellow), G (for green), N (for neutral)

** Will be factory assembled on a returned hook handle. Please contact your a1-medical dealer for details.



Haefliger Phaco cleaver, 45° angled,
right

HS-0340



Phaco Manipulator, finger shaped,
curved tip

HS-0350



Nagahara Phaco-Chopper, 1.25mm tip, offset blade,
angled 60°, right

HS-0360



Nagahara Phaco-Chopper, 1.5mm tip, offset blade,
angled 60°, right

HS-0370

For Ceramic Coating please add "C"



Lund Iris Hook,
14.0mm tip, 45° angled

HS-0380



Dardenne Iris Hook,
10.0mm, angled shaft

HS-0390



Miyoshi Nucleus Divider,
straight

HS-0400



Miyoshi Nucleus Divider,
left

HS-0410



Miyoshi Nucleus Divider,
right

HS-0420



Neuhann Nucleus Rotator,
10.0mm angled tip, blunt forked tip

HS-0430



Neuhann Nucleus Spatula / Phaco Instrument,
curved spatulated tip with 45° shaft

HS-0440



Manipulator, Spatula 0.50mm x 0.25mm, round, Ø 0.25mm,
angled

HS-0450



Barraquer Iris Spatula, round smooth, Ø 0.5mm, angled, 12.0mm

HS-0460



Barraquer Iris Spatula, round smooth, Ø 2.5mm, angled, 12.0mm

HS-0470



Hirschman Iris Hook, blunt, angled

HS-0480



Osher Y-Hook Lens Manipulator, straight

HS-0490



Osher Y-Hook Lens Manipulator, angled

HS-0500



Deitz ICL Manipulator, blunt forked tip Y-Type, angled

HS-0510



Maloney IOL Manipulator, valuted shank

HS-0520



Koch Stop & Chop Manipulator, 1.5mm tip, for „Stop & Chop technique“.

HS-0521

For Ceramic Coating please add "C"



Manipulator, Spatula 1.0 x 0.25mm,
angled

HS-0530



Knolle Nucleus Spatula, 0.5mm wide,
straight

HS-0540



Hunkeler IOL Manipulator, with Ø 0.3mm semi-ball tip,
angled

HS-0550



Hirschman Lens Spatula, flat, 0.5mm wide, multiple noches

HS-0560



Koch Nucleus Spatula, sharp & smooth tip,
gently curved

HS-0570



Koch Nucleus Spatula, duckbill tip,
gently curved

HS-0571

For Ceramic Coating please add "C"



Kuglen Iris Hook and Lens Manipulator, push-pull model, angled

HS-0580

Kuglen Iris Hook, delicate tip

HS-0581



Rosen Phaco Splitter, wedge shaped blunt tip, inferior edges, angled 10.0mm

HS-0590




Rosen Phaco Splitter, wedge shaped blunt tip, inferior edges, angled 10.0mm, 60°

HS-0600



Pisacano Nucleus Rotator, vertical flat disc

HS-0610



Pisacano Nucleus Rotator, horizontal flat disc

HS-0620



Pisacano Nucleus Rotator, 11.0mm angled, olive round tip

HS-0621



Bechert Nucleus Rotator, blunt forked tip, Y-horizontal, 7.0mm angled

HS-0630



Drysdale Nucleus Manipulator, paddle shaped tip

HS-0640

For Ceramic Coating please add "C"



Lester IOL-Manipulator, Ø 0.25mm,
straight

HS-0650



Lester IOL-Manipulator, Ø 0.25mm,
10.0mm angled

HS-0660



Rentsch IOL Manipulator, blunt,
straight

HS-0670




Sinsky II Lens Manipulating Hook, Ø 0.2mm, blunt tip,
straight

HS-0680




Sinsky II Lens Manipulating Hook, Ø 0.2mm, blunt tip,
angled

HS-0690



Sinsky Lens Manipulating Hook, Ø 0.2mm, blunt tip,
angled, 180° reversed

HS-0691



Sinsky Lens Manipulating Hook, Ø 0.20mm blunt tip,
straight

HS-0700



Sinsky Lens Manipulating Hook, Ø 0.25mm blunt tip,
angled

HS-0710

For Ceramic Coating please add "C"



Fenzl Lens Manipulating Hook, Ø 0.15mm tip, straight

HS-0720



Fenzl Lens Manipulating Hook, Ø 0.15mm tip, angled

HS-0730



Maltzman-Fenzl Lens Manipulator Hook, Ø 0.15mm, V-shaped tip, straight

HS-0740



Maltzman-Fenzl Lens Manipulator Hook, Ø 0.15mm, V-shaped tip, angled

HS-0750



Kagawa Phaco-Chopper, offset blade, angled 90°, right


HS-0760

For Ceramic Coating please add "C"



Graether Collar Button Micro Iris Hook,
straight

HS-0770



Graether Collar Button Micro Iris Hook,
angled

HS-0780



Barraquer Iris Spatula, with flat disc,
angled 8.0mm

HS-0790



Barraquer Iris Spatula, with flat disc,
angled 10.0mm

HS-0800



Barraquer Iris Spatula, with flat disc,
angled 12.0mm

HS-0810



Barraquer Iris Spatula, with flat disc,
angled 15.0mm

HS-0820

For Ceramic Coating please add "C"



Bonn Micro Iris Hook, extra delicate, blunt tip, straight

HS-0830



Lieberman Micro Finger Nucleus Manipulator, for right hand

HS-0840



Lieberman Micro Finger Nucleus Manipulator, for left hand

HS-0850



Gimbel Nucleus Spatula, round, smooth, Ø 0.25mm, angled 30 °

HS-0860



Fukasaku Snapper Hook, with blunt tip & front notch

HS-0870




Fukasaku Snapper Hook, angled

HS-0880



LASIK Spatula rounded anterior & flat posterior surface with dissecting tip

HS-0890



Akahoshi Nucleus Sustainer, used with F-2850

HS-0895

For Ceramic Coating please add "C"



Shepard Iris Hook, blunt, 0.5mm hook, straight

HS-0900



Bonn Micro Iris Hook, blunt, 0.3mm hook, extra delicate

HS-0910



Kuglen Iris Hook & Lens Manipulator, push-pull model, straight

HS-0920



Kuglen Lens Manipulator push-pull, "Clover Leaf" model, straight

HS-0930



Kuglen Iris Hook & Lens Manipulator, push-pull model, angled 10.0mm

HS-0940



Jaffe Lens Manipulating Hook, Ø 0.15mm blunt tip, straight

HS-0950



Jaffe Lens Manipulating Hook, Ø 0.15mm blunt tip, angled

HS-0960

For Ceramic Coating please add "C"



Sinskey II Lens Manipulating Hook, Ø 0.2mm, blunt tip, straight

HS-0970




Sinskey II Lens Manipulating Hook, Ø 0.2mm, blunt tip, angled

HS-0980



Sinskey Lens Manipulating Hook, Ø 0.25mm, blunt tip, straight

HS-0990



Jaffe Lens Spatula, smooth, 0.5mm wide blade, angled 10.0mm

HS-1000



Obstbaum Synechia spatula, 0.25mm wide blade, 10.0mm angled

HS-1010



Clayman Nucleus Rotator, 0.5mm wide, flat Spatula, angled

HS-1020



DALK Corneal Dissector, 12.0mm long spatula, flat posterior surface, blunt, beveled tip

HS-1022

For Ceramic Coating please add "C"



Rowen Phaco Spatula & Bechert Nucleus Rotator

HS-1030



Nucleus Rotator & Dividing Hook

HS-1040



Nucleus Rotator (push-pull) & Dividing Hook

HS-1050



Pfaffl Femto Flaplifter

HS-1051



LASIK Flap Lifter
Combines a modified Sinskey Hook & a bi-convex spatula with beveled notches

HS-1053

For Ceramic Coating please add "C"



Castroviejo Spatula, double ended,
10.0mm & 15.0mm long

HS-1060



Chang Phaco Chopper, a combination of a modified MicroFinger,
and a pointed „Quick Chop“ tip,
for left hand

HS-1090



Chang Phaco Chopper, a combination of a modified MicroFinger,
and a pointed „Quick Chop“ tip,
for right hand

HS-1100



Masket Phaco Spatula, a combination of the Rosen Phaco Splitter
and a 12.0mm long round spatula with disc shaped tip

HS-1110

For Ceramic Coating please add "C"



Culler Iris Spatula,
1.0mm wide

HS-1120

Culler Iris Spatula,
2.0mm wide

HS-1130



Bangerter Iris Spatula,
1.0 x 0.25mm

HS-1140

Bangerter Iris Spatula,
1,5 x 0.25mm

HS-1150

Bangerter Iris Spatula,
2,0 x 0.25mm

HS-1160

For Ceramic Coating please add "C"



Giessen Orbital Spatula

HS-1170



Hamburg Orbital Spatula,
18.0mm x 12.0mm

HS-1180



Schepens Orbital Spatula

HS-1250

For Ceramic Coating please add "C"



Arruga Globe Retractor, cup shaped,
length 14.5cm

HS-1260



Wells Enucleation Spoon

HS-1270



Culler Lens Spoon,
sterling silver tip

HS-1280

For Ceramic Coating please add "C"



Paton Spatula & Spoon

HS-1290



Lewis Lens Loop,
small

HS-1300



Lewis Lens Loop,
large

HS-1310



New Orleans Lens Loop

HS-1320



Wilder Lens Loop,
15.0mm x 4.0mm smooth loop

HS-1330



Snellen Lens Loop,
15.0mm x 4.0mm, serrated

HS-1340



Arlt Lens Loop,
serrated, 3.0mm wide

HS-1350



Keener-Arlt Lens Loop, serrated, 3.0mm wide,
strong curved

HS-1360



Kansas Nucleus Loop, serrated, 3.0mm wide

HS-1370

For Ceramic Coating please add "C"



Evisceration Spoon, bunge small,
6.0mm HS-1460

Evisceration Spoon, bunge small,
8.0mm HS-1470

Evisceration Spoon, bunge large,
10.0mm HS-1480

Evisceration Spoon, bunge large,
12.0mm HS-1490

For Ceramic Coating please add "C"



Meyerhoefer Chalazions Curette,
size 00, Ø 1.0mm HS-1500

Meyerhoefer Chalazions Curette,
size 0, Ø 1.5mm HS-1510

Meyerhoefer Chalazions Curette,
size 1, Ø 1.7mm HS-1520

Meyerhoefer Chalazions Curette,
size 2, Ø 2.0mm HS-1530


Meyerhoefer Chalazions Curette,
size 3, Ø 2.5mm HS-1540

Meyerhoefer Chalazions Curette,
size 3.5, Ø 3.0mm HS-1545

Meyerhoefer Chalazions Curette,
size 4, Ø 3.5mm HS-1550

Meyerhoefer Chalazions Curette,
size 5, Ø 4.0mm HS-1560

For Ceramic Coating please add "C"




MARKING
PATTERN

Bores Axis Marker,
used with Mendez Degree Gauge

HS-1570

Toric IOL Marker
two 2.55mm radial marks extending to 11.0mm

HS-1571




MARKING
PATTERN

Mendez Degree Gauge 12.0mm,
0 - 180° in 10° increments

*HS-1580


For Ceramic Coating please add "C"

*Not available in Ceramic Coating



Lu-Mendez Fixation Ring,
with blunt teeth, 11.0mm

*HS-1590



Fine-Thornton Fixation Ring,
Ø 16.0mm, 16 blunt teeth, 15°

HS-1600

Hofmann-Thornton Fixation Ring,
Ø 13.0mm, 16 blunt teeth, 15°

HS-1610

For Ceramic Coating please add "C"

*Not available in Ceramic Coating



Globe Fixation Ring


HS-1620



Lindstrom Incision Marker,
used with Mendez gauge, 7.0mm optic zone

*HS-1630

For Ceramic Coating please add "C"
* Not available in Ceramic Coating



MARKING
PATTERN

Anis Radial Marker,
8 blades

HS-1640

Anis Radial Marker,
12 blades

HS-1650

Anis Radial Marker,
16 blades

HS-1660



MARKING
PATTERN

Radial Marker

12 blades

HS-1663

16 blades

HS-1664



MARKING
PATTERN

Dulaney LASIK Corneal Marker

HS-1670

For Ceramic Coating please add "C"



Dhayet LASIK Corneal Marker

HS-1690



Mendez LASIK Corneal Marker

*HS-1700

For Ceramic Coating please add "C"
* Not available in Ceramic Coating



Hoffer Zone Marker,
3.0mm

HS-1710

Hoffer Zone Marker,
3.25mm

HS-1720

Hoffer Zone Marker,
3.5mm

HS-1730

Hoffer Zone Marker,
3.75mm

HS-1740

Hoffer Zone Marker,
4.0mm

HS-1750

Hoffer Zone Marker,
4.25mm

HS-1760

Hoffer Zone Marker,
4.5mm

HS-1770

Hoffer Zone Marker,
4.75mm

HS-1780

For Ceramic Coating please add "C"



Hoffer Zone Marker,
5.0mm

HS-1790

Hoffer Zone Marker,
5.5mm

HS-1800

Hoffer Zone Marker,
6.0mm

HS-1810

Hoffer Zone Marker,
7.0mm

HS-1820

Hoffer Zone Marker,
8.0mm

HS-1830

Hoffer Zone Marker,
9.0mm

HS-1840



Bores Optic Zone Marker,
3.25mm

HS-1850

Bores Optic Zone Marker,
3.5mm

HS-1860

For Ceramic Coating please add "C"



Bore's Optic Zone Marker,
3.75mm

HS-1870

Bore's Optic Zone Marker,
4.0mm

HS-1880

Bore's Optic Zone Marker,
4.25mm

HS-1890

Bore's Optic Zone Marker,
4.5mm

HS-1900

Bore's Optic Zone Marker,
4.75mm

HS-1910

Bore's Optic Zone Marker,
5.0mm

HS-1920

Bore's Optic Zone Marker,
6.0mm

HS-1930

Bore's Optic Zone Marker,
7.0mm

HS-1940

For Ceramic Coating please add "C"



Borek Optic Zone Marker,
8.0mm

HS-1950

Borek Optic Zone Marker,
9.0mm

HS-1960



Lasek Epithelial Trepine,
8.0mm

HS-1970

Lasek Epithelial Trepine,
9.0mm

HS-1980

Lasek Epithelial Trepine,
10.0mm

HS-1990

For Ceramic Coating please add "C"



Camellin Style, LASEK Alcohol Well,
Ø 8.5mm

HS-2000


Camellin Style, LASEK Alcohol Well,
Ø 9.5mm

HS-2010



Buratto LASIK Flap Protector

HS-2020



MARKING
PATTERN

Fechtner Trabeculectomy Marker

HS-2021

For Ceramic Coating please add "C"



Membran-Peeler, blunt,
length 13.0cm

HS-2030



Membran-Peeler, blunt,
angled 110°,
length 13.0cm

HS-2040



Membran-Peeler, blunt,
angled 130°,
length 13.0cm

HS-2050



Membrane Hook,
length 13.0cm

HS-2060



Membrane Spatula,
length 13.0cm

HS-2070



Membrane Scratcher,
length 13.0cm

HS-2080

For Ceramic Coating please add "C"



Wadsworth Todd Kauter,
Length 14.0cm

*HS-2090



Eckardt Pressure Plate, 23G/ 0.6mm,
for Vitrectomy System, left hand, used with F-1540,
length 9.5cm

HS-2100

Eckardt Pressure Plate, 23G/ 0.6mm,
for Vitrectomy System, right hand, used with F-1540,
length 9.5cm

HS-2110

For Ceramic Coating please add "C"
* Not available in Ceramic Coating



Maloney Intra-Operative Keratometer,
titanium

HS-2120




Marking Pen, sterile,
„single use“ (Box of 5 psc.)

HS-2130

FORCEPS
MIC FORCEPS
TITANIUM FORCEPES
BIPOLAR FORCEPS
ERGO GRIP FORCEPS
TC FORCEPS




GENERAL INFORMATION




Pplatform Forceps

The usable area of the platform with guaranteed parallel closure. In this way a secure grip on the sutures can be ensured.




Pierse Type Forceps

These forceps are ideal for tissue fixation, manipulation of limbal area, corneal lip and minimize trauma to the tissue.




Ring Forceps

These ring forceps have been designed for perfect grasping soft tissue & handling of the conjunctiva. Additionally it provides an optimal counter-pressor function.



Tothed Forceps

These delicate suture forceps are ideal for fixation & manipulating of conjunctiva, limbal area or other delicate tissue & grasping fine suture tying.



Capsulorhexis Forceps

These forceps with triangular jaws, small platform and ultra fine grasping tips are optimal to perform the capsulorhexis.

...sharp tip is used as a cystotome.

...blunt tip is used to grip, lift and pull the iris.

ULTRA LIGHT ERGRO GRIP

...Perfect balance for the surgeon's use

Every surgical operation requires good skill and sensitivity. Operations with longer instruments under the microscope, in particular require thoroughness and precision work. Long operations are especially challenging for the surgeon. In order to facilitate the work of the surgeon, a1 medical has developed an ergonomically shaped and light-weight line of instruments.

The "Ultra-Light Ergo Grip". At a1 medical, our micro instruments are manufactured only from the finest stainless steel; their protruding ribs give you a perfect grip, and they can be easily cleaned. Their ergonomic design is not only pleasing to the eye but also offers perfect balance in the hand. Particularly noteworthy is the price-performance ratio as compared to titanium. Instruments made of titanium usually cost almost twice as much as stainless steel instruments.

Due to their excellent design, "Ultra-Light Ergo Grip" instruments weigh hardly any more than those made of titanium and are thus a good cost alternative.

Of course, our "Ultra-Light Ergoline" is also available with features such as TC coating and ceramic coating.

ADVANTAGES AT A GLANCE

- provides perfect balance during surgery
- ergonomic & attractive design
- excellent grip
- high quality stainless steels
- easy to maintain & sterilize
- saves costs





With guide pin & stop pin, Ø 8.0mm, length 11.5cm



Ergo Grip Tennant Tying Forceps, straight,
with 6.0mm platform, for 9.0-11.0 sutures

F-4800




Ergo Grip Tennant Tying Forceps, curved,
with 6.0mm tying platform, for 9.0-11.0 sutures,

F-4805




Ergo Grip Kelman-McPherson Forceps, angled,
7.5mm long smooth jaws
(ideal for implanting IOL)

F-4810




Ergo Grip McPherson Suturing Forceps, angled,
with 5.0mm tying platform

F-4820




Ergo Grip Gaskin Fragment Forceps, angled 11.0mm,
smooth jaws

F-4840




Ergo Grip Bechert Forceps,
angled, 10.0mm

F-4850



Ergo Grip Faulkner Lens Holding Forceps, angled,
with 11.0mm cone shaped jaws,
rounded surfaces & edges

F-4860



Ergo Grip Livernois Lens Forceps, angled,
with 5.0mm cone-shaped jaws,
flat inside surfaces & rounded edges

F-4870

For Ceramic Coating please add "C"



With guide pin & stop pin, Ø 8.0mm, length 11.5cm

| | | |
|--|---|---------------|
| | <p>Ergo Grip Steinert Direct Action Forceps, angled, with 7.0mm long cone-shaped jaws</p> | <p>F-4880</p> |
| | <p>Ergo Grip Spaleck Acrylic IOL Implantation Forceps, angled, 15.0mm</p> | <p>F-4890</p> |
| | <p>Ergo Grip Buratto Lens Folding Forceps, convex jaws</p> | <p>F-4900</p> |
| | <p>Ergo Grip Utrata Capsulorhexis Forceps, angled, triangular shaped tips</p> | <p>F-4910</p> |
| | <p>Ergo Grip Lehner Capsulorhexis Forceps, curved, triangular shaped tips</p> | <p>F-4920</p> |
| | <p>Ergo Grip Jaffe Capsulorhexis Forceps, angled, blunt, olive shaped tips</p> | <p>F-4930</p> |
| | <p>Ergo Grip Nevyas Capsulorhexis Forceps, angled, cystotome shaped tips</p> | <p>F-4940</p> |
| | <p>Ergo Grip Bonn Forceps, straight, with 6.0mm tying platform, 1x2 teeth, 0.12mm 90° angle</p> | <p>F-4950</p> |

For Ceramic Coating please add "C"

TC COATING

TC-Coating



When coated with tungsten carbide, tips of instruments last significantly longer since tungsten carbide is much harder than stainless steel. This allows an atraumatic and save handling of delicate tissue, vessels and fine needles.



With guide pin & stop pin, Ø 8.0mm, length 11.5cm




Ergo GripTennant Tying Forceps, straight, with 6.0mm tying platform, for 9.0-11.0 sutures

F-4800TC



Ergo Grip Tennant Tying Forceps, curved, with 6.0mm tying platform, for 9.0-11.0 sutures

F-4805TC



Ergo Grip Bonn Forceps, straight, with 5.0mm tying platform, 1x2 teeth, 0.12mm 90° angle


F-4950TC

For Ceramic Coating please add "C"




| | | |
|--|--|---------|
| | Harms Colibri Forceps, with tying platform, 1x2 teeth, 0.12mm, 45° angle, length 7.5cm | F-0010 |
| | Same in titanium | F-0010T |
| | Barraquer Corneal Utility Forceps, with tying platform, 1x2 teeth, 0.12mm, 90° angle, length 7.5cm | F-0020 |
| | Same in titanium | F-0020T |
| | Barraquer Corneal Utility Forceps, with tying platform, 1x2 teeth, 0.3mm, 90° angle, length 7.5cm | F-0030 |
| | Same in titanium | F-0030T |
| | Barraquer Corneal Utility Forceps, with tying platform, 1x2 teeth, 0.4mm, 90° angle, length 7.5cm | F-0040 |
| | Same in titanium | F-0040T |
| | Barraquer Corneal Utility Forceps, 1x2 teeth, 0.12mm, 90° angle, length 7.5cm | F-0050 |
| | Same in titanium | F-0050T |
| | Iris Botvin Forceps, with 1x2 teeth, double curved, length 7.5cm | F-0060 |
| | Same in titanium | F-0060T |

For Ceramic Coating please add "C"




Micro Colibri Forceps, with tying platform,
1x2 teeth, 0.12mm, 45° angle,
length 7.5cm

F-0070



Micro Colibri Forceps, with tying platform,
1x2 teeth, 0.3mm, 45° angle,
length 7.5cm

F-0080




Troutman Corneal Fixation Forceps,
1x2 teeth, 0.12mm, 45° angle,
length 7,5.cm

F-0100



Barraquer Corneal Utility Forceps, with tying platform,
1x2 teeth, 0.12mm, 90° angle,
length 7.5cm

F-0120




Barraquer Corneal Utility Forceps,
1x2 teeth, 0.12m, 90° angle,
length 7.5cm

F-0140




Suture Forceps (Zurich model) straight, with tying platform,
1x2 teeth, 0.12mm, 45° angle,
length 7.5cm

F-0160




Corneal Colibri Forceps, with tying platform,
1x2 teeth, 0.12mm, 90° angle,
length 7.5cm

F-0230




Barraquer Forceps for removing sutures, with tying platform,
length 7.5cm

F-0250




Polack Double Corneal Forceps,
with 2x 1x2 teeth, 0.12mm,
length 7.5cm

F-0260



Fixation Colibri Forceps, with tying platform,
1x2 teeth, 0.12mm, 45° angle,
length 7.5cm


F-0270



Corneal Colibri Utility Forceps, with tying platform,
1x2 teeth, 0.3mm, 45° angle,
length 7.5cm


F-0280

For Ceramic Coating please add "C"




Gill Iris Forceps,
criss-cross serrated tips,
length 7.5cm

F-0410




Micro Suture Forceps (Bonn model), straight,
5.0mm tying platform
length 7.5cm

F-0440



Iris Forceps (Bonn model), straight, fine serrated,
length 7.5cm

F-0450



Yoshitomi-Neuhann Lens Holding Forceps „diamond dusted tips“
with one flat & one convex surface, short model,
length 7.5cm

F-0490



Bonn Suturing Forceps, straight, with tying platform, 1x2 teeth, 0.12mm, 90° angle, length 7.5cm

F-0510



Iris Forceps Bonn model, straight, with 1x2 teeth, 0.12mm, 90° angle, length 7.5cm

F-0520




Corneal Tying Forceps, straight, with extra fine tying platform, length 8.5cm

F-0560



Same in titanium

F-0560T




Corneal Tying Forceps, curved, with extra fine tying platform, length 8.5cm

F-0570

Same in titanium

F-0570T



Corneal Tying Forceps, angled, with extra fine 7.0mm tying platform, length 8.5cm

F-0580

Same in titanium

F-0580T

For Ceramic Coating please add "C"



Corneal Fixation Forceps, straight, with tying platform, 1x1 Pierse teeth, 0.1mm, length 8.5cm

F-0590



Same in titanium

F-0590T



Corneal Fixation Forceps, straight, with tying platform, 1x1 Pierse teeth, 0.2mm, length 8.5cm

F-0600



Same in titanium

F-0600T



Corneal Fixation Forceps, straight, 1x1 Pierse teeth, 0.2mm, length 8.5cm

F-0610



Same in titanium

F-0610T



Corneal Fixation Forceps, straight, 1x1 Pierse teeth, 0.3mm, length 8.5cm

F-0620



Same in titanium

F-0620T



Corneal Suturing Forceps, straight, with tying platform, 1x2 teeth, 0.12mm, 90° angle, length 8.5cm

F-0650



Same in titanium

F-0650T



Corneal Suturing Forceps, curved, with tying platform, 1x2 teeth, 0.12mm, 45° angle, length 8.5cm


F-0680




Same in titanium

F-0680T

For Ceramic Coating please add "C"



| | | |
|--|---|----------------|
| | <p>Corneal Fixation Forceps, straight, with tying platform, 1x1 Pierse teeth, 0.2mm, length 8.5cm</p> | <p>F-0740</p> |
| | <p>Same in titanium</p> | <p>F-0740T</p> |
| | <p>Corneal Fixation Forceps, angled, with 1x1 Pierse teeth, 0.2mm, length 8.5cm</p> | <p>F-0770</p> |
| | <p>Same in titanium</p> | <p>F-0770T</p> |
| | <p>Colibri Corneal Fixation Forceps, with tying platform, 1x1 Pierse teeth, 0.1mm, length 8.5cm</p> | <p>F-0790</p> |
| | <p>Same in titanium</p> | <p>F-0790T</p> |
| | <p>Colibri Corneal Utility Forceps, straight, with tying platform, 1x2 teeth, 0.12mm, 90° angle, length 8.5cm</p> | <p>F-0810</p> |
| | <p>Same in titanium</p> | <p>F-0810T</p> |
| | <p>Colibri Corneal Utility Forceps, with tying platform, 1x2 teeth, 0.12mm, 45° angle, length 8.5cm</p> | <p>F-0820</p> |
| | <p>Same in titanium</p> | <p>F-0820T</p> |



| | | |
|---|---|---------------|
|  | <p>Paufique Suture Forceps, with 5.0mm tying platform, standard jaws, length 8.5cm</p> | <p>F-0830</p> |
|  | <p>Paufique Suture Forceps, with 5.0mm tying platform, very delicate jaws, length 8.5cm</p> | <p>F-0840</p> |
|  | <p>Suture Forceps (St. Martin model), with 5.0mm tying platform, 1x2 teeth, 0.12mm, 90° angle, length 8.5cm</p> | <p>F-0850</p> |
|  | <p>Suture Forceps (St. Martin model), with 5.0mm tying platform, 1x2 teeth, 0.3mm, 90° angle, length 8.5cm</p> | <p>F-0860</p> |
|  | <p>Suture Forceps (Paufique model), with 5.0mm tying platform, 1x2 teeth, 0.12mm, 45° angle, length 8.5cm</p> | <p>F-0870</p> |
|  | <p>Paufique Suture Forceps, with 5.0mm tying platform, 1x2 teeth, 0.3mm, 45° angle, length 8.5cm</p> | <p>F-0880</p> |
|  | <p>Manhattan Eye & Ear Suturing Forceps, with 5.0mm tying platform, 1x2 teeth, 0.5mm, 45° angle, length 8.5cm</p> | <p>F-0890</p> |

For Ceramic Coating please add "C"




Paufigue Suture Forceps, straight, with 5.0mm tying platform, 1x2 teeth, 0.7mm, 45° angle, length 8.5cm

F-0900



Sauer Suturing Forceps, straight, 1x2 teeth, 0.5mm, 45° angle, length 8.5cm

F-0930



Suture Forceps (Tuebingen model), with 5.0mm tying platform very delicate, length 8.5cm

F-0940

For Ceramic Coating please add "C"




McPherson Tying Forceps straight, with 5.0mm tying platform, length 8.5cm

F-0950




McPherson Tying Forceps angled, with 5.0mm tying platform, length 8.5cm

F-0960




Kelman-McPherson Tying Forceps, 7.5mm, angled, with 5.0mm platform, length 8.5cm

F-0970




Bechert-McPherson Tying Forceps, 10.0mm, angled, with 5.0mm tying platform, length 8.5cm

F-0980



Sheets McPherson IOL Forceps, 10.0mm, angled, with 5.0mm tying platform, length 8.5cm

F-0990



Utrata Capsulorhexis Forceps, angled, triangular shaped tips, length 8.5cm


F-1000



Lehner Capsulorhexis Forceps, angled, triangular shaped tips, length 8.5cm

F-1010

For Ceramic Coating please add "C"




Gimbel Capsulorhexis Forceps, angled, only 3.0mm opening, triangular shaped tips, length 8.5cm

F-1020



Lehner Capsulorhexis Forceps, curved, criss-cross serrated, triangular shaped tips, length 8.5cm

F-1030



McPherson Corneal Forceps, straight, with 5.0mm tying platform, 1x2 teeth, 0.12mm, 90° angle, length 8.5cm

F-1040




Bonn Tying Forceps, straight, with 5.0mm tying platform, for 9.0-11.0 sutures, length 9.5cm

F-1050




Bonn Tying Forceps, curved, with 5.0mm tying platform, for 9.0-11.0 sutures, length 9.5cm

F-1060



Bonn Tying Forceps, angled, with 5.0mm tying platform, for 9.0-11.0 sutures, length 9.5cm


F-1070



Bonn Forceps straight, long model, with 5.0mm tying platform, 1x2 teeth, 0.12mm, 90° angle, length 9.5cm


F-1090

For Ceramic Coating please add "C"




Thorpe Conjunctiva Fixation Forceps, straight, with tying platform, with 2x3 teeth, 0.6mm, 90° angle, length 9.5cm

F-1120




Lester Fixation Forceps, straight, with 2x3 teeth, 0.6mm, 90° angle, length 9.5cm

F-1130




Castroviejo Tying Forceps, straight, with 5.0mm tying platform, length 11.0cm

F-1140



Kelman McPherson Tying Forceps, 7.5mm angled, with tying platform, length 11.0cm

F-1150



Utrata Capsulorhexis Forceps, angled, triangular shaped tips, length 11.0cm



F-1160

For Ceramic Coating please add "C"




| | | |
|--|--|---------------|
| | <p>Castroviejo Suturing Forceps, with 6.0mm tying platform, 1x2 teeth, 0.12, 45° angle, length 11.0cm</p> | <p>F-1190</p> |
| | <p>Castroviejo Suturing Forceps, with 6.0mm tying platform, 1x2 teeth, 0.3mm, 45° angle, length 11.0cm</p> | <p>F-1200</p> |
| | <p>Castroviejo Suturing Forceps, with 6.0mm tying platform, 1x2 teeth, 0.5mm, 45° angle, length 11.0cm</p> | <p>F-1210</p> |
| | <p>Castroviejo Suturing Forceps, with 6.0mm tying platform, 1x2 teeth, 0.9mm, 45° angle, length 11.0cm</p> | <p>F-1220</p> |
| | <p>Snider-Castroviejo Sutures Forceps, with 6.0mm tying platform, 1x2 teeth, 0.12mm, 45° angle, with lock, length 11.0cm</p> | <p>F-1240</p> |
| | <p>Castroviejo Suturing Forceps, 1x2 teeth, 0.3mm, 45° angle, length 11.0cm</p> | <p>F-1250</p> |
| | <p>Castroviejo Suturing Forceps, 1x2 teeth, 0.5mm, 45° angle, length 11.0cm</p> | <p>F-1260</p> |
| | <p>McPherson Corneal Forceps, with 5.0mm tying platform, 1x2 teeth, 0.4mm, 90° angle, length 11.0cm</p> | <p>F-1290</p> |

For Ceramic Coating please add "C"




Castroviejo Suturing Forceps, angled, with 6.0mm tying platform, 1x2 teeth, 0.12mm, 45° angle, length 11.0cm

F-1300




Castroviejo Suturing Forceps angled, with 6.0mm tying platform, 1x2 teeth, 0.3mm, 45° angle, length 11.0cm

F-1310



Castroviejo Suturing Forceps angled, with 6.0mm tying platform, 1x2 teeth, 0.5mm, 45° angle, length 11.0cm

F-1320




Colibri Forceps, with tying platform, 1x2 teeth, 0.12mm, 90° angle, length 11.0cm

F-1380




Stern-Castroviejo Fixation Forceps, with 6.0mm tying platform, 1x2 teeth, 0.5mm, 45° angle, with lock, length 11.0cm

F-1400



Moody Fixation Forceps, right, 1x2 teeth, 0.5mm, 45° angle, with lock, length 11.0cm


F-1410



Moody Fixation Forceps, left, 1x2 teeth, 0.5mm, angle, with lock, length 11.0cm


F-1420

For Ceramic Coating please add "C"




Harm's Tying Forceps, straight, with 6.0mm tying platform, length 11.0cm

F-1510




Pierce Corneal Forceps, straight, 1x1 teeth 0.10mm, length 10.5cm

F-1520




Harm's Tying Forceps, curved, with 6.0mm tying platform, length 11.0cm

F-1530



Removal Cannula Forceps for 23G, length 11.0cm

F-1540



Dodick Nucleus Cracker, cross-action, longitudinally serrated paddle shaped jaws, length 11.0cm

F-1550



Castroviejo Suturing Forceps, straight, with tying platform, 1x2 teeth, 0.12mm, 45° angle, length 11.0cm

F-1560

Same in titanium

F-1560T



Castroviejo Suturing Forceps, straight, with tying platform, 1x2 teeth, 0.3mm, 45° angle, length 11.0cm

F-1570

Same in titanium

F-1570T



Castroviejo Suturing Forceps, straight, with tying platform, 1x2 teeth, 0.5 mm, 45° angle, length 11.0cm

F-1580

Same in titanium

F-1580T



Castroviejo Suturing Forceps, with 6.0mm tying platform, 1x2 teeth, 0.9mm, 45° angle, length 11.0cm

F-1590


Same in titanium

F-1590T



Stern-Castroviejo Fixation Forceps, 1x2 teeth, 0.5mm, 45° angle, with lock, length 11.0cm


F-1610



Moody Fixation Forceps, curved, with 1x2 teeth, 0.5mm, 45° angle, with lock, left, length 11.0cm


F-1620

For Ceramic Coating please add "C"



Moody Fixation Forceps, curved,
1x2 teeth, 0.5mm 45°angle, with lock, right,
length 11.0cm

F-1630



Colibri Forceps,
1x2 teeth, 0.12mm, 90° angle,
length 11.0cm

F-1670



McPherson Suture Forceps straight, with 5.0mm tying platform,
length 10.5cm

F-1680



Same in titanium

F-1680T




Helveston Tying Forceps,
standard smooth jaws,
length 10.5cm

F-1690



Livernois St Lens Forceps straight, with 0.5mm cone-shaped jaws,
flat inside surface & rounded edges,
length 10.5cm

F-1700




Kelman-McPherson Forceps, angled,
very delicate 7.5mm long smooth jaws,
length 10.5cm

F-1710

Same in titanium


F-1710T

For Ceramic Coating please add "C"




McPherson Suturing Forceps, angled,
with 5.0mm tying platform,
length 10.5cm

F-1720




Same in titanium

F-1720T



Gaskin Fragment Forceps, 9.0mm angled,
with very delicate smooth jaws
length 10.5cm

F-1730



Gaskin Fragment Forceps, 11.0mm angled,
smooth jaws,
length 10.5cm

F-1740



Bechert Forceps, 10.0mm angled,
length 10.5cm

F-1750



Faulkner Lens Holding Forceps, angled,
11.0mm cone shaped jaws with, rounded surfaces & edges,
length 10.5cm


F-1760



Livernois Lens Forceps, angled,
5.0mm cone-shaped jaws with flat inside surfaces & rounded edges,
length 10.5cm

F-1770

For Ceramic Coating please add "C"




Steinert Direct Action Forceps, angled, with 7.0mm long cone-shaped jaws, length 10.5cm

F-1780



Haberle-McPherson Lens Holding Forceps, 7.5mm angled, smooth jaws, diagonal groove inside tips, length 10.5cm

F-1790




Spaleck Acrylic IOL Implantation Forceps, angled, 15.0mm, length 10.5cm

F-1800



Buratto Lens Folding Forceps, convex jaws, length 10.5cm

F-1810



Utrata Capsulorhexis Forceps, angled, triangular shaped tips, length 10.5cm

F-1820

Same in titanium

F-1820T




Lehner Capsulorhexis Forceps, curved, triangular shaped tip, length 10.5cm

F-1830

Same in titanium

F-1830T



Jaffe Capsulorhexis Forceps, angled, blunt, olive shaped tips, length 10.5cm

F-1840



| | | |
|--|--|----------------|
| | <p>Nevyas Capsulorhexis Forceps, angled, cystotome shaped tips, length 10.5cm</p> | <p>F-1850</p> |
| | <p>LASIK Flap Forceps, curved disc shaped & serrated jaws, length 10.5cm</p> | <p>F-1860</p> |
| | <p>Kurstin Flap Stretching Forceps, length 10.5cm</p> | <p>F-1870</p> |
| | <p>Bonn Forceps, straight, with tying platform, 1x2 teeth, 0.12mm, 90° angle, length 10.5cm</p> | <p>F-1910</p> |
| | <p>Same in titanium</p> | <p>F-1910T</p> |
| | <p>Bonn Forceps, straight, with tying platform, 1x2 teeth, 0.3mm, 90° angle, length 10.5cm</p> | <p>F-1920</p> |
| | <p>Same in titanium</p> | <p>F-1920T</p> |
| | <p>McPherson Corneal Forceps, straight, with 5.0mm tying platform, 1x2 teeth, 0.4mm 90° angle, length 10.5cm</p> | <p>F-1930</p> |
| | <p>Same in titanium</p> | <p>F-1930T</p> |
| | <p>Kelman-McPherson Forceps, angled 7.5mm, with tying platform, 1x2 teeth, 0.3mm 90° angle, length 10.5cm</p> | <p>F-1940</p> |
| | <p>Same in titanium</p> | <p>F-1940T</p> |
| | <p>Fechtner Conjunctiva Forceps Ø 0.8mm ring, length 10.5cm</p> | <p>F-2421</p> |

For Ceramic Coating please add "C"




Arita Meibomian Gland Compressor,
length 10.5cm

F-1971



Endothelial / Descemetorhexis Forceps, extra-delicate,
reversed triangular tips, angled 75°,
length 10.5cm

F-1975



Jaffe Forceps, straight, delicate jaw,
with 6.0mm tying platform,
length 11.0cm

F-1990




Jaffe Forceps, curved, delicate jaw,
with 6.0mm tying platform,
length 11.0cm

F-2000

For Ceramic Coating please add "C"



MIC Forceps for small incisions




MIC Utrata Capsulorhexis Forceps, triangular grasping tips, extremely thin shanks, for 2.2mm incision, controlled opening, length 10.5cm

F-1960




MIC Corydon Capsulorhexis Forceps, cystotome shaped tips, extremely thin shanks, for 2.2mm incision, controlled opening, length 10.5cm

F-1970




MIC GIANNETTI Capsulorhexis Forceps, curved, for 2.0mm incision, controlled opening, length 10.5cm

F-1976



MIC Utrata Capsulorhexis Forceps, triangular shaped tips, 14.0mm, angled, very delicate shanks, length 8.5cm

F-1003



MIC Jaffe Capsulorhexis Forceps, olive shaped blunt tips, 14.0mm, angled, very delicate shanks, length 8.5cm

F-1004

For Ceramic Coating please add "C"



MIC Forceps for small incisions

- | | | |
|--|---|---------------|
| | <p>MIC Corydon Capsulorhexis Forceps, curved, 12.0mm, very delicate shanks, cystotome shaped tips, for 2.0mm incision, controlled opening, Ø 6.0mm, length 11.0cm</p> | <p>F-3190</p> |
| | <p>MIC Corydon Capsulorhexis Forceps, curved, 12.0mm, very delicate shanks, cystotome shaped tips, for 1.8mm incision, controlled opening, Ø 6.0mm, length 11.0cm</p> | <p>F-3191</p> |
| | <p>MIC GIANNETTI Capsulorhexis Forceps, curved, for 2.0mm incision, controlled opening, Ø 6.0mm, length 11.0cm</p> | <p>F-3200</p> |



- | | | |
|--|--|---------------|
| | <p>MIC Utrata Capsulorhexis Forceps, 13.0mm, angled, extra delicate shanks, triangular shaped tips, Ø 6.0mm, length 10.5cm</p> | <p>F-3151</p> |
|--|--|---------------|

For Ceramic Coating please add "C"



Capsulorhexis Forceps with markings on jaw



Corydon Capsulorhexis Forceps 12.0mm angled, cystotome tips,
marking 2.5.mm & 5.0mm, Ø 8.0mm,
length 11.0cm

F-2981




MIC Forceps for small incisions



MIC Lehner Capsulorhexis Forceps, 12.0mm, angled,
for 2.2mm incision, triangular tips, delicate shaft, controlled opening,
Ø 8.0mm, length 11.0cm

F-2982



GIANNETTI MIC Capsulorhexis Forceps, curved,
for 2.0mm incision, controlled opening,
Ø 8.0mm,length 11.0cm

F-2983



| | | |
|--|---|---------------|
| | <p>Harms Tying Forceps (Tuebingen model), straight, with 6.0mm tying platform, length 10.0cm</p> | <p>F-2460</p> |
| | <p>McPherson Suture Forceps, straight, with 5.0mm tying platform, length 10.0cm</p> | <p>F-2470</p> |
| | <p>McGregor Conjunctiva Fixation Forceps, with 5.0mm tying platform, cross-serrated tips, length 10.0cm</p> | <p>F-2480</p> |
| | <p>Harms Tying Forceps (Tuebingen model), curved, with 6.0mm tying platform, length 10.0cm</p> | <p>F-2490</p> |
| | <p>Harms Tying Forceps (Tuebingen model), angled, with 6.0mm tying platform, length 10.0cm</p> | <p>F-2500</p> |
| | <p>Suture Forceps Kelman-McPherson, 7.5mm angled, length 10.0cm</p> | <p>F-2510</p> |
| | <p>Utrata Capsulorhexis Forceps, 12.0mm angled, triangular shaped tips, length 10.0cm</p> | <p>F-2520</p> |

For Ceramic Coating please add "C"



Iris Forceps Bonn model, straight, with 6.0mm tying platform, 1x2 teeth, 0.12mm, 90° angle, length 10.0cm

F-2540



Bonn Iris Forceps, 1x2 teeth, 0.12mm, 90° angle, length 10.0 cm

F-2550



Castroviejo Suturing Forceps, straight, with 6.0mm tying platform, 1x2 teeth, 0.12mm, 45° angle, length 10.0cm

F-2560



Castroviejo Suturing Forceps, straight, with 6.0mm tying platform, 1x2 teeth, 0.3mm, 45° angle, length 10.0cm

F-2570




Castroviejo Suturing Forceps, straight, with 6.0mm tying platform, 1x2 teeth, 0.5mm, 45° angle, length 10.0cm

F-2580



Castroviejo Suturing Forceps, straight, with tying platform, 1x2 teeth, 0.9mm, 45° angle, length 10.0cm

F-2590




Troutman Tying Forceps, straight, with extra fine 6.0mm tying platform, for 8-0 to 11-0 sutures, length 10.0cm F-2800



Troutman Tying Forceps, curved, with extra fine 6.0mm tying platform, for 8-0 to 11-0 sutures, length 10.0cm F-2810



Troutman Suturing Forceps, straight, with tying platform, 1x2 teeth, 0.12mm, 90° angle, length 10.0cm F-2820



Akahoshi Phaco Pre-Chopper, with 4.5mm delicate jaws & sharp tips, length 10.0cm F-2850




Ernest-McDonald II IOL Inserting Forceps, 5.0mm straight jaws, length 10.0cm F-2860

For Ceramic Coating please add "C"




Calf III IOL Forceps,
folder for foldable IOL's,
length 10.0cm

F-2830



Ernest-Calf Nucleus Cracking Forceps (cross-action),
extremely thin paddle,
length 10.0cm

F-2840



Livernois Pickup & Folding Forceps, angled,
L-shaped jaws with concave inside surfaces,
length 10.5cm

F-2880

For Ceramic Coating please add "C"



With guide pin & stop pin, Ø 8.0mm, length 11.0cm

| | | |
|---|--|---------|
|  | Tennant Tying Forceps, straight, with 6.0mm platform, for 9.0-11.0 sutures | F-2890 |
|  | Same in titanium | F-2890T |
|  | Tennant Tying Forceps, curved, with 6.0mm tying platform, for 9.0-11.0 sutures | F-2900 |
|  | Same in titanium | F-2900T |
|  | Kelman-McPherson Tying Forceps, angled, with 8.0mm tying platform for 9.0-11.0 sutures | F-2910 |
|  | Same in titanium | F-2910T |
|  | Gaskin Fragment Forceps, angled, 11.0mm, ultra thin smooth jaws | F-2920 |
|  | Same in titanium | F-2920T |
|  | Faulkner Lens Holding Forceps, 11.0mm, angled | F-2930 |
|  | Utrata Capsulorhexis Forceps, angled, 12.0mm, triangular tips | F-2940 |
|  | Same in titanium | F-2940T |
|  | Jaffe Capsulorhexis Forceps, angled, 12.0mm, blunt shaped tips | F-2960 |
|  | Same in titanium | F-2960T |

For Ceramic Coating please add "C"



With guide pin & stop pin, Ø 8.0mm, length 11.0cm



Lehner Utrata Capsulorhexis Forceps, angled, 12.0mm, triangular shaped tips

F-2970



Same in titanium

F-2970T



Bonn Forceps straight, with 6.0mm tying platform, 1x2 teeth, 0.12mm

F-2990



Same in titanium

F-2990T



Castroviejo Suture Forceps, straight, with tying platform, 1x2 teeth, 0.12mm 45° angle

F-3000




Same in titanium

F-3000

For Ceramic Coating please add "C"



With guide pin & stop pin, Ø 6.0mm, length 11.0cm




Tennant Tying Forceps straight, with 6.0mm tying platform, for 9.0-11.0 sutures

F-3100




Tennant Tying Forceps curved, with 6.0mm tying platform, for 9.0-11.0 sutures

F-3110



Kelman-McPherson Tying Forceps, angled, with 8.0mm tying platform, for 9.0-11.0 sutures

F-3120




Gaskin Fragment Forceps, angled, 11.0mm, with ultra thin smooth jaws

F-3130




Utrata Capsulorhexis Forceps, angled, 12.0mm, triangular shaped tips

F-3140



Lehner Utrata Capsulorhexis Forceps, angled, 12.0mm, triangular shaped tips

F-3150



Jaffe Utrata Capsulorhexis Forceps, angled, 12.0mm, blunt shaped tips



F-3160



Bonn Forceps straight, with tying platform, 1x2 teeth, 0.12mm, 90° angle

F-3170

For Ceramic Coating please add "C"



Tennant Tying Forceps, straight, 6.0mm tying platform for 9-0 to 1-0 sutures, length 10.0cm

F-3210




Tennant Tying Forceps, curved, with 6.0mm tying platform, for 9-0 to 1-0 sutures, length 10.0cm

F-3220



Utrata Capsulorhexis Forceps, angled, 12.0mm, triangular tips, length 10.0cm

F-3250



Lehner Capsulorhexis Forceps, curved, 12.0mm, triangular shaped tips, length 10.0cm

F-3260



Tennant Suturing Forceps, straight, with tying platform, 1x2 teeth, 0.12mm 90° angle, length 10.0cm

F-3270




Corneoscleral Forceps, straight, with tying platform, 1x2 teeth, 0.12mm, 45° angle, length 10.0cm

F-3280



Colibri Forceps, curved, with tying platform, 1x2 teeth, 0.12mm, 45° angle, length 10.0cm

F-3290



Colibri Forceps, straight, with tying platform, 1x2 teeth, 0.12mm, 45° angle, length 10.0cm

F-3300

For Ceramic Coating please add "C"



Tennant Tying Forceps, straight, with 6.0mm tying platform, for 9.0-1.0 sutures, length 10.5cm

F-3310




Tennant Tying Forceps, curved, with 6.0mm tying platform, for 9.0-1.0 sutures, length 10.5cm

F-3320



Utrata Capsulorhexis Forceps, angled, 12.0mm, triangular shaped tips, length 10.5cm

F-3330



Jaffe Capsulorhexis Forceps, angled, 12.0mm, blunt shaped tips, length 10.5cm

F-3331



O'Gawa Tying Forceps, straight, with 5.5mm tying platform, length 10.5cm

F-3340



O'Gawa Tying Forceps, curved, with 5.5mm tying platform, length 10.5cm

F-3350



Utrata Capsulorhexis Forceps, angled, 13.0mm, triangular shaped tips, length 10.5cm


F-3360



Lehner Capsulorhexis Forceps, curved, 13.0mm, triangular shaped tips, length 10.5cm

F-3370

For Ceramic Coating please add "C"



Designed for the „Phaco Chopper Technique“



Salvitti Akahoshi Combo Pre-Chopper,
length 12.5cm

F-3390



Akahoshi Combo Pre-Chopper,
length 12.5cm

F-3400




Akahoshi Combo Pre-Chopper, Universal,
length 12.5cm

F-3410




Akahoshi Combo Pre-Chopper, cross-action,
length 12.5cm

F-3420




Bishop Harmon Fixation Forceps, serrated, 0.5mm tips, length 8.5cm

F-3530




Bishop Harmon Fixation Forceps, serrated, 0.8mm tips, length 8.5cm

F-3540




Bishop Harmon Fixation Forceps, serrated, 1.0mm tips, length 8.5cm

F-3550




Bishop-Harmon Dressing Forceps, criss-cross serrated, 0.5mm tips, length 8.5cm

F-3560




Bishop-Harmon Dressing Forceps, criss-cross serrated, 0.8mm tips, length 8.5cm

F-3570



Bishop-Harmon Forceps(Bonn model), with 5.0mm tying platform, 1x2 teeth, 0.12mm, 90° angle, length 8.5cm

F-3590



Bishop-Harmon Forceps (Bonn model), with 5.0mm tying platform, 1x2 teeth, 0.3mm, 90° angle, length 8.5cm

F-3600

For Ceramic Coating please add "C"



Bishop-Harmon Forceps (Bonn model), with 5.0mm tying platform,
1x2 teeth, 0.4mm, 90° angle,
length 8.5cm

F-3610



Bishop-Harmon Forceps (Bonn model), with 5.0mm tying platform,
1x2 teeth, 0.7mm, 90° angle,
length 8.5cm

F-3620



Bishop-Harmon Tissue Forceps,
with 1x2 teeth, 0.12mm, 90° angle,
length 8.5cm

F-3630




Bishop-Harmon Tissue Forceps,
with 1x2 teeth, 0.3mm, 90° angle,
length 8.5cm

F-3640



Bishop-Harmon Tissue Forceps,
with 1x2 teeth, 0.5mm, 90° angle,
length 8.5cm

F-3650



Bishop-Harmon Tissue Forceps,
with 1x2 teeth, 0.6mm, 90° angle,
length 8.5cm

F-3660



Bishop-Harmon Tissue Forceps,
with 1x2 teeth, 0.8mm, 90°, angle,
length 8.5cm


F-3670

For Ceramic Coating please add "C"




Lester Conjunctiva Fixation Forceps, straight,
with 1x2 teeth, 0.6mm, 90° angle,
length 9.5cm

F-3690




Lester Conjunctiva Fixation Forceps, straight,
with 2x3 teeth, 0.6mm, 90° angle,
length 9.5cm

F-3700




Thorpe Conjunctiva Fixation Forceps, with 5.0mm tying platform,
with 2x3 teeth, 0.6mm, 90° angle,
length 9.5cm

F-3710



Bonaccolto Utility Forceps, 1.2mm wide,
length 11.0cm

F-3720



Bonaccolto Utility Forceps, 1.7mm wide,
length 11.0cm

F-3730



Elschnig Skleral Fixation Forceps, straight,
with 1x2 teeth, 0.6mm, 45° angle,
length 11.0cm


F-3750



Elschnig Skleral Fixation Forceps, straight,
1x2 teeth, 0.8mm, 45° angle,
length 11.0cm


F-3760

For Ceramic Coating please add "C"




Remky Superior Rectus Forceps,
length 11.0cm

F-3780




Lister Conjunctiva Forceps,
length 11.0cm

F-3810



Bonaccolto Forceps (French model), straight,
serrated jaws, 0.8mm,
length 11.0cm

F-3820



Bonaccolto Forceps (French model), curved,
serrated jaws, 1.0mm,
length 11.0cm

F-3830



Jayle Forceps, straight, with 7.0mm tying platform,
1x2 teeth, 0.3mm, 45° angle,
length 11.0cm


F-3840



Jayle Forceps, straight, with 7.0mm tying platform,
1x2 teeth, 0.8mm, 45° angle,
length 11.0cm


F-3860

For Ceramic Coating please add "C"




Dressing Forceps, serrated, straight, 0.5mm, length 10.0cm

F-3870




Dressing Forceps, serrated, straight, 0.8mm, length 10.0cm

F-3880




Dressing Forceps, serrated, straight, 1.0mm, length 10.0cm

F-3890




Dressing Forceps, serrated, straight, 1.2mm, length 10.0cm

F-3900




Dressing Forceps, serrated, curved, 0.5mm, length 10.0cm

F-3910



Dressing Forceps, serrated, strong curved, 0.6mm, length 10.0cm


F-3920



Iris Forceps serrated, curved, 0.8mm, length 10.0cm


F-3930

For Ceramic Coating please add "C"




Dressing Forceps, serrated, strong curved, 0.8mm, length 10.0cm

F-3940




Dressing Forceps, serrated, curved, 1.0mm, length 10.0cm

F-3950




Dressing Forceps, serrated, strong curved, 1.0mm, length 10.0cm

F-3960




Dressing Forceps, serrated, curved, 1.2 mm, length 10.0cm

F-3970




Nugent Utility Forceps, curved shafts, 45° angle, with smooth 10.0mm jaws, length 10.0cm

F-3980




Iris Forceps, straight, with 1x2 teeth, 0.7mm, 90° angle, length 10.0cm

F-4000




Iris Forceps, straight, with 1x2 teeth, 1.0mm, 90° angle, length 10.0cm

F-4010



Iris Forceps, curved, with 1x2 teeth, 0.8mm, 90° angle, length 10.0cm


F-4020



Iris Forceps, curved, with 1x2 teeth, 1.0mm, 90° angle, length 10.0cm


F-4030

For Ceramic Coating please add "C"




Barraquer Fixation Forceps,
length 10.0cm

F-4040



Graefe Fixation Forceps, without lock,
length 11.0cm

F-4050



Graefe Fixation Forceps, with lock,
length 11.0cm

F-4060



Barraquer Cilia Forceps,
length 10.5cm

F-4070

For Ceramic Coating please add "C"



Barraquer Cilia Forceps, slightly angled with 6.0mm tying platform,
length 11.0cm

F-1980




Cilia Forceps, 2.0mm,
length 9.0cm

F-4080

Cilia Forceps, 3.0mm,
length 9.0cm

F-4090




Gradle Cilia Forceps,
smooth jaws,
length 9.5cm

F-4100

Gradle Cilia Forceps,
with positive gripping jaws,
length 9.5cm


F-4110

For Ceramic Coating please add "C"



Jeweler's Forceps #5, straight,
extra delicate

F-4120



Jeweler's Forceps #4, straight,
extra delicate pointed tips

F-4130




Jeweler's Forceps #7, curved,
extra delicate, pointed,
length 11.5 cm

F-4140



Jeweler's Forceps, curved,
extra pointed, 0.4mm,
length 11.5cm


F-4150



Jeweler's Forceps, angled,
extra delicate, pointed serrated tips,
length 11.5cm


F-4160

For Ceramic Coating please add "C"



Moorfields Conjunctiva Forceps,
standard model,
length 11.0cm

F-4170



Moorfield Conjunctiva Forceps,
fine model,
length 11.0cm

F-4180



Watzke Sleeve Spreader,
length 11.0cm

F-4200

For Ceramic Coating please add "C"



Shaaf Cilia Foreign Body Forceps,
length 9.5cm

F-4190



Ayer Chalazion Forceps,
round, Ø 8.0mm,
length 9.2cm

F-4210



Lambert Chalazion Forceps,
round, tip width: 11.0mm,
length 9.2cm

F-4220

Lambert Chalazion Forceps,
round, tip width: 15.0mm,
length 9.2cm

F-4230

For Ceramic Coating please add "C"



Desmarres Chalazion Forceps, 20.0mm,
length 9.2cm

F-4240

Desmarres Chalazion Forceps, 26.0mm,
length 9.2cm

F-4250

Desmarres Chalazion Forceps, 31.0mm,
length 9.2cm

F-4260



Erhardt Chalazion Forceps,
oval, Ø 20.0mm,
length 9.2cm

F-4231



Hunt Chalazion Forceps,
round, Ø 12.0mm,
length 9.7cm


F-4270

For Ceramic Coating please add "C"



Baird Chalazion Forceps, oval,
11.0 x 8.0mm, small,
length 9.5cm

F-4280



Francis Chalazion Forceps oval,
14.0 x 11.0mm I.D., medium,
length 9.5cm

F-4290



Wies Chalazion Forceps,
length 9.5cm

F-4300

For Ceramic Coating please add "C"



Desmarres cross-action Forceps, oval,
14.0mm x 11.0mm,
Length 9.5cm

F-4310



Lordan Chalazion Forceps,
length 9.5cm

F-4320



Berke Ptosis Forceps,
with slide lock, longitudinally grooved jaws,
20.0mm, length 10.5cm

F-4330

Berke Ptosis Forceps,
with slide lock, longitudinally grooved jaws,
27.0mm, length 10.5cm

F-4340

For Ceramic Coating please add "C"



Snellen Entropium Forceps, right,
with small solid blades, 28.0mm,
length 9.5cm

F-4350

Snellen Entropium Forceps, left,
with small solid blades, 28.0mm,
length 9.5cm

F-4360

Snellen Entropium Forceps, right,
with solid blades, 35.0mm,
length 9.5cm

F-4370

Snellen Entropium Forceps, left,
with solid blades, 35.0mm,
length 9.5cm

F-4380



Downs Lid Forceps,
length 9.8cm

F-4395

For Ceramic Coating please add "C"



Putterman Forceps, Lid Clamp,
with 6 pins, serrated jaws, sliding lock,
length 10.0cm

F-4390



Jameson Muscle Forceps, right,
with slide lock, 4 teeth & holes,
length 10.0cm

F-4400

Jameson Muscle Forceps, left,
with slide lock, 4 teeth & holes,
length 10.0cm

F-4405

Jameson Muscle Forceps, right,
with slide lock, 6 teeth & holes,
length 10.0cm

F-4410


Jameson Muscle Forceps, left,
with slide lock, 6 teeth & holes,
length 10.0cm

F-4415



Ocular Landers Lens Forceps,
length 10,3cm

F-4420



Scleral Plug Forceps,
length 10.3cm

F-4430

Designed to hold scleral plugs during insertion and removal,
cross-action handle



Bangerter Muscle Clamp,
length 10.5cm

F-4440

For Ceramic Coating please add "C"



Dieffenbach Bulldog Clamp, serrifine straight,
length 3.5cm F-4450


Dieffenbach Bulldog Clamp, serrifine curved,
length 3.5cm F-4460

Dieffenbach Bulldog Clamp, serrifine straight,
length 5.0cm F-4470

Dieffenbach Bulldog Clamp, serrifine curved,
length 5.0cm F-4480


Dieffenbach Bulldog Clamp, serrifine straight, with hollow rivet,
length 3.5cm F-4485

Dieffenbach Bulldog Clamp, serrifine curved, with hollow rivet,
length 3.5cm F-4487




Blaydes Implantation Forceps, without lock, very delicate narrow jaws angled 45°, length 11.0cm

F-4490



Clayman Lens Holding Forceps, without lock, length 11.0cm

F-4500




Clayman Lens Holding Forceps, without lock, length 11.0cm

F-4510




Kraff Lens Holding Forceps, grasp edge lens optic, length 11.0cm

F-4520



Bechert Lens Holding Forceps, gently curved, length 11.0cm


F-4530



Shepard Lens Forceps, without lock, length 11.0cm


F-4540

For Ceramic Coating please add "C"




Zaldivar Implantation Forceps, without lock,
length 11.0cm

F-4550



Tennant-Maumenee Colibri Forceps, with 5.0mm tying platform,
1x2 teeth, 0.12mm, 90° angle,
length 11.0cm

F-4560



Tennant-Maumenee Colibri Forceps, with 5.0mm tying platform,
1x2 teeth, 0.12mm, 45° angle,
length 11.0cm



F-4570



Inamura Capsulorhexis cross-action Forceps, titanium,
13.0mm curved, triangular shaped tips,
length 10.5cm


F-4575T

For Ceramic Coating please add "C"




Bechert Lens holding Forceps,
gently curved, without lock,
length 11.0cm

F-4580




Kratz Lens Holding Forceps,
gently curved, without lock,
length 11.0cm

F-4590




Blaydes Implantations Forceps,
without lock, very delicate narrow jaws, angled 45°,
length 11.0cm

F-4600



Shepard Lens Holding Forceps,
without lock, gently curved,
length 11.0cm

F-4610



Clayman Lens Holding Forceps,
without lock,
length 11.0cm

F-4620



Sinsky Lens Holding Forceps,
without lock,
length 11.0cm

F-4630

For Ceramic Coating please add "C"



Hartmann Mosquito Forceps, straight,
length 9.5cm

F-4640

Hartmann Mosquito Forceps, curved,
length 9.5cm

F-4650

Halsted Mosquito Forceps, straight,
length 14.5cm

F-4660

Halsted Mosquito Forceps, curved,
length 14.5cm

F-4670



Baby Jones Towel Clamp,
length 5.5cm

F-4680

For Ceramic Coating please add "C"



Backhaus towel clamp,
length 9.0cm

F-4690

Backhaus towel clamp,
length 13.5cm

F-4700



Bipolar Forceps straight,
slightly blunted tips, 0.25mm,
length 11.0cm

F-4710

Bipolar Forceps angled,
slightly blunted tips, 0.25mm,
length 11.0cm

F-4720

Bipolar Cord, autoclavable,
5.0m

F-4730

For Ceramic Coating please add "C"



McPherson Bipolar Forceps, straight,
slightly blunted tips, 0.25mm, isolated,
length 9.5cm

F-4740

McPherson Bipolar Forceps, angled,
slightly blunted tips, uninsulated,
length 9.5cm

F-4750

Bipolar Cord, autoclavable,
length 5.0m

F-4730

NEEDLE HOLDERS
ERGOGRIP NEEDLE HOLDERS
BLADE BREAKERS
TITANIUM NEEDLE HOLDERS
TC NEEDLE HOLDERS



ULTRA LIGHT ERGRO GRIP

...Perfect balance for doctors' use

Every surgical operation requires skill and good sensitivity. Operations with longer instruments under the microscope in particular require thoroughness and precision work. Long operations are especially challenging for the surgeon. In order to facilitate the work of the surgeon, a1 medical has developed an ergonomically shaped and light-weight line of instruments.

The "Ultra-Light Ergo Grip". At a1 medical, our micro instruments are manufactured only from the finest stainless steel; their protruding ribs give you a perfect grip, and they can be easily cleaned. Their ergonomic design is not only pleasing to the eye but also offers perfect balance in the hand. Particularly noteworthy is the price-performance ratio as compared to titanium. Instruments made of titanium usually cost almost twice as much as stainless steel instruments.

Due to their excellent design, "Ultra-Light Ergo Grip" instruments weigh hardly any more than those made of titanium and are thus a good cost alternative.

Of course, our "Ultra-Light Ergoline" is also available with features such as TC coating and ceramic coating.

ADVANTAGES AT A GLANCE

- provides perfect balance during surgery
- ergonomic & attractive design
- excellent grip
- high quality stainless steel
- easy to maintain & sterilize
- saves costs



For Ceramic Coating please add "C"



Ergo Grip Barraquer Needle Holder, straight, 10.0 x 1.0mm, smooth jaw, with lock, length 13.5cm N-0800

Ergo Grip Barraquer Needle Holder, straight, 10.0 x 1.0mm, smooth jaw, without lock, length 13.5cm N-0810

Ergo Grip Barraquer Needle Holder, curved, 10.0 x 1.0mm, smooth jaw, with lock, length 13.5cm N-0820

Ergo Grip Barraquer Needle Holder, curved, 10.0 x 1.0mm, smooth jaw, without lock, length 13.5cm N-0830



Ergo Grip Barraquer Needle Holder, straight, 10.0 x 0.6mm, N-0840
smooth jaw, with lock, length 13.5cm

Ergo Grip Barraquer Needle Holder, straight, 10.0 x 0.6mm, N-0850
smooth jaw, without lock, length 13.5cm

Ergo Grip Barraquer Needle Holder, curved, 10.0 x 0.6mm, N-0860
smooth jaw, with lock, length 13.5cm

Ergo Grip Barraquer Needle Holder, curved, 10.0 x 0.6mm, N-0870
smooth jaw, without lock, length 13.5cm

Ergo Grip Barraquer Needle Holder, curved, 10.0 x 0.5mm, N-0880
smooth jaw, with lock, length 13.5cm

Ergo Grip Barraquer Needle Holder, curved, 10.0 x 0.5mm, N-0890
smooth jaw, without lock, length 13.5cm

For Ceramic Coating please add "C"

TC-COATING

TC-Coating



When coated with tungsten carbide, tips of instruments last significantly longer since tungsten carbide is much harder than stainless steel. This allows a save handling of fine needles.



Ergo Grip Barraquer Needle Holder with TC-Coating, straight, 10.0 x 1.0mm, with lock, length 13.5cm

N-0800TC



Ergo Grip Barraquer Needle Holder, with TC-Coating, straight, 10.0 x 1.0mm, without lock, length 13.5cm

N-0810TC



Ergo Grip Barraquer Needle Holder with TC-Coating, curved, 10.0 x 1.0mm, with lock, length 13.5cm

N-0820TC



Ergo Grip Barraquer Needle Holder with TC-Coating, curved, 10.0 x 1.0mm, without lock, length 13.5cm

N-0830TC

For Ceramic Coating please add "C"



McPherson Needle Holder, straight, 10.0mm x 0.5mm, smooth jaw, with lock, 11.0cm N-0010

McPherson Needle Holder, straight, 10.0mm x 0.5mm, smooth jaw, without lock, 11.0cm N-0020

McPherson Needle Holder, curved, 10.0mm x 0.5mm, smooth jaw, with lock, 11.0cm N-0030

McPherson Needle Holder, curved, 10.0mm x 0.5mm, smooth jaw, without lock, 11.0cm N-0040



McPherson Needle Holder, straight, 10.0mm x 0.5mm, smooth jaw, with lock, 10.0cm N-0050

McPherson Needle Holder, straight, 10.0mm x 0.5mm, smooth jaw, without lock, 10.0cm N-0060

McPherson Needle Holder, curved, 10.0mm x 0.5mm, smooth jaw, with lock, 10.0cm N-0070

McPherson Needle Holder, curved, 10.0mm x 0.5mm, smooth jaw, without lock, 10.0cm N-0080

For Ceramic Coating please add "C"



Castroviejo Needle Holder, straight, 10.0mm x 1.0mm,
smooth jaw, with lock, length 14.0cm N-0090

Castroviejo Needle Holder, straight, 10.0mm x 1.0mm,
smooth jaw, without lock, length 14.0cm N-0100

Castroviejo Needle Holder, curved, 10.0mm x 1.0mm,
smooth jaw, with lock, length 14.0cm N-0110

Castroviejo Needle Holder, curved, 10.0mm x 1.0mm,
smooth jaw, without lock, length 14.0cm N-0120

Castroviejo Needle Holder, straight, 10.0mm x 0.75mm,
smooth jaw, with lock, length 14.0cm N-0130

Castroviejo Needle Holder, straight, 10.0mm x 0.75mm,
smooth jaw, without lock, length 14.0cm N-0140

Castroviejo Needle Holder, curved, 10.0mm x 0.75mm,
smooth jaw, with lock, length 14.0cm N-0150

Castroviejo Needle Holder, curved, 10.0mm x 0.75mm,
smooth jaw, without lock, length 14.0cm N-0160

For Ceramic Coating please add "C"



Kalt Needle Holder, straight, 10.0 x 1.0mm,
length 14.0cm

N-0250

Kalt Needle Holder, curved, 10.0 x 1.0mm,
length 14.0cm

N-0260

Kalt Needle Holder, straight, 12.0 x 1.8mm,
length 14.0cm

N-0270

For Ceramic Coating please add "C"



Barraquer Needle Holder (Baby Model), straight, 9.0 x 0.6mm, smooth jaw, with lock, length 10.0cm N-0280


Barraquer Needle Holder (Baby Model), straight, 9.0 x 0.6mm, smooth jaw, without lock, length 10.0cm N-0290

Barraquer Needle Holder (Baby Model), curved, 9.0 x 0.6mm, smooth jaw, with lock, length 10.0cm N-0300

Barraquer Needle Holder (Baby Model), curved, 9.0 x 0.6mm, smooth jaw, without lock, length 10.0cm N-0310

Barraquer Needle Holder, curved, 9.0 x 0.5mm, smooth jaw, with lock, length 12.0cm N-0320

Barraquer Needle Holder, curved, 9.0 x 0.5mm, smooth jaw, without lock, length 12.0cm N-0330



Barraquer Needle Holder, curved, 9.0 x 0.75mm,
smooth jaw, with lock, length 12.0cm N-0340


Barraquer Needle Holder, curved, 9.0 x 0.75mm,
smooth jaw, without lock, length 12.0cm N-0350

Barraquer Needle Holder, straight, 10.0 x 1.0mm,
smooth jaw, with lock, length 13.5cm N-0360

Barraquer Needle Holder, straight, 10.0 x 1.0mm,
smooth jaw, without lock, length 13.5cm N-0370

Barraquer Needle Holder, curved, 10.0 x 1.0mm,
smooth jaw, with lock, length 13.5cm N-0380

Barraquer Needle Holder curved, 10.0 x 1.0mm,
smooth jaw, without lock, length 13.5cm N-0390



Barraquer Needle Holder, straight, 10.0 x 0.6mm,
with lock, length 13.5cm N-0400

Barraquer Needle Holder, straight, 10.0 x 0.6mm,
without lock, length 13.5cm N-0410

Barraquer Needle Holder, curved, 10.0 x 0.6mm,
smooth jaw, with lock, length 13.5cm N-0420

Barraquer Needle Holder, curved, 10.0 x 0.6mm,
smooth jaw, without lock, length 13.5cm N-0430

Barraquer Needle Holder, curved, 10.0 x 0.5mm,
smooth jaw, with lock, length 13.5cm N-0440

Barraquer Needle Holder, curved, 10.0 x 0.5mm,
smooth jaw, without lock, length 13.5cm N-0450

Needle Holder in titanium upon request!

For Ceramic Coating please add "C"



Anis Needle Holder, straight, 7.0 x 0.4mm,
smooth jaw, with lock, length 12.0cm N-0460

Anis Needle Holder, straight, 7.0 x 0.4mm,
smooth jaw, without lock, length 12.0cm N-0470

Anis Needle Holder, curved, 7.0 x 0.4mm,
smooth jaw, with lock, length 12.0cm N-0480

Anis Needle Holder, curved, 7.0 x 0.4mm,
smooth jaw, without lock, length 12.0cm N-0490



Jakobi Needle Holder, straight, 7.0 x 0.5mm,
smooth jaw, with lock, length 12.0cm

N-0500

Jakobi Needle Holder, straight, 7.0 x 0.5mm,
smooth jaw, without lock, length 12.0cm

N-0510

Jakobi Needle Holder, curved, 7.0 x 0.5mm,
smooth jaw, with lock, length 12.0cm

N-0520

Jakobi Needle Holder, curved, 7.0 x 0.5mm,
smooth jaw, without lock, length 12.0cm

N-0530



Barraquer Needle Holder, French style, straight, 10.0 x 0.6mm
smooth jaw, without lock, length 12.0cm N-0540

Barraquer Needle Holder, French style, straight, 10.0 x 0.6mm
smooth jaw, with lock, length 12.0cm N-0550

Barraquer Needle Holder, French style, curved, 10.0 x 0.6mm
smooth jaw, without lock, length 12.0cm N-0560

Barraquer Needle Holder, French style, curved, 10.0 x 0.6mm
smooth jaw, with lock, length 12.0cm N-0570

Barraquer Needle Holder, French style, curved, 10.0 x 0.4mm
smooth jaw, without lock, length 13.5cm N-0620

Barraquer Needle Holder, French style, curved, 10.0 x 0.6mm
smooth, without lock, length 13.5cm N-0630

Barraquer Needle Holder, French style, curved, 10.0 x 0.6mm
smooth jaw, with lock, length 13.5cm N-0640

Barraquer Needle holder, French style, curved, 10.0 x 1.0mm
smooth jaw, without lock, length 13.5cm N-0650

For Ceramic Coating please add "C"



Castroviejo Blade Breaker and Holder,
with sliding lock,
length 13.5cm

N-0660



Swiss Model Blade Breaker & Holder, straight,
flat inside surface, length 12.0cm

N-0670

Swiss Model Blade Breaker & Holder, straight,
concave inside surface, length 12.0cm

N-0680

For Ceramic Coating please add "C"



Troutman-Chris Blade Breaker & Holder,
with concave/ convex inside surface,
length 9.5cm

*N-0690



Derf Needle Holder straight,
serrated jaws with groove,
length 12.5cm

N-0700

For Ceramic Coating please add "C"

* Not available in Ceramic Coating

SCISSORS
SUPER-CUT SCISSORS
RING SCISSORS
MICRO SCISSORS



GENERAL INFORMATION




Sharp scissors (Adventitia)

For fine trimming and secure suture cutting.



Blunt scissors (Dissection)

Blunt blades avoid possible damages during dissection



Super Cut Scissors

A razor sharp cutting edge on one blade and a micro-serration on the opposite blade assure a perfectly precise and non traumatic cut with minimal force.



Vannas blade

The vannas scissors are also applied in vascular surgery because due to their more delicate, smaller and shorter blades.



Castroviejo Universal Corneal Scissors, curved, blunt tips & small blades, cutting length 5.0mm, length 11.5cm S-0010

Castroviejo Universal Corneal Scissors, curved, blunt tips & medium blades, cutting length 6.5mm, length 11.5cm S-0020

Castroviejo Universal Corneal Scissors, curved, blunt tips & large blades, cutting length 8.0mm, length 11.5cm S-0030



Westcott Tenotomy Suture Scissors, curved, blunt tips & medium blades, cutting length 9.0mm, length 11.5cm S-0040

McPherson-Westcott Suture Scissors, curved, sharp tips & small blades, cutting length 6.0mm, length 11.5cm S-0050

McPherson-Westcott Suture Scissors, curved, blunt tips & small blades, cutting length 6.0mm, length 11.5cm S-0060

For Ceramic Coating please add "C"



Jaffe Stitch Scissors, curved,
sharp tips & medium blades, cutting length 8.0mm,
length 11.5cm

S-0070



Fine Stitch Scissors, straight,
sharp tips & medium blades, cutting length 11.0mm,
length 11.5cm

S-0080

Fine Stitch Scissors, curved,
sharp tips & medium blades, cutting length 11.0mm,
length 11.5cm

S-0090

For Ceramic Coating please add "C"



Vannas Scissors, straight,
sharp tips & 10.0mm blades, cutting length 8.0mm,
length 11.5cm

S-0100

Vannas Scissors, curved,
sharp tips & 10.0mm blades, cutting length 8.0mm,
length 11.5cm

S-0110



Castroviejo Corneoscleral Scissors, right,
small blades, cutting length 5.0mm,
length 10.0cm

S-0120

Castroviejo Corneoscleral Scissors, left,
small blades, cutting length 5.0mm,
length 10.0cm

S-0130

Castroviejo Corneoscleral Scissors, right,
medium blades, cutting length 7.0mm,
length 10.0cm

S-0140

Castroviejo Corneoscleral Scissors, left,
medium blades, cutting length 7.0mm,
length 10.0cm

S-0150

For Ceramic Coating please add "C"



Castroviejo Universal Corneal Scissors, curved,
blunt tips & small blades, cutting length 6.0mm,
length 10.0cm S-0160

Castroviejo Universal Corneal Scissors, strong curved,
blunt tips & small blades, cutting length 5.0mm,
length 10.0cm S-0161

Castroviejo Universal Corneal Scissors, curved,
blunt tips & medium blades, cutting length 10.0mm,
length 10.0cm S-0170

Castroviejo Universal Corneal Scissors, curved,
blunt tips & large blades, cutting length 10.0mm,
length 10.0cm S-0180



Castroviejo Keratoplasty Scissors, angled to side,
blunt tips & medium blades, cutting length 10.0mm,
length 10.0cm S-0190

For Ceramic Coating please add "C"



McPherson-Westcott Conjunctival Scissors, curved,
blunt tips & small blades, cutting length 5.0mm,
length 10.5cm S-0200

Westcott Tenotomy Scissors, curved,
blunt tips & medium blades, cutting length 8.0mm,
length 11.5cm S-0210

Westcott Tenotomy Scissors, curved right,
blunt tips & standard blades, cutting length 11.0mm,
length 11.5cm S-0382

Westcott Tenotomy Scissors, curved right,
blunt tips & standard blades, cutting length 14.0mm,
length 11.5cm S-0220

Westcott Tenotomy Scissors, curved left,
blunt tips & standard blades, cutting length 9.0mm,
length 11.5cm S-0230

Westcott Tenotomy Scissors, curved,
blunt tips & small blades, cutting length 5.0mm,
length 10.0cm S-0240

Westcott Tenotomy Scissors, curved,
pointed tips & standard blades, cutting length 14.0mm,
length 11.5cm S-0250

For Ceramic Coating please add "C"



Iris Scissors, straight, sharp tips & small blades
cutting length 6.0mm,
length 10.0cm

S-0260

Iris Scissors, curved, sharp tips & small blades,
cutting length 6.0mm,
length 10.0cm

S-0270



Noyes Iris Scissors, straight, sharp tips & small blades,
cutting length 10.0mm,
length 11.0cm

S-0280

Noyes Iris Scissors, lightly curved, sharp tips & small blades,
cutting length 10.0mm,
length 11.0cm

S-0290

Noyes Iris Scissors, straight, sharp tips & standard blades,
cutting length 12.5mm,
length 11.0cm

S-0300

Noyes Iris Scissors, lightly curved, sharp tips & standard blades,
cutting length 12.5mm,
length 11.0cm

S-0310

For Ceramic Coating please add "C"



Westcott Tenotomy Scissors, curved,
blunt tips & 13.0mm blades, cutting length 7.5mm,
length 14.0cm

S-0380

Westcott Tenotomy Scissors, curved,
blunt tips & 22.0mm blades, cutting length 11.0mm,
length 14.0cm

S-0390



Troutman-Katzin Corneal Transplant Scissors, strong curved, right,
pointed tips & small blades, cutting length 6.0mm,
length 10.0cm

S-0400

Troutman-Katzin Corneal Transplant Scissors, strong curved, left,
pointed tips & small blades, cutting length 6.0mm
length 10.0cm

S-0410

For Ceramic Coating please add "C"



Vannas Capsulotomy Scissors, straight,
sharp tips, 5.0mm blade,
length 8.5cm S-0420

Vannas Capsulotomy Scissors, straight,
blunt tips, 5.0mm blade,
length 8.5cm S-0430

Vannas Capsulotomy Scissors, curved,
sharp tips, 5.0mm blade,
length 8.5cm S-0440

Vannas Capsulotomy Scissors, curved,
blunt tips, 5.0mm blade,
length 8.5cm S-0450

Gills-Vannas Capsulotomy Scissors, straight,
sharp tips, 7.0mm blade,
length 8.5cm S-0460

Gills-Vannas Capsulotomy Scissors, curved,
sharp tips, 7.0mm blade,
length 8.5cm S-0470

Vannas Capsulotomy Scissors, angled forward,
sharp tips, 5.0mm blade,
length 8.5cm S-0480

For Ceramic Coating please add "C"



Gills-Vannas Capsulotomy Scissors, angled forward,
sharp tips, 7.0mm blade,
length 8.5cm

S-0490

Vannas Capsulotomy Scissors, angled to side,
sharp tips, 5.0mm blade,
length 8.5cm

S-0500

Gills-Welsh Capsulotomy Scissors, angled forward,
sharps tips, 11.0mm blade,
length 8.5cm

S-0501

Uribe-Stern Capsulotomy Scissors, angled forward,
concave curved blade,
length 8.5cm

S-0505



DeWecker Iridectomy Scissors, angled,
sharp tips & 7.0mm blade,
length 11.0cm

S-0510

DeWecker Iridectomy Scissors, angled,
blunt tips & 7.0mm blade,
length 11.0cm

S-0520

For Ceramic Coating please add "C"



Vannas Scissors, sharp, curved, cutting length 6.0mm,
Ø 6.0mm handle, length 11.5cm

S-0415

Vannas Scissors, blunt, curved, cutting length 6.0mm,
Ø 6.0mm handle, length 11.5cm

S-0416



Barraquer Iris Scissors,
blunt tips & 7.0mm blades

S-0530

Barraquer Iris Scissors,
pointed tips & 7.0mm blades

S-0540

Barraquer Iris Scissors,
pointed tips & 11.0mm blades

S-0550

Barraquer Iris Scissors,
blunt tips & 11.0mm blades

S-0560



Eye Scissors, straight, pointed tips, small blades,
14.0mm blades,
length 9.0cm S-0570

Eye Scissors, curved, pointed tips, small blades,
14.0mm blades,
length 9.0cm S-0580

Eye Scissors, curved, blunt tips, small blades,
14.0mm blades,
length 9.0cm S-0581

Eye Scissors, straight, pointed tips, large blades,
23.0mm blades,
length 9.0cm S-0590

Eye Scissors, curved, pointed tips, large blades,
24.0mm blades,
length 9.0cm S-0600

Eye Scissors, straight, pointed tips, standard model,
23.0mm blades,
length 10.5cm S-0610

Eye Scissors, straight, blunt tips, standard model,
23.0mm blades,
length 10.5cm S-0611

Eye Scissors, curved, pointed tips, standard model,
24.0mm blades,
length 10.5cm S-0620

For Ceramic Coating please add "C"



Eye Scissors, curved, blunt tips, standard model,
24.0mm blades,
length 10.5cm

S-0621

Eye Scissors, straight, pointed tips, standard model,
23.5mm blades,
length 9.5cm

S-0630

Eye Scissors, curved, pointed tips, big rings,
22.5mm blades,
length 9.5cm

S-0640

Eye Scissors, curved, blunt tips, big rings,
22,5mm blades,
length 9.5cm

S-0641



Stitch Scissors, ribbon type, curved,
pointed tips, 12.50mm blades,
length 11.0cm

S-0650

For Ceramic Coating please add "C"



Stitch Scissors, lightly curved,
pointed tips, 12.5mm blades,
length 9.5cm

S-0660




Stevens Tenotomy Scissors, straight,
blunt tips, 14.0mm blades,
length 11.5cm

S-0670

Stevens Tenotomy Scissors, curved,
blunt tips, 13.0mm blades,
length 11.5cm

S-0680

For Ceramic Coating please add "C"



Stevens Tenotomy Scissors, ribbon type, straight,
blunt tips, 13.0mm blades,
length 9.5cm

S-0690

Stevens Tenotomy Scissors, ribbon type, curved,
blunt tips, 12.0mm blades,
length 9.5cm

S-0700



Strabismus Scissors, ribbon type, straight,
blunt tips, 30.0mm blades,
length 10.5cm

S-0710

Strabismus Scissors, ribbon type, curved,
blunt tips, 29.0mm blades,
length 10.5cm

S-0720

For Ceramic Coating please add "C"



Knapp Strabismus Scissors, straight,
blunt-blunt tips, 31.0mm blades,
length 10.3cm

S-0730

Knapp Strabismus Scissors, curved,
blunt-blunt tips, 30.0mm blades,
length 10.3cm

S-0740



Enucleation Scissors, medium curved,
blunt tips, 43.5mm blades,
length 13.5cm

S-0750

Enucleation Scissors, strong curved,
blunt tips, 40.2mm blades,
length 13.5cm

S-0760

For Ceramic Coating please add "C"



Utility Scissors, straight,
sharp/ blunt tips,
length 13.5cm

S-0770



Metzenbaum Scissors, curved,
blunt tips, 48.0mm blades,
length 14.5cm

S-0780



Spencer Ligature Scissors, straight,
with one notched tip, 21.0mm blades,
length 9.0cm

S-0790

For Ceramic Coating please add "C"



Super-Cut Scissors

Perfectly precise and non-traumatic cutting through any kind of tissue with minimal force requires a SuperCut Scissor. SuperCut Scissors have

- a razor sharp cutting edge on one blade
- and a standard edge combined with a micro-serration on the opposite blade.

The micro-serration prevents tissue from slipping away during the cutting process.

Super-Cut Stevens Tentotomy Scissors, straight,
blunt tips, 14.0mm blades,
length 11.5cm

S-0670SC

Super-Cut Stevens Tentotomy Scissors, curved,
blunt tips, 13.0mm blades,
length 11.5cm

S-0680SC



Super-Cut Metzenbaum Scissors, curved,
blunt tips, 48.0mm, length 14.5cm

S-0780SC

For Ceramic Coating please add "C"



Super-Cut Scissors

Perfectly precise and non-traumatic cutting through any kind of tissue with minimal force requires a SuperCut Scissor. SuperCut Scissors have

- a razor sharp cutting edge on one blade
- and a standard edge combined with a micro-serration on the opposite blade.

The micro-serration prevents tissue from slipping away during the cutting process.

Super-Cut Westcott Tenotomy Scissors, curved,
blunt tips, standard blade, cutting length 14.0mm,
length 11.5cm

S-0220SC

Super-Cut Westcott Tenotomy Scissors, curved,
pointed tips, standard blades, cutting length 14.0mm,
length 11.5cm

S-0250SC

VITRECTOMY INSTRUMENTS



Titanium handle & stainless steel tip

| | | | |
|---|---|-----|-------------|
|  | Vitrectomy Forceps, smooth, with grip straight | 23G | VR-0010-23T |
|  | Vitrectomy Forceps, serrated tip, straight | 23G | VR-0030-23T |
|  | Vitrectomy Forceps, serrated tip, straight | 23G | VR-0050-23T |
|  | Vitrectomy Forceps, with notch, 0.3mm, straight | 23G | VR-0090-23T |
|  | Vitrectomy Forceps, grasping tip | 23G | VR-0110-23T |
|  | Vitrectomy Forceps, serrated, grasping, 60° | 23G | VR-0170-23T |
|  | Endgripping Forceps, straight, extra fine, 0.2mm x 0.5mm | 23G | VR-0205-23T |



Titanium handle & stainless steel tip



Vitrectomy scissors,
straight

23G

VR-0220 -23T



Vitrectomy Scissors,
angled 60°

23G


VR-0240-23T



Vitrectomy Scissors,
with scoop tips, 60°


23G

VR-0250-23T



Vitrectomy Foreign Body Forceps, 20G,
three 1.9mm long arms, slightly angled inwards,
titanium handle & stainless steel tip

VR-0270



Vitrectomy Foreign Body Forceps, 20G,
titanium handle & stainless steel tip

VR-0280



Lens cutting scissors with 4mm stainless steel blades to cut and explant silicone or acrylic lenses through a 1.8mm incision. Both blades have small serrations inside to prevent slippage of lens during cutting and removal procedure.

Lens cutting Scissors, straight blades, 15mm shaft, 18G

VR-1600

Lens cutting Scissors, curved blades, 15mm shaft, 18G

VR-1610



Lens holding forceps for removal of acrylic or silicone IOL's, with multiple grooved jaws, for securely gripping lens during explantation procedure

Grasping Forceps, 20G

VR-1620







Detachable System

Choose the instrument head you require and attach to your handle. Each head will be supplied in a protection cap. This guarantees an easier cleaning & storage.

All instruments are suitable to use with 23G trocar system.

Titanium handle & stainless steel tip

| | | | | |
|---|---|-----|---------------------|----------------|
| | Titanium Handle only, used with tip variations. | | Handle only | VR-5000T |
|  | Micro Gripping Forceps, serrated, 1.5mm tip, 30.0mm shaft. Excellent for peeling and manipulating thicker membrane. | 23G | Complete instrument | VR-0300-23TDT* |
| | Tip for VR-0300-23TDT | | Tip only | VR-0300-23-TIP |
|  | ILM Forceps, 2.0mm tip length, the jaw surface maximizes the gripping ability during ILM peeling. | 23G | Complete instrument | VR-0310-23TDT* |
| | Tip for VR-0310-23TDT | | Tip only | VR-0310-23-TIP |
|  | Straight Scissors, sharp, tip length 2.5mm. Ideal for cutting in the iris plane and peripheral membranes. | 23G | Complete instrument | VR-0400-23TDT* |
| | Tip for VR-0400-23TDT | | Tip only | VR-0400-23-TIP |
|  | Curved Scissors, sharp, tip length 1.25mm. Ideal for cutting parallel to the retina and cutting membranes adherent to the retina surface. | 23G | Complete instrument | VR-0410-23TDT* |
| | Tip for VR-0410-23TDT | | Tip only | VR-0410-23-TIP |

*DT= Detachable

CANNULAS

LACRIMAL INSTRUMENTS

BIMANUAL I/A INSTRUMENTS

DETACHABLE I/A INSTRUMENTS



Kerrison Bone Nibbling Rongeur,
size # 0, 3.0mm

L-0010

Kerrison Bone Nibbling Rongeur,
size # 1, 3.5mm

L-0020

Kerrison Bone Nibbling Rongeur,
size # 2, 4.0mm

L-0030

Kerrison Bone Nibbling Rongeur,
size # 3, 4.5mm

L-0040



Citelli Bone Nibbling Rongeur,
1.5mm wide

L-0050

Citelli Bone Nibbling Rongeur,
2.0mm wide

L-0060

Citelli Bone Nibbling Rongeur,
3.0mm wide

L-0070



Beyer Rongeur, single action, curved

L-0080



Nasal Speculum,
infant

L-0090

Nasal Speculum,
child

L-0100

Nasal Speculum,
adult

L-0110



Castroviejo Lacrimal Dilator, double-ended,
fine & medium taper,
length 14.0cm

L-0120



Wilder Lacrimal Dilator,
size # 1, fine taper L-0130

Wilder Lacrimal Dilator,
size # 2, medium taper L-0140

Wilder Lacrimal Dilator,
size # 3, heavy taper L-0145

Ruedemann Infant Dilator,
length 8.0cm L-0150



Freer Periosteal Elevator,
length 18.0cm L-0160

For Ceramic Coating please add "C"



Partsch Mallet, lead filled,
stainless steel head,
length 18.0cm

L-0170



Joseph Periosteal Elevator,
gently curved tip

L-0180

For Ceramic Coating please add "C"



Partsch Lacrimal Sac Chisel,
length 13.5cm

3.0mm wide L-0200

4.0mm wide L-0210

5.0mm wide L-0220

7.0mm wide L-0230

For Ceramic Coating please add "C"



Partsch Lacrimal Sac Gouge,
length 13.5cm

3.0mm wide

L-0250

4.0mm wide

L-0260

5.0mm wide

L-0270

7.0mm wide

L-0280



West Lacrimal Chisel,
straight

L-0290

For Ceramic Coating please add "C"



Freer Lacrimal Chisel,
curved, 6.0mm,
length 16.0cm

L-0300



Harms Trabeculotomy Probe,
9.0mm long pointed tips, right

L-0310

Harms Trabeculotomy Probe,
9.0mm long pointed tips, left

L-0320

For Ceramic Coating please add "C"



Bowman Lacrimal Probe, Stainless steel, set of 6 L-0330

size 0000 & 000 L-0340

size 00 & 0 L-0350

size 1 & 2 L-0360

size 3 & 4 L-0370

size 5 & 6 L-0380

size 7 & 8 L-0390

For Ceramic Coating please add "C"



Kellnar Retrograde Probe,
right

*L-0400

Kellnar Retrograde Probe,
right

*L-0410



Catalano Lacrimal Intubations Set,
0.7mm x 14.0cm

*L-0420

Silicone tubing

*L-0430

*Not available in Ceramic Coating



Crawford Lacrimal Intubation Set,
0.6mm x 14.0cm

*L-0440

Silicone tubing

*L-0450



Crawford Retrieving Hook,
length 15.0cm

L-0460

For Ceramic Coating please add "C"
* Not available in Ceramic Coating



Pigtail Probe with suture holes,
length 14.8cm

L-0470



Simcoe Cortex Extractor, 0.3mm side port,
23G

L-0480



Simcoe Cortex Extractor, 0.3mm side port,
23G, textured tip

L-0490



Simcoe Cortex Extractor, 0.4mm side port,
21G

L-0500

*Not available in Ceramic Coating



Anis Cortex Aspirating Cannula, straight,
0.4mm, 23G, textured tip

L-0520




Anis Cortex Aspirating Cannula, angled,
0.4mm, 23G, textured tip

L-0530



Vidaurri LASIK Cannula, 8 irrigation ports,
25G, vertical

L-0540



Vidaurri LASIK Cannula, 8 irrigation ports,
25G, horizontal

L-0550



Gills Aspirating Cannula,
5.0mm tip, 25G

L-0560



Gills Aspirating Cannula,
7.0mm tip, 25G

L-0570



Gills Aspirating Cannula,
8.0mm tip, 25G

L-0580



Gills Aspirating Cannula,
10.0mm tip, 25G

L-0590



Lacrimal Cannula, straight, tapered,
23G

L-0600



Lacrimal Cannula, curved, tapered,
23G

L-0610



Lacrimal Cannula, angled, 10.0mm conical

L-0620




Lacrimal Cannula, strong curved, conical

L-0630



McIntyre Lacrimal Cannula, straight, 0.3mm side port, 23G

L-0640



McIntyre Lacrimal Cannula, curved, 0.3mm side port, 23G

L-0650



Anel Lacrimal Cannula, straight,
23G

L-0660



Anel Lacrimal Cannula, curved,
23G

L-0670



Bangerter Lacrimal Probe Cannula,
40.0mm, 21G

L-0680



Bangerter Lacrimal Probe Cannula,
60.0mm, 21G

L-0690



Rycroft Air Injection Cannula, angled, 5.0mm, 27G

L-0700



Rycroft Air Injection Cannula, angled, 5.0mm, 30G

L-0710



Air Injection Cannula, angled, 7.0mm, 27G

L-0720



Air Injection Cannula, angled, 7.0mm, 30G

L-0730




Rainin Injection Cannula, angled, 7.0mm, 27G

L-0740



Hydrodissection Cannula, angled, 8.0mm, flat, front opening 25G

L-0750



Seeley Hydrodissection Cannula, angled, 12.0mm, flat, front opening 26G

L-0760



Pearce Hydrodissection Cannula, J-shaped, straight, front opening 25G & 22.0mm tube

L-0780




Pearce Hydrodissection Cannula, J-shaped, 8.0mm, angled, front opening 25G, 22.0mm tube

L-0790



Frankfurt Model Hydrodissections Cannula, 27G

L-0800



Sauter Hydrodissection Cannula, vertically & flat tip,
front opening 27G

L-0810



McIntyre-Binkhorst Cannula, right, blunt tip,
26G

L-0820




McIntyre-Binkhorst Cannula, straight, blunt tip,
26G

L-0830



McIntyre-Binkhorst Cannula, left, blunt tip,
26G

L-0840



Bishop-Harmon Anterior Chamber Irrigator, complete,
(Cannula, adaptor and silicone bulb),
19G

L-0850

silicone bulb & adaptor

L-0860

adaptor

L-0870

Cannula,
19G

L-0880

Cannula,
21G

L-0890



| | | |
|--|--|--------|
| | McIntyre Anterior Chamber Irrigating Cannula, straight, 15.0mm long tube, blunt tip, 30G | L-0900 |
| | McIntyre Anterior Chamber Irrigating Cannula, straight, 15.0mm long tube, blunt tip, 26G | L-0910 |
| | McIntyre Anterior Chamber Irrigating Cannula, straight, 15.0mm long tube, blunt tip, 23G | L-0920 |
| | McIntyre Anterior Chamber Irrigating Cannula, angled, 15.0mm long tube, blunt tip, 30G | L-0930 |
| | McIntyre Anterior Chamber Irrigating Cannula, angled, 15.0mm long tube, blunt tip, 26G | L-0940 |
| | McIntyre Anterior Chamber Irrigating Cannula, angled, 15.0mm long tube, blunt tip, 23G | L-0950 |



Knolle Anterior Chamber Irrigating Cannula, angled, 2.0mm, front opening, 27G

L-0960



Knolle Anterior Chamber Irrigating Cannula, angled, 2.0mm, front opening, 23G

L-0970



Heslin Anterior Chamber Maintaining Cannula, 25G

L-0980



Heslin Anterior Chamber Maintaining Cannula, 23G

L-0990




Alfonso-McIntyre
Irrigating Nucleus Spoon, single 0.3mm,
with front opening port 2.5mm wide

L-1000




Alfonso-McIntyre
Irrigating Nucleus Spoon, single 0.3mm,
with front opening port 4.0mm wide

L-1001



Jensen Capsule Polisher, front opening, sandblasted,
23G

L-1010



Jensen Capsule Polisher, front opening, sandblasted,
25G

L-1020



Jensen Capsule Polisher, front opening, polished,,
25G

L-1021



Jensen Posterior Capsule Polisher, olive tip, front opening,
diamond dusted, 27G

L-1030



Jensen Posterior Capsule Polisher, olive tip,
front opening sandblasted, 23G

L-1040



Gills Welsh Olive Tip Cannula , 25G

L-1041



Kratz Capsule Polisher,
22G

L-1050



Rubman-Katzin Capsule Polisher,
22G

L-1060



Capsule Polishing Cannula with bulbous tip,
18G

L-1070



MICS Capsule Polishing Cannula, 23G,
textured disc-shaped tip
with 0.5mm posterior facing port

L-1075



Knolle-Pearce Irrigating Vectis Cannula,
with three front opening ports,
23G Loop

L-1080



Sheets Irrigating Vectis Cannula,
with three front opening ports, serrated,
20G Loop

L-1090



Cranz Eye Douche

L-1100



Atkinson Retrobulbar Needle,
23G & 35.0mm long

L-1110



Fukasaku Anaesthesia Cannula, gently curved,
24G

L-1120



Buratto LASIK Cannula, on front and two side ports,
25G

L-1130



Buratto LASIK Cannula, on front and two side ports,
23G

L-1140



Viscoelastic Aspirating Cannula, angled, 9.0mm, 23G

L-1150



Viscoelastic Aspirating Cannula, angled, 9.0mm, 21G

L-1160



Viscoelastic Aspirating Cannula, angled, 10.0mm, 23G

L-1170



Viscoelastic Aspirating Cannula, angled, 10.0mm, 21G

L-1180



Viscoelastic Aspirating Cannula, angled, 10.0mm, 22G

L-1200



Viscoelastic Aspirating Cannula, angled, 10.0mm, 19G

L-1210



Viscoelastic Aspirating Cannula, angled, 11.0mm, flattened tip 23G

L-1220



Viscoelastic Aspirating Cannula, angled, 11.0mm, flattened tip 22G

L-1230



Viscoelastic Aspirating Cannula, angled, 11.0mm, flattened tip 21G

L-1240



Kratz Cystotome Cannula, 10.0mm, angled tube, with sharp points, bent downward, 22G

L-1250




McIntyre Cystotome Cannula, 10.0mm, angled tube, 23G

L-1260




Kelman Cannula, 30.0mm, straight tube,
with sharp point & sharp cutting edge,
22G

L-1270



Randolph Cyclodialysis Cannula, angled, 13.0mm
19G

L-1280



Randolph Cyclodialysis Cannula, angled, 13.0mm
21G

L-1290



Lippmann Micro Cystotome, with silicone protection,
22G

L-1300



DALK Cannula, 27G,
flattened tip, 0.2mm posterior port

L-1304



Detachable I/A Instrument complete

Instruments complete = titanium handle + sealing ring + stainless steel tip



Irrigation Handpiece,
2.0 x 0.5mm, 21G, polished

L-1311



Irrigation Handpiece,
2.0 x 0.3mm, 23G, rough

L-1321



Aspiration Handpiece,
0.3mm, 21G, polished

L-1331



Aspiration Handpiece,
0.5mm, 21G, rough

L-1341



Brauweiler Aspiration Handpiece,
0.5mm, 21G, rough

L-1351



Brauweiler Irrigation Handpiece,
0.5mm, 21G, rough

L-1361



Brauweiler Aspirating Handpiece,
0.5mm, 21G, rough

L-1371



Koch Aspirating Handpiece,
0.5mm, 21G, rough

L-1381



Tip only, for detachable I/A system



Tip for Irrigation Handpiece L-1311
2.0 x 0.5mm, 21G, polished

L-1311-TIP



Tip for Irrigation Handpiece, L-1321,
2.0 x 0.3mm, 23G, rough

L-1321-TIP




Tip for Aspiration Handpiece, L-1331,
0.3mm, 21G, polished

L-1331-TIP



Tip for Aspiration Handpiece, L-1341,
0.5mm, 21G, rough

L-1341-TIP



Tip for Aspirating Handpiece, L-1351,
0.5mm, 21G, rough

L-1351-TIP



Tip for Brauweiler Irrigation Handpiece, L-1361,
0.5mm, 21G, rough

L-1361-TIP



Tip for Brauweiler Irrigation Handpiece, L-1371,
0.5mm, 21G, rough

L-1371-TIP



Tip for Aspirating Handpiece, L-1381,
0.5mm, 21G, rough

L-1381-TIP



Titanium handles only,
for detachable I/A system

Aspiration Handpiece only,
(includes 2x sealing rings)

L-1300-A



Irrigation Handpiece only,
(includes 2x sealing rings)

L-1300-I



Gills Aspirating Cannula,
23G, with silicone tubing

L-1400




Simcoe Irrigating/ Aspirating Cannula,
23G

L-1410



Simcoe Irrigating/ Aspirating Cannula,
reversed, 23G

L-1420



Gess Irrigating-Aspirating Cannula, curved right,
with end opening aspiration port,
23G

L-1430



Gess Irrigating-Aspirating Cannula, curved left,
with end opening aspiration port,
23G

L-1440



Infusions Plug, 20G,
45° bevelled tip, 2.5mm long tip

L-1450



Infusions Plug, 20G,
45° bevelled tip, 3.5mm long tip

L-1460



Infusions Plug, 20G,
45° bevelled tip, 6.0mm long tip

L-1470



Infusions Plug, 20G,
45° bevelled tip, 7.0mm long tip

L-1480



Monomanual Phaco Safe-Lock I/A Handpiece,
with male/female luer slip fittings, compatible with
all I/A and phaco systems

L-1700



Safe - Lock I/A tip, 16G, straight,
for standard aspiration procedures with 0.35mm aspirating port

L-1710



Safe - Lock I/A tip 16G, straight,
for removal of viscoelastic fluids with 0.5mm aspirating port

L-1720



Safe - Lock I/A tip, 16G, angled 45°,
for peripheral cortex removal

L-1730



Safe - Lock I/A tip, 16G, angled 90°,
for peripheral and superior cortex removal

L-1740



Safe - Lock I/A tip, 16G, curved,
for initial cortex removal

L-1750



O-Ring Replacement Set,
for I/A tips, 5 each rings 2x1, 3x1
and 5x1mm, (pkg/15)

L-1760



Safe - Lock I/A SI-Sleeve tip, 16G, straight,
without soft style si-sleeve

L-1770



Safe - Lock I/A SI-Sleeve tip, 20G, straight,
without soft style si-sleeve

L-1780



Safe - Lock I/A SI-Sleeve tip, 45°, 16G, straight,
without soft style si-sleeve

L-1790



Silicone Irrigation Sleeve,
non-sterile (pack of 5).
used for L-1770 - L-1790

L-1800



Luer-Lock Adaptor,
2x Luer male

L-1820



Luer-Lock Adaptor,
2x female connectors

L-1830

MARKING INSTRUMENTS



Castroviejo Caliper, 20.0mm,
straight

*M-0010

Castroviejo Caliper, 20.0mm,
curved

*M-0020

Castroviejo Caliper, 20.0mm,
with 0.5mm scale, straight

*M-0030



Jameson Caliper, adjustable,
0-80mm scale in 1.0mm increments

*M-0040

* Not available in Ceramic Coating



Osher Incision Caliper,
1.0 to 10.0mm in increments of 1.0mm

*M-0050



Ruler, stainless steel

*M-0060



Braunstein Measuring Instrument,
3.0mm/ 3.5mm

*M-0070

Braunstein Measuring Instrument,
3.0mm/ 4.0mm

*M-0080

Braunstein Measuring Instrument,
3.5mm/ 4.0mm

*M-0090

* Not available in Ceramic Coating



Schiøtz Tonometer, autoclavable, straight scale,
incl. certificate

MIS-0011

Schiøtz Tonometer, autoclavable, curved scale,
with metrological controll report "add M to the article number"

MIS-0020
MIS-0020M



Honan Balloon,
Intraocular Pressure Reducer, complete

MIS-0030

MISCELLANEOUS



Carter Sphere Introducer,
length 12.5cm

MIS-0010



Steri Oil Pen for lubrication of instruments with articulations,
12ml, 0,4 fl. oz., 7.0mm tip

MIS-0770



Scleral Eye Shield (pair),
large, stainless steel

MIS-0142

Scleral Eye Shield (pair),
small, stainless steel

MIS-0143



Metal Bole, Ø 58.0mm,
height 40.0mm

MIS-0150



Lens Injector for use with Type A, B & C cartridges,
single-hand IOL Implantation, titanium,
length 20.0cm

I-0010

Lens Injector for use with Type D cartridges,
single-hand IOL Implantation, titanium,
length 20.0cm

I-0020




Lens Injector for use with Type A, B & C cartridges,
screw mechanism, IOL Implantation, titanium,
length 16.0cm

I-0030

Lens Injector for use with Type D cartridges,
screw mechanism, IOL Implantation, titanium
length 16.0cm

I-0040



Sterilizing Case, autoclavable plastic,
base, lid & silicone mat,
outside dimension 38x68x19mm

MIS-0160

Sterilizing Case, autoclavable plastic,
base deep, lid & silicone mat,
outside dimension 38x68x25mm

MIS-0170

Sterilizing Case, autoclavable plastic,
base deep, lid & silicone mat,
outside dimension 38x68x32mm

MIS-0180

Sterilizing Case, autoclavable plastic,
base & silicone mat,
outside dimension 200x110x25mm

MIS-0190

Sterilizing Case, autoclavable plastic,
base, silicone mat,
outside dimension 265x160x25mm


MIS-0200

Sterilizing Case, autoclavable plastic,
base, double & silicone mat,
outside dimension 265x160x50mm

MIS-0210

Sterilizing Case, autoclavable plastic,
base lid & silicone mat,
outside dimension 64x150x19mm

MIS-0220



Sterilization Case, autoclavable plastic,
medium base, lid & silicone mat ,
outside dimension 64 x150x32mm

MIS-0230

Sterilization Case, autoclavable plastic,
medium base, lid & silicone mat,
outside dimension 64x190x19mm

MIS-0240

Sterilization Case, autoclavable plastic,
medium base, lid & silicone mat,
outside dimension 102x62x19mm

MIS-0250

Sterilization Case, autoclavable plastic,
medium base, lid & silicone mat,
outside dimension 102x190x38mm

MIS-0260

Sterilization Case, autoclavable plastic,
medium base, lid & silicone mat,
outside dimension 152x254x38mm


MIS-0270

Phaco Tray, autoclavable plastic,
base, lid & 2 bars / 1 Slot,
outside dimension 100x19x38mm

MIS-0280

Phaco Tray, autoclavable plastic,
base, lid & 2 bars / 1 Slot
outside dimension 150x25x38mm

MIS-0290



Phaco Trays, Base, autoclavable plastic,
lid, two silicone bars / 3 Slot,
outside dimension 150x25x38mm

MIS-0300

Knife Tray, autoclavable plastic,
1 slot-base, lid, two silicone bars,
outside dimension 64x152x19mm

MIS-0310

Knife Tray, autoclavable plastic,
2 slot-base, lid, two silicone bars,
outside dimension 64x152x19mm

MIS-0320

Knife Tray, autoclavable plastic,
3 slot-base, lid, two silicone bars,
outside dimension 64x152x19mm

MIS-0330

Knife Tray, autoclavable plastic,
5 slot-base, lid, 2 silicone bars,
outside dimension 102x162x19mm

MIS-0340



Sterilisation Case KSI 1, stainless steel, with lid & bottom perforated, 1 Silicone mat, textil filter, outside dimension 220x165x55mm MIS-0350

Sterilisation Case KSI 1, stainless steel, with lid & bottom perforated, 1 Silicone mat, textil filter, outside dimension 300x220x55mm MIS-0360

Sterilisation Case KSI 1, stainless steel, with lid & bottom perforated, 1 Silicone mat, textil filter, outside dimension 285x280x55mm MIS-0370

Sterilisation Case KSI 1, stainless steel, with lid & bottom perforated, 1 Silicone mat, textil filter, outside dimension 400x290x55mm MIS-0380

Double Sterilisation Case KSI 1, stainless steel, with lid & bottom perforated, 2 Silicone mat, textil filter, outside dimension 220x165x75mm MIS-0390

Double Sterilisation Case KSI 1, stainless steel, with lid & bottom perforated, 2 Silicone mat, textil filter, outside dimension 300x220x75mm MIS-0400

Double Sterilisation Case KSI 1, stainless steel, with lid & bottom perforated, 2 Silicone mat, textil filter, outside dimension 285x280x75mm MIS-0410

Sterilisation Case KSI 1, stainless steel, with lid & bottom perforated, double, 2 Silicone mat & filter, outside dimension 400x290x75mm MIS-0420



The aluminium sterilisation containers are specially designed to accommodate 1/1, 1/2 & 1/3 Trays assuring complete sterilisation, as well as safe in-house transportation and storage.

Sterilisation Case aluminium, complete with quick-release, MIS-0430
with positive action spring mechanism on the lid, bottom perforated,
outside dimension 310x190x40

Sterilisation Case aluminium, complete with quick-release, MIS-0440
with positive action spring mechanism on the lid, bottom perforated,
outside dimension 310x190x65

Sterilisation Case aluminium, complete with quick-release, MIS-0450
with positive action spring mechanism on the lid, bottom perforated,
outside dimension 310x190x80



Paper Filters,with chemical indicator,
disposable (pack of 100) for single use.
For MIS-0430; MIS-0440; MIS-0450.

MIS-0460

Textile Filters, reusable (set of 2) for frequent use.
For MIS-0430; MIS-0440; MIS-0450.

MIS-0470

PTFE Teflon Filters, long-life (set of 2) for permanent use.
For MIS-0430; MIS-0440; MIS-0450.

MIS-0480

All filters are suitable for all steam sterilisation methods.




Identification Plates ALUMINIUM (50.0 x 18.0mm),
can be engraved to your specification;
max. 16 letters

MIS-0490



Indicator lable,
disposable (pack of 100)

MIS-0500



Plastic Safety Seals,
disposable (pack of 100)

MIS-0510



Wash Tray stainless steel, complete,
with 2 silicone instrument supports
and 1 instrument fixation bar,
1/1 size: 275x178x24mm
for 16 hand instruments

MIS-0520

Wash Tray stainless steel, complete,
with 2 silicone instrument supports
and 1 instrument fixation bar,
1/2 size: 178x140x24mm
for 8 hand instruments

MIS-0530


Wash Tray stainless steel, 24 mm height,
complete with 2 silicone instrument supports
and 1 instrument fixation bar,
1/3 size: 178x90x24mm
for 5 hand instruments

MIS-0540



Silicone Instrument Support, 265.0mm long rubber strip,
with 16-slots to fit 1/1 size tray MIS-0520.
Can be cut to fit 1/2 MIS-0530 or 1/3 MIS-0540 size trays, too.

MIS-0550



Silicone Finger Mats, extra soft, light blue material,
with narrow spaced fine nubbins,
275x178x24mm


MIS-0560

Silicone Finger Mats, extra soft, light blue material,
with narrow spaced fine nubbins,
178x140x24mm

MIS-0570

Silicone Finger Mats, extra soft, light blue material,
with narrow spaced fine nubbins,
178x90x24mm

MIS-0580




Sterilisation tray with lid, stainless steel, with silicone mats,
storage of up to 16 instruments,
outer dimensions: 261x219x26mm

MIS-0614



Double-level sterilization tray, stainless steel, with silicone mats,
storage of up to 50 instruments,
outer dimensions: 325x265x 42mm

MIS-0615



Sterilization tray with lid, stainless steel, electropolished,
270x175x30 mm; 4 lower and 2 upper instrument supports,
for use in autoclave

MIS-0611



Upper Instrument Support element,
for 10 Instruments

MIS-0612



Lower instrument support element,
for 10 Instruments,

MIS-0613

CONTENT

A

| | |
|------|----|
| A-1 | 51 |
| A-2 | 51 |
| A-3 | 51 |
| A-3 | 55 |
| A-4 | 51 |
| A-5 | 51 |
| A-6 | 52 |
| A-7 | 52 |
| A-8 | 52 |
| A-9 | 52 |
| A-10 | 52 |
| A-11 | 52 |
| A-12 | 53 |
| A-13 | 53 |
| A-14 | 53 |
| A-15 | 53 |
| A-16 | 53 |
| A-17 | 53 |
| A-18 | 53 |
| A-19 | 53 |
| A-20 | 54 |
| A-21 | 54 |
| A-22 | 54 |
| A-23 | 54 |
| A-24 | 54 |
| A-25 | 54 |
| A-26 | 55 |
| A-27 | 55 |
| A-28 | 55 |
| A-29 | 55 |
| A-30 | 55 |
| A-31 | 55 |
| A-32 | 55 |
| A-33 | 55 |
| A-34 | 55 |
| A-35 | 55 |
| A-36 | 55 |
| A-39 | 52 |
| A-45 | 54 |

C

| | |
|--------|----|
| C-0010 | 35 |
|--------|----|

| | |
|--------|----|
| C-0020 | 35 |
| C-0030 | 35 |
| C-0040 | 35 |
| C-0050 | 36 |
| C-0060 | 36 |
| C-0070 | 36 |
| C-0080 | 36 |
| C-0090 | 36 |
| C-0100 | 36 |
| C-0110 | 36 |
| C-0120 | 37 |
| C-0130 | 37 |
| C-0140 | 37 |
| C-0150 | 37 |
| C-0160 | 38 |
| C-0170 | 38 |
| C-0180 | 38 |
| C-0190 | 38 |
| C-0200 | 38 |
| C-0210 | 38 |
| C-0220 | 38 |
| C-0230 | 39 |
| C-0240 | 39 |
| C-0250 | 39 |
| C-0260 | 39 |
| C-0280 | 40 |
| C-0290 | 40 |
| C-0300 | 40 |
| C-0310 | 40 |
| C-0320 | 40 |
| C-0330 | 41 |
| C-0340 | 41 |
| C-0350 | 41 |
| C-0360 | 42 |
| C-0361 | 42 |
| C-0362 | 42 |
| C-0370 | 43 |
| C-0380 | 43 |
| C-0390 | 43 |
| C-0400 | 43 |
| C-0410 | 43 |
| C-0420 | 43 |
| C-0430 | 44 |
| C-0440 | 44 |
| C-0450 | 44 |
| C-0460 | 44 |
| C-0470 | 44 |
| C-0480 | 45 |
| C-0490 | 45 |
| C-0500 | 45 |

| | |
|---------------|----|
| C-0510 | 45 |
| C-0520 | 45 |
| C-0530 | 45 |
| C-0550 | 46 |
| C-0560 + size | 46 |
| C-0570 + size | 46 |
| C-0670 | 46 |
| C-0680 | 47 |
| C-0690 | 47 |
| C-0700 | 47 |
| C-0710 | 47 |
| C-0711 | 47 |
| C-0712 | 47 |
| C-0714 | 47 |
| C-0720 | 48 |
| C-0730 | 48 |
| C-0740 | 48 |
| C-0750 | 49 |
| C-0760 | 49 |

ES

| | |
|---------|----|
| ES-0010 | 11 |
| ES-0020 | 11 |
| ES-0030 | 11 |
| ES-0050 | 12 |
| ES-0060 | 12 |
| ES-0070 | 12 |
| ES-0080 | 12 |
| ES-0090 | 13 |
| ES-0100 | 13 |
| ES-0122 | 13 |
| ES-0130 | 14 |
| ES-0140 | 14 |
| ES-0150 | 14 |
| ES-0160 | 14 |
| ES-0170 | 15 |
| ES-0180 | 15 |
| ES-0190 | 15 |
| ES-0200 | 15 |
| ES-0210 | 15 |
| ES-0220 | 16 |
| ES-0230 | 16 |
| ES-0240 | 16 |
| ES-0250 | 16 |
| ES-0260 | 17 |
| ES-0270 | 17 |
| ES-0280 | 17 |

| | | | | | |
|---------|----|---------|----|---------|-----|
| ES-0290 | 17 | ES-0820 | 30 | F-0570T | 98 |
| ES-0300 | 18 | ES-0830 | 30 | F-0580 | 98 |
| ES-0310 | 18 | ES-0840 | 30 | F-0580T | 98 |
| ES-0311 | 11 | ES-0850 | 31 | F-0590 | 99 |
| ES-0320 | 18 | ES-0860 | 31 | F-0590T | 99 |
| ES-0330 | 18 | ES-0870 | 31 | F-0600 | 99 |
| ES-0340 | 18 | ES-0880 | 31 | F-0600T | 99 |
| ES-0341 | 18 | ES-0890 | 31 | F-0610 | 99 |
| ES-0350 | 19 | ES-0900 | 31 | F-0610T | 99 |
| ES-0360 | 19 | ES-0910 | 31 | F-0620 | 99 |
| ES-0370 | 19 | ES-0920 | 32 | F-0620T | 99 |
| ES-0371 | 19 | ES-0930 | 32 | F-0650 | 99 |
| ES-0380 | 19 | ES-0940 | 32 | F-0650T | 99 |
| ES-0390 | 20 | ES-0950 | 32 | F-0680 | 99 |
| ES-0400 | 20 | ES-0960 | 33 | F-0680T | 99 |
| ES-0410 | 20 | ES-1050 | 33 | F-0740 | 100 |
| ES-0420 | 20 | ES-1060 | 33 | F-0740T | 100 |
| ES-0430 | 21 | | | F-0770 | 100 |
| ES-0440 | 21 | | | F-0770T | 100 |
| ES-0450 | 21 | | | F-0790 | 100 |
| ES-0460 | 21 | | | F-0790T | 100 |
| ES-0500 | 22 | | | F-0810 | 100 |
| ES-0510 | 22 | | | F-0810T | 100 |
| ES-0520 | 22 | | | F-0820 | 100 |
| ES-0530 | 22 | F-0010 | 94 | F-0820T | 100 |
| ES-0540 | 23 | F-0010T | 94 | F-0830 | 101 |
| ES-0550 | 23 | F-0020 | 94 | F-0840 | 101 |
| ES-0560 | 23 | F-0020T | 94 | F-0850 | 101 |
| ES-0570 | 23 | F-0030 | 94 | F-0860 | 101 |
| ES-0580 | 23 | F-0030T | 94 | F-0870 | 101 |
| ES-0590 | 24 | F-0040 | 94 | F-0880 | 101 |
| ES-0600 | 24 | F-0040T | 94 | F-0890 | 101 |
| ES-0610 | 24 | F-0050 | 94 | F-0900 | 102 |
| ES-0620 | 24 | F-0050T | 94 | F-0930 | 102 |
| ES-0630 | 24 | F-0060 | 94 | F-0940 | 102 |
| ES-0640 | 25 | F-0060T | 94 | F-0950 | 103 |
| ES-0650 | 25 | F-0070 | 95 | F-0960 | 103 |
| ES-0660 | 26 | F-0080 | 95 | F-0970 | 103 |
| ES-0670 | 26 | F-0100 | 95 | F-0980 | 103 |
| ES-0680 | 26 | F-0120 | 95 | F-0990 | 103 |
| ES-0690 | 26 | F-0140 | 95 | F-1000 | 103 |
| ES-0700 | 27 | F-0160 | 95 | F-1003 | 115 |
| ES-0710 | 27 | F-0230 | 96 | F-1004 | 115 |
| ES-0715 | 27 | F-0250 | 96 | F-1010 | 103 |
| ES-0720 | 28 | F-0260 | 96 | F-1020 | 104 |
| ES-0725 | 28 | F-0270 | 96 | F-1030 | 104 |
| ES-0730 | 28 | F-0280 | 96 | F-1040 | 104 |
| ES-0740 | 29 | F-0410 | 97 | F-1050 | 104 |
| ES-0750 | 29 | F-0440 | 97 | F-1060 | 104 |
| ES-0760 | 29 | F-0450 | 97 | F-1070 | 104 |
| ES-0770 | 29 | F-0490 | 97 | F-1090 | 104 |
| ES-0770 | 29 | F-0510 | 98 | F-1120 | 105 |
| ES-0780 | 29 | F-0520 | 98 | F-1130 | 105 |
| ES-0790 | 30 | F-0560 | 98 | F-1140 | 105 |
| ES-0800 | 30 | F-0560T | 98 | F-1150 | 105 |
| ES-0810 | 30 | F-0570 | 98 | | |

F

| | | | | | |
|---------|-----|---------|-----|---------|-----|
| F-1160 | 105 | F-1850 | 113 | F-2970T | 123 |
| F-1190 | 106 | F-1860 | 113 | F-2980 | 123 |
| F-1200 | 106 | F-1870 | 113 | F-2981 | 117 |
| F-1210 | 106 | F-1910 | 113 | F-2982 | 117 |
| F-1220 | 106 | F-1910T | 113 | F-2983 | 117 |
| F-1240 | 106 | F-1920 | 113 | F-2990 | 123 |
| F-1250 | 106 | F-1920T | 113 | F-2990T | 123 |
| F-1260 | 106 | F-1930 | 113 | F-3000 | 123 |
| F-1290 | 106 | F-1930T | 113 | F-3000 | 123 |
| F-1300 | 107 | F-1940 | 113 | F-3100 | 124 |
| F-1310 | 107 | F-1940T | 113 | F-3110 | 124 |
| F-1320 | 107 | F-1960 | 115 | F-3120 | 124 |
| F-1380 | 107 | F-1970 | 115 | F-3130 | 124 |
| F-1400 | 107 | F-1971 | 114 | F-3140 | 124 |
| F-1410 | 107 | F-1975 | 114 | F-3150 | 124 |
| F-1420 | 107 | F-1976 | 115 | F-3151 | 116 |
| F-1510 | 108 | F-1980 | 135 | F-3160 | 124 |
| F-1520 | 108 | F-1990 | 114 | F-3170 | 124 |
| F-1530 | 108 | F-2000 | 114 | F-3190 | 116 |
| F-1540 | 108 | F-2421 | 113 | F-3191 | 116 |
| F-1550 | 108 | F-2460 | 118 | F-3200 | 116 |
| F-1560 | 109 | F-2470 | 118 | F-3210 | 125 |
| F-1560T | 109 | F-2480 | 118 | F-3220 | 125 |
| F-1570 | 109 | F-2490 | 118 | F-3250 | 125 |
| F-1570T | 109 | F-2500 | 118 | F-3260 | 125 |
| F-1580 | 109 | F-2510 | 118 | F-3270 | 125 |
| F-1580T | 109 | F-2520 | 118 | F-3280 | 125 |
| F-1590 | 109 | F-2540 | 119 | F-3290 | 125 |
| F-1590T | 109 | F-2550 | 119 | F-3300 | 125 |
| F-1610 | 109 | F-2560 | 119 | F-3310 | 126 |
| F-1620 | 109 | F-2570 | 119 | F-3320 | 126 |
| F-1630 | 110 | F-2580 | 119 | F-3330 | 126 |
| F-1670 | 110 | F-2590 | 119 | F-3331 | 126 |
| F-1680 | 110 | F-2800 | 120 | F-3340 | 126 |
| F-1680T | 110 | F-2810 | 120 | F-3350 | 126 |
| F-1690 | 110 | F-2820 | 120 | F-3360 | 126 |
| F-1700 | 110 | F-2830 | 121 | F-3370 | 126 |
| F-1710 | 110 | F-2840 | 121 | F-3390 | 127 |
| F-1710T | 110 | F-2850 | 120 | F-3400 | 127 |
| F-1720 | 111 | F-2860 | 120 | F-3410 | 127 |
| F-1720T | 111 | F-2880 | 121 | F-3420 | 127 |
| F-1730 | 111 | F-2890 | 122 | F-3530 | 128 |
| F-1740 | 111 | F-2890T | 122 | F-3540 | 128 |
| F-1750 | 111 | F-2900 | 122 | F-3550 | 128 |
| F-1760 | 111 | F-2900T | 122 | F-3560 | 128 |
| F-1770 | 111 | F-2910 | 122 | F-3570 | 128 |
| F-1780 | 112 | F-2910T | 122 | F-3590 | 128 |
| F-1790 | 112 | F-2920 | 122 | F-3600 | 128 |
| F-1800 | 112 | F-2920T | 122 | F-3610 | 129 |
| F-1810 | 112 | F-2930 | 122 | F-3620 | 129 |
| F-1820 | 112 | F-2940 | 122 | F-3630 | 129 |
| F-1820T | 112 | F-2940T | 122 | F-3640 | 129 |
| F-1830 | 112 | F-2960 | 122 | F-3650 | 129 |
| F-1830T | 112 | F-2960T | 122 | F-3660 | 129 |
| F-1840 | 112 | F-2970 | 123 | F-3670 | 129 |

| | | | | | |
|--------|-----|---------|-----|-----------|----|
| F-3690 | 130 | F-4290 | 140 | F-4800TC | 93 |
| F-3700 | 130 | F-4300 | 140 | F-4805 | 91 |
| F-3710 | 130 | F-4310 | 141 | F-4805TC | 93 |
| F-3720 | 130 | F-4320 | 141 | F-4810 | 91 |
| F-3730 | 130 | F-4330 | 141 | F-4820 | 91 |
| F-3750 | 130 | F-4340 | 141 | F-4840 | 91 |
| F-3760 | 130 | F-4350 | 142 | F-4850 | 91 |
| F-3780 | 131 | F-4360 | 142 | F-4860 | 91 |
| F-3810 | 131 | F-4370 | 142 | F-4870 | 91 |
| F-3820 | 131 | F-4380 | 142 | F-4880 | 92 |
| F-3830 | 131 | F-4390 | 143 | F-4890 | 92 |
| F-3840 | 131 | F-4395 | 142 | F-4900 | 92 |
| F-3860 | 131 | F-4400 | 143 | F-4910 | 92 |
| F-3870 | 132 | F-4405 | 143 | F-4930 | 92 |
| F-3880 | 132 | F-4410 | 143 | F-4940 | 92 |
| F-3890 | 132 | F-4415 | 143 | F-4950 | 92 |
| F-3900 | 132 | F-4420 | 144 | F-4950TC | 93 |
| F-3910 | 132 | F-4430 | 144 | | |
| F-3920 | 132 | F-4440 | 144 | | |
| F-3930 | 132 | F-4450 | 145 | | |
| F-3940 | 133 | F-4460 | 145 | | |
| F-3950 | 133 | F-4470 | 145 | | |
| F-3960 | 133 | F-4480 | 145 | | |
| F-3970 | 133 | F-4485 | 145 | | |
| F-3980 | 133 | F-4487 | 145 | | |
| F-4000 | 133 | F-4490 | 146 | | |
| F-4010 | 133 | F-4500 | 146 | | |
| F-4020 | 133 | F-4510 | 146 | | |
| F-4030 | 133 | F-4520 | 146 | | |
| F-4040 | 134 | F-4530 | 146 | | |
| F-4050 | 134 | F-4540 | 146 | | |
| F-4060 | 134 | F-4550 | 147 | | |
| F-4070 | 134 | F-4560 | 147 | | |
| F-4080 | 135 | F-4570 | 147 | | |
| F-4090 | 135 | F-4575T | 147 | | |
| F-4100 | 135 | F-4580 | 148 | | |
| F-4110 | 135 | F-4590 | 148 | | |
| F-4120 | 136 | F-4600 | 148 | | |
| F-4130 | 136 | F-4610 | 148 | | |
| F-4140 | 136 | F-4620 | 148 | | |
| F-4150 | 136 | F-4630 | 148 | | |
| F-4160 | 136 | F-4640 | 149 | | |
| F-4170 | 137 | F-4650 | 149 | | |
| F-4180 | 137 | F-4660 | 149 | | |
| F-4190 | 138 | F-4670 | 149 | | |
| F-4200 | 137 | F-4680 | 149 | | |
| F-4210 | 138 | F-4690 | 150 | | |
| F-4220 | 138 | F-4700 | 150 | | |
| F-4230 | 138 | F-4710 | 150 | | |
| F-4231 | 139 | F-4720 | 150 | | |
| F-4240 | 139 | F-4730 | 150 | | |
| F-4250 | 139 | F-4730 | 151 | | |
| F-4260 | 139 | F-4740 | 151 | | |
| F-4270 | 139 | F-4750 | 151 | | |
| F-4280 | 140 | F-4800 | 91 | | |
| | | | | HS | |
| | | | | HS-0000TS | 51 |
| | | | | HS-0001TS | 51 |
| | | | | HS-0002TS | 55 |
| | | | | HS-0010TS | 51 |
| | | | | HS-0020TS | 51 |
| | | | | HS-0030TS | 51 |
| | | | | HS-0040TS | 51 |
| | | | | HS-0050TS | 52 |
| | | | | HS-0060TS | 52 |
| | | | | HS-0070TS | 52 |
| | | | | HS-0080TS | 52 |
| | | | | HS-0090TS | 52 |
| | | | | HS-0091TS | 52 |
| | | | | HS-0110TS | 52 |
| | | | | HS-0120TS | 53 |
| | | | | HS-0130TS | 53 |
| | | | | HS-0140TS | 53 |
| | | | | HS-0150TS | 53 |
| | | | | HS-0160TS | 53 |
| | | | | HS-0170TS | 53 |
| | | | | HS-0180TS | 53 |
| | | | | HS-0190TS | 53 |
| | | | | HS-0200TS | 54 |
| | | | | HS-0210TS | 54 |
| | | | | HS-0220TS | 54 |
| | | | | HS-0230TS | 54 |
| | | | | HS-0240TS | 54 |
| | | | | HS-0250TS | 54 |
| | | | | HS-0256TS | 54 |
| | | | | HS-0260TS | 55 |
| | | | | HS-0270TS | 55 |
| | | | | HS-0280TS | 55 |
| | | | | HS-0310TS | 55 |

| | | | | | |
|-----------|----|---------|----|----------|----|
| HS-0320TS | 55 | HS-0820 | 63 | HS-1490 | 73 |
| HS-0330TS | 55 | HS-0830 | 64 | HS-1500 | 74 |
| HS-0340 | 56 | HS-0840 | 64 | HS-1510 | 74 |
| HS-0350 | 56 | HS-0850 | 64 | HS-1520 | 74 |
| HS-0360 | 56 | HS-0860 | 64 | HS-1530 | 74 |
| HS-0370 | 56 | HS-0870 | 64 | HS-1540 | 74 |
| HS-0380 | 57 | HS-0880 | 64 | HS-1545 | 74 |
| HS-0390 | 57 | HS-0890 | 64 | HS-1550 | 74 |
| HS-0400 | 57 | HS-0895 | 64 | HS-1560 | 74 |
| HS-0410 | 57 | HS-0900 | 65 | HS-1570 | 75 |
| HS-0420 | 57 | HS-0910 | 65 | HS-1571 | 75 |
| HS-0430 | 57 | HS-0920 | 65 | *HS-1580 | 75 |
| HS-0440 | 57 | HS-0930 | 65 | *HS-1590 | 76 |
| HS-0450 | 57 | HS-0940 | 65 | HS-1600 | 76 |
| HS-0460 | 58 | HS-0950 | 65 | HS-1610 | 76 |
| HS-0470 | 58 | HS-0960 | 65 | HS-1620 | 77 |
| HS-0480 | 58 | HS-0970 | 66 | *HS-1630 | 77 |
| HS-0490 | 58 | HS-0980 | 66 | HS-1640 | 78 |
| HS-0500 | 58 | HS-0990 | 66 | HS-1650 | 78 |
| HS-0510 | 58 | HS-1000 | 66 | HS-1660 | 78 |
| HS-0520 | 58 | HS-1010 | 66 | HS-1663 | 78 |
| HS-0521 | 58 | HS-1020 | 66 | HS-1664 | 78 |
| HS-0530 | 59 | HS-1022 | 66 | HS-1670 | 78 |
| HS-0540 | 59 | HS-1030 | 67 | HS-1690 | 79 |
| HS-0550 | 59 | HS-1040 | 67 | *HS-1700 | 79 |
| HS-0560 | 59 | HS-1050 | 67 | HS-1710 | 80 |
| HS-0570 | 59 | HS-1051 | 67 | HS-1720 | 80 |
| HS-0571 | 59 | HS-1053 | 67 | HS-1730 | 80 |
| HS-0580 | 60 | HS-1060 | 68 | HS-1740 | 80 |
| HS-0581 | 60 | HS-1090 | 68 | HS-1750 | 80 |
| HS-0590 | 60 | HS-1100 | 68 | HS-1760 | 80 |
| HS-0600 | 60 | HS-1110 | 68 | HS-1770 | 80 |
| HS-0610 | 60 | HS-1120 | 69 | HS-1780 | 80 |
| HS-0620 | 60 | HS-1130 | 69 | HS-1790 | 81 |
| HS-0621 | 60 | HS-1140 | 69 | HS-1800 | 81 |
| HS-0630 | 60 | HS-1150 | 69 | HS-1810 | 81 |
| HS-0640 | 60 | HS-1160 | 69 | HS-1820 | 81 |
| HS-0650 | 61 | HS-1170 | 70 | HS-1830 | 81 |
| HS-0660 | 61 | HS-1180 | 70 | HS-1840 | 81 |
| HS-0670 | 61 | HS-1250 | 70 | HS-1850 | 81 |
| HS-0680 | 61 | HS-1260 | 71 | HS-1860 | 81 |
| HS-0690 | 61 | HS-1270 | 71 | HS-1870 | 82 |
| HS-0691 | 61 | HS-1280 | 71 | HS-1880 | 82 |
| HS-0700 | 61 | HS-1290 | 72 | HS-1890 | 82 |
| HS-0710 | 61 | HS-1300 | 72 | HS-1900 | 82 |
| HS-0720 | 62 | HS-1310 | 72 | HS-1910 | 82 |
| HS-0730 | 62 | HS-1320 | 72 | HS-1920 | 82 |
| HS-0740 | 62 | HS-1330 | 72 | HS-1930 | 82 |
| HS-0750 | 62 | HS-1340 | 72 | HS-1940 | 82 |
| HS-0760 | 62 | HS-1350 | 72 | HS-1950 | 83 |
| HS-0770 | 63 | HS-1360 | 72 | HS-1960 | 83 |
| HS-0780 | 63 | HS-1370 | 72 | HS-1970 | 83 |
| HS-0790 | 63 | HS-1460 | 73 | HS-1980 | 83 |
| HS-0800 | 63 | HS-1470 | 73 | HS-1990 | 83 |
| HS-0810 | 63 | HS-1480 | 73 | HS-2000 | 84 |

| | | | | | |
|----------|-----|---------|-----|------------|-----|
| HS-2010 | 84 | L-0270 | 201 | L-0840 | 213 |
| HS-2020 | 84 | L-0280 | 201 | L-0850 | 214 |
| HS-2021 | 84 | L-0290 | 201 | L-0860 | 214 |
| HS-2030 | 85 | L-0300 | 202 | L-0870 | 214 |
| HS-2040 | 85 | L-0310 | 202 | L-0880 | 214 |
| HS-2050 | 85 | L-0320 | 202 | L-0890 | 214 |
| HS-2060 | 85 | L-0330 | 203 | L-0900 | 215 |
| HS-2070 | 85 | L-0340 | 203 | L-0910 | 215 |
| HS-2080 | 85 | L-0350 | 203 | L-0920 | 215 |
| *HS-2090 | 86 | L-0360 | 203 | L-0930 | 215 |
| HS-2100 | 86 | L-0370 | 203 | L-0940 | 215 |
| HS-2110 | 86 | L-0380 | 203 | L-0950 | 215 |
| HS-2120 | 87 | L-0390 | 203 | L-0960 | 216 |
| HS-2130 | 87 | *L-0400 | 204 | L-0970 | 216 |
| | | *L-0410 | 204 | L-0980 | 216 |
| | | *L-0420 | 204 | L-0990 | 216 |
| | | *L-0430 | 204 | L-1000 | 217 |
| | | *L-0440 | 205 | L-1001 | 217 |
| | | *L-0450 | 205 | L-1010 | 217 |
| | | L-0460 | 205 | L-1020 | 217 |
| I-0010 | 241 | L-0470 | 206 | L-1021 | 217 |
| I-0020 | 241 | L-0480 | 206 | L-1030 | 217 |
| I-0030 | 241 | L-0490 | 206 | L-1040 | 217 |
| I-0040 | 241 | L-0500 | 206 | L-1041 | 218 |
| | | L-0520 | 207 | L-1050 | 218 |
| | | L-0530 | 207 | L-1060 | 218 |
| | | L-0540 | 207 | L-1070 | 218 |
| | | L-0550 | 207 | L-1075 | 219 |
| | | L-0560 | 208 | L-1080 | 219 |
| | | L-0570 | 208 | L-1090 | 219 |
| | | L-0580 | 208 | L-1100 | 219 |
| L-0010 | 195 | L-0590 | 208 | L-1110 | 220 |
| L-0020 | 195 | L-0600 | 208 | L-1120 | 220 |
| L-0030 | 195 | L-0610 | 208 | L-1130 | 220 |
| L-0040 | 195 | L-0620 | 209 | L-1140 | 220 |
| L-0050 | 196 | L-0630 | 209 | L-1150 | 221 |
| L-0060 | 196 | L-0640 | 209 | L-1160 | 221 |
| L-0070 | 196 | L-0650 | 209 | L-1170 | 221 |
| L-0080 | 196 | L-0660 | 210 | L-1180 | 221 |
| L-0090 | 197 | L-0670 | 210 | L-1200 | 221 |
| L-0100 | 197 | L-0680 | 210 | L-1210 | 221 |
| L-0110 | 197 | L-0690 | 210 | L-1220 | 222 |
| L-0120 | 197 | L-0700 | 211 | L-1230 | 222 |
| L-0130 | 198 | L-0710 | 211 | L-1240 | 222 |
| L-0140 | 198 | L-0720 | 211 | L-1250 | 222 |
| L-0145 | 198 | L-0730 | 211 | L-1260 | 222 |
| L-0150 | 198 | L-0740 | 211 | L-1270 | 223 |
| L-0160 | 198 | L-0750 | 211 | L-1280 | 223 |
| L-0170 | 199 | L-0760 | 212 | L-1290 | 223 |
| L-0180 | 199 | L-0780 | 212 | L-1300 | 224 |
| L-0200 | 200 | L-0790 | 212 | L-1300-A | 227 |
| L-0210 | 200 | L-0800 | 212 | L-1300-I | 227 |
| L-0220 | 200 | L-0810 | 213 | L-1304 | 224 |
| L-0230 | 200 | L-0820 | 213 | L-1311 | 225 |
| L-0250 | 201 | L-0830 | 213 | L-1311-TIP | 226 |
| L-0260 | 201 | | | | |

| | |
|------------|-----|
| L-1321 | 225 |
| L-1321-TIP | 226 |
| L-1331 | 225 |
| L-1331-TIP | 226 |
| L-1341 | 225 |
| L-1341-TIP | 226 |
| L-1351 | 225 |
| L-1351-TIP | 226 |
| L-1361 | 225 |
| L-1361-TIP | 226 |
| L-1371 | 225 |
| L-1371-TIP | 226 |
| L-1381 | 225 |
| L-1381-TIP | 226 |
| L-1400 | 228 |
| L-1410 | 228 |
| L-1420 | 228 |
| L-1430 | 229 |
| L-1440 | 229 |
| L-1450 | 229 |
| L-1460 | 229 |
| L-1470 | 229 |
| L-1480 | 229 |
| L-1700 | 230 |
| L-1710 | 230 |
| L-1720 | 230 |
| L-1730 | 230 |
| L-1740 | 230 |
| L-1750 | 230 |
| L-1760 | 230 |
| L-1770 | 231 |
| L-1780 | 231 |
| L-1790 | 231 |
| L-1800 | 231 |
| L-1820 | 232 |
| L-1830 | 232 |

M

| | |
|--------|-----|
| M-0010 | 235 |
| M-0020 | 235 |
| M-0030 | 235 |
| M-0040 | 235 |
| M-0050 | 236 |
| M-0060 | 236 |
| M-0070 | 236 |
| M-0080 | 236 |
| M-0090 | 236 |

MIS

| | |
|-----------|-----|
| MIS-0010 | 239 |
| MIS-0011 | 237 |
| MIS-0020 | 237 |
| MIS-0020M | 237 |
| MIS-0030 | 237 |
| MIS-0142 | 240 |
| MIS-0143 | 240 |
| MIS-0150 | 240 |
| MIS-0160 | 242 |
| MIS-0170 | 242 |
| MIS-0180 | 242 |
| MIS-0190 | 242 |
| MIS-0200 | 242 |
| MIS-0210 | 242 |
| MIS-0220 | 242 |
| MIS-0230 | 243 |
| MIS-0240 | 243 |
| MIS-0250 | 243 |
| MIS-0260 | 243 |
| MIS-0270 | 243 |
| MIS-0280 | 243 |
| MIS-0290 | 243 |
| MIS-0300 | 244 |
| MIS-0310 | 244 |
| MIS-0320 | 244 |
| MIS-0330 | 244 |
| MIS-0340 | 244 |
| MIS-0350 | 245 |
| MIS-0360 | 245 |
| MIS-0370 | 245 |
| MIS-0380 | 245 |
| MIS-0390 | 245 |
| MIS-0400 | 245 |
| MIS-0410 | 245 |
| MIS-0420 | 245 |
| MIS-0430 | 246 |
| MIS-0440 | 246 |
| MIS-0450 | 246 |
| MIS-0460 | 247 |
| MIS-0470 | 247 |
| MIS-0480 | 247 |
| MIS-0490 | 247 |
| MIS-0500 | 248 |
| MIS-0510 | 248 |
| MIS-0520 | 249 |
| MIS-0530 | 249 |
| MIS-0540 | 249 |
| MIS-0550 | 250 |
| MIS-0560 | 250 |
| MIS-0570 | 250 |
| MIS-0580 | 250 |

| | |
|----------|-----|
| MIS-0611 | 252 |
| MIS-0612 | 252 |
| MIS-0613 | 252 |
| MIS-0614 | 251 |
| MIS-0615 | 251 |
| MIS-0770 | 239 |

N

| | |
|--------|-----|
| N-0010 | 157 |
| N-0020 | 157 |
| N-0030 | 157 |
| N-0040 | 157 |
| N-0050 | 157 |
| N-0060 | 157 |
| N-0070 | 157 |
| N-0080 | 157 |
| N-0090 | 158 |
| N-0100 | 158 |
| N-0110 | 158 |
| N-0120 | 158 |
| N-0130 | 158 |
| N-0140 | 158 |
| N-0150 | 158 |
| N-0160 | 158 |
| N-0250 | 159 |
| N-0260 | 159 |
| N-0270 | 159 |
| N-0280 | 160 |
| N-0290 | 160 |
| N-0300 | 160 |
| N-0310 | 160 |
| N-0320 | 160 |
| N-0330 | 160 |
| N-0340 | 161 |
| N-0350 | 161 |
| N-0360 | 161 |
| N-0370 | 161 |
| N-0380 | 161 |
| N-0390 | 161 |
| N-0400 | 162 |
| N-0410 | 162 |
| N-0420 | 162 |
| N-0430 | 162 |
| N-0440 | 162 |
| N-0450 | 162 |
| N-0460 | 163 |
| N-0470 | 163 |
| N-0480 | 163 |
| N-0490 | 163 |
| N-0500 | 164 |
| N-0510 | 164 |

| | | | | | |
|----------|-----|----------|-----|----------------|-----|
| N-0520 | 164 | S-0210 | 174 | S-0690 | 183 |
| N-0530 | 164 | S-0220 | 174 | S-0700 | 183 |
| N-0540 | 165 | S-0220SC | 187 | S-0710 | 183 |
| N-0550 | 165 | S-0230 | 174 | S-0720 | 183 |
| N-0560 | 165 | S-0240 | 174 | S-0730 | 184 |
| N-0570 | 165 | S-0250 | 174 | S-0740 | 184 |
| N-0620 | 165 | S-0250SC | 187 | S-0750 | 184 |
| N-0630 | 165 | S-0260 | 175 | S-0760 | 184 |
| N-0640 | 165 | S-0270 | 175 | S-0770 | 185 |
| N-0650 | 165 | S-0280 | 175 | S-0780 | 185 |
| N-0660 | 166 | S-0290 | 175 | S-0780SC | 186 |
| N-0670 | 166 | S-0300 | 175 | S-0790 | 185 |
| N-0680 | 166 | S-0310 | 175 | | |
| *N-0690 | 167 | S-0380 | 176 | | |
| N-0700 | 167 | S-0382 | 174 | | |
| N-0800 | 154 | S-0390 | 176 | | |
| N-0800TC | 156 | S-0400 | 176 | | |
| N-0810 | 154 | S-0410 | 176 | | |
| N-0810TC | 156 | S-0415 | 179 | | |
| N-0820 | 154 | S-0416 | 179 | VR-0010-23T | 189 |
| N-0820TC | 156 | S-0420 | 177 | VR-0030-23T | 189 |
| N-0830 | 154 | S-0430 | 177 | VR-0050-23T | 189 |
| N-0830TC | 156 | S-0440 | 177 | VR-0090-23T | 189 |
| N-0840 | 155 | S-0450 | 177 | VR-0110-23T | 189 |
| N-0850 | 155 | S-0460 | 177 | VR-0170-23T | 189 |
| N-0860 | 155 | S-0470 | 177 | VR-0205-23T | 189 |
| N-0870 | 155 | S-0480 | 177 | VR-0220 -23T | 190 |
| N-0880 | 155 | S-0490 | 178 | VR-0240-23T | 190 |
| N-0890 | 155 | S-0500 | 178 | VR-0250-23T | 190 |
| | | S-0501 | 178 | VR-0270 | 190 |
| | | S-0505 | 178 | VR-0280 | 190 |
| | | S-0510 | 178 | VR-0300-23TDT | 192 |
| | | S-0520 | 178 | VR-0300-23-TIP | 192 |
| | | S-0530 | 179 | VR-0310-23TDT | 192 |
| | | S-0540 | 179 | VR-0310-23-TIP | 192 |
| S-0010 | 170 | S-0550 | 179 | VR-0400-23TDT | 192 |
| S-0020 | 170 | S-0560 | 179 | VR-0400-23-TIP | 192 |
| S-0030 | 170 | S-0570 | 180 | VR-0410-23TDT | 192 |
| S-0040 | 170 | S-0580 | 180 | VR-0410-23-TIP | 192 |
| S-0050 | 170 | S-0581 | 180 | VR-1600 | 191 |
| S-0060 | 170 | S-0590 | 180 | VR-1610 | 191 |
| S-0070 | 171 | S-0600 | 180 | VR-1620 | 191 |
| S-0080 | 171 | S-0610 | 180 | VR-5000T | 192 |
| S-0090 | 171 | S-0611 | 180 | | |
| S-0100 | 172 | S-0620 | 180 | | |
| S-0110 | 172 | S-0621 | 181 | | |
| S-0120 | 172 | S-0630 | 181 | | |
| S-0130 | 172 | S-0640 | 181 | | |
| S-0140 | 172 | S-0641 | 181 | | |
| S-0150 | 172 | S-0650 | 181 | | |
| S-0160 | 173 | S-0660 | 182 | | |
| S-0161 | 173 | S-0670 | 182 | | |
| S-0170 | 173 | S-0670SC | 186 | | |
| S-0180 | 173 | S-0680 | 182 | | |
| S-0190 | 173 | S-0680SC | 186 | | |
| S-0200 | 174 | | | | |

VR

S

ALPHABETICAL INDEX

A

Absorbent Strips 32
Adjustable Eye Speculum 18
Anaesthesia Cannula 220
Anis Needle Holder 163
Anterior Chamber Irrigating Cannula 215, 216
Anterior Chamber Irrigator 214
Anterior Chamber Maintaining Cannula 216
Aspirating Cannula 207, 208, 221, 222, 228, 229
Aspiration Handpiece 225, 227
Aspiration Handpiece only 227
Atkinson Retrobulbar Needle 220
Axenfeld Hook 29
Axis Marker 75
Ayer Chalazion Forceps 138

B

Baby Jones 149
Backhaus towel clamp 150
Baird Chalazion 140
Baird Chalazion Forceps 140
Bangerter Muscle Clamp 144
Bard-Parker Handle 46, 47
Barraquer Cilia Forceps 134, 135
Barraquer Fixation Forceps 134
Barraquer Iris Scissors 179
Barraquer Iris Spatula 54, 58, 63
Barraquer Needle Holder 154, 155, 156, 161, 162
Barraquer Wire Speculum 11, 12, 13
Barrett Eye Speculum 11
Bechert Lens Holding Forceps 146
Bechert Nucleus Rotator 52, 55, 60, 67
Berke Ptosis Forceps 141
Beyer Rongeur 196
Bimanual I/A 6, 194
Bipolar Cord 150, 151
Bipolar Forceps 6, 150, 151
Bishop-Harmon Anterior Chamber Irrigator 214
Bishop Harmon Fixation Forceps 128
Bishop-Harmon Tissue Forceps 129
Blade Breaker & Holder 166, 167
Blaydes Implantation Forceps 146
Bonaccolto Utility Forceps 130

Bonn Forceps 92, 93, 104, 113, 123, 124
Bores Optic Zone 81, 82, 83
Bores Optic Zone Marker 81, 82, 83
Braunstein Measuring Instrument 236
Brauweiler Aspirating Handpiece 225
Breaker and Holder 166
Brown Triple 20
Bunge 73
Burrs 6, 34

C

Caliper 235, 236
Cannula 108, 207, 208, 209, 210, 211, 212, 213, 214, 215, 216, 218, 219, 220, 221, 222, 223, 224, 228, 229
Cannulas 6, 194
Capsule Polisher 217, 218
Capsulorhexis Forceps 92, 103, 104, 105, 112, 113, 115, 116, 117, 118, 122, 123, 124, 125, 126, 147
Capsulotomy Scissors 177, 178
Carter Sphere Introducer 239
Case 243, 245, 246
Castroviejo Breaker and Holder 166
Castroviejo Caliper 235
Castroviejo Corneoscleral Scissors 172
Castroviejo Eye Speculum 21
Castroviejo Keratoplasty Scissors 173
Castroviejo Lacrimal Dilator 197
Castroviejo Needle Holder 158
Castroviejo Spatula 68
Castroviejo Suturing Forceps 106, 107, 109, 119
Castroviejo Tying Forceps 105
Castroviejo Universal Corneal Scissors 170, 173
Cataract Knife 37
Chalazion Forceps 138, 139, 140, 141
Chalazions Curette 74
Chang Phaco Chopper 55, 68
Chisel 200, 201, 202
Chisels 6, 10
Chuck Handle 47
Cilia Forceps 134, 135
Citelli Bone Nibbling Rongeur 196
Clayman Lens Holding Forceps 146, 148
Colibri Forceps 94, 95, 96, 107, 110, 125, 147
Combo Pre-Chopper 127
Conjunctiva Fixation Forceps 105, 118, 130
Conjunctiva Forceps 113, 131, 137

Conjunctival Scissors 174
Cook Eye Speculum 16
Corneal Dissector 36, 66
Corneal Fixation Forceps 95, 99, 100
Corneal Knife 35, 38
Corneal Rust Ring Remover 40
Corneal Scarifier 36
Corneal Scarifiers 34
Corneal Transplant Scissors 176
Corneal Utility Forceps 94, 95, 100
Corneoscleral Punch 48
Corneoscleral Scissors 172
Cranz Eye Douche 219
Crawford Retrieving Hook 205
Crescent Diamond Knife 42
Cyclodialysis Cannula 223

D

DALK Cannula 224
DALK Corneal Dissector 66
Depressor and Fixation 33
Derf Needle Holder 167
Desmarres Chalazion Forceps 139
Desmarres cross-action Forceps 141
Desmarres Lid Retractor 24
DeWecker Iridectomy Scissors 178
Diamond step knife 41
Dieffenbach 145
Direct Action Forceps 92, 112
Discission Knife 37
Discission Needle 39
Dissector 34, 36, 66
Double Ended Strabismus Hook 32
double hook 55
Double-level tray 251
Downs Lid Forceps 142
Dressing Forceps 128, 132, 133
Drysdale Nucleus Manipulator 52, 60

E

Elevator 198, 199
Elschnig Skleral Fixation Forceps 130
Endgripping Forceps 189
Endothelial / Descemetorhexis Forceps 114
Entropium Forceps 142
Enucleation Scissors 184
Enucleation Spoon 71
Ergo Grip 91
Ergo Grip Barraquer Needle Holder 154, 155, 156

Erhardt Chalazion Forceps 139
Eye Scissors 180, 181
Eye Shield 240
Eye Speculum 11, 13, 14, 16, 17, 18, 19, 20, 21

F

Fascia Needle 10, 25
Femto Flaplifter 67
Fine Stitch Scissors 171
Fishkind Eye Speculum 21
Fixation Forceps 95, 99, 100, 105, 107, 109, 110, 118, 128, 130, 134
Fixation Hook 28
Fixation Ring 22, 76, 77
Fixation Rings 6, 10
Flap Stretching Forceps 113
Flieringa Scleral Fixation Ring 22
Foreign Body Forceps 138, 190
Foreign Body Needle 38, 39
Foreign Body Spud 39
Fragment Forceps 91, 111, 122, 124
Francis Chalazion Forceps 140
Freer Lacrimal Chisel 202
Freer Periosteal Elevator 198

G

Gaskin Fragment Forceps 91, 111, 122, 124
Gass Retinal Detachment Hook 31
Gass Scleral Marker 23
GIANNETTI MIC Capsulorhexis Forceps 115, 117
Gills-Vannas Capsulotomy Scissors 177, 178
Gills-Welsh Capsulotomy Scissors 178
Gills Welsh Olive Tip Cannula 218
Ginsberg Eye Speculum 19
Globe Retractor 71
Gradle Cilia Forceps 135
Graefe model 134
Graefe Muscle Hook 30
Graether Collar Button Micro Iris Hook 63
Grasping Forceps 191
Green Muscle Hook 31

H

Haefliger Phaco cleaver 51, 56
Halsted Mosquito Forceps 149
Harms Trabeculotomy Probe 202

Hartmann Mosquito Forceps 149
Helveston Muscle Hook 30
Hokey Knife 35
Hooks 6, 10, 50
Hunt Chalazion Forceps 139
Hydrodissection Cannula 211, 212, 213

I

ICL Manipulator 58
ILM Forceps 192
Implantation Forceps 92, 112, 146, 147
Inamura Capsulorhexis Forceps 147
Incision Marker 77
Infant Lid Speculum 15
Infusions Plug 229
Injection Cannula 211
Instruments 6, 57, 90, 153, 194, 234, 252
Intraocular Pressure Reducer 237
IOL Cutter 6, 49
IOL Cutters 34
IOL Implantation Forceps 92, 112
IOL Inserting Forceps 120
IOL-Manipulator 52, 53, 61
Iridectomy Scissors 178
Iris Botvin Forceps 94
Iris Forceps 97, 98, 119, 132, 133
Iris Hook 51, 57, 58, 60, 63, 64, 65
Iris Knife 38
Iris Retractor 33
Iris Scissors 175, 179
Iris Spatula 54, 58, 63, 69
Irrigating-Aspirating Cannula 229
Irrigating Nucleus Spoon 217
Irrigating Vectis Cannula 219
Irrigation Handpiece 225, 227
Irrigation Handpiece only 227

J

Jaeger Lid Plate 25
Jaeger Strabismus Hook 32
Jaffe Capsulorhexis Forceps 92, 112, 115, 122
Jaffe Forceps 114
Jaffe Stitch Scissors 171
Jaffe Wire Lid Retractor 24
Jakobi Needle Holder 164
Jameson Caliper 235
Jameson Muscle Forceps 143
Jameson Muscle Hook 31
Jayle Forceps 131

Jensen Capsule Polisher 217
Jeweler's Forceps 136
Joseph Periosteal Elevator 199

K

Kalt Needle Holder 159
Kauter 86
Kelman Cannula 223
Kelman-McPherson Forceps 91, 110, 113
Keratometer 6, 87
Keratoplasty Scissors 173
Kerrison Bone Nibbling Rongeur 195
Knapp Strabismus Scissors 184
Knife Tray 244
Knives 6, 34
Knolle-Pearce Irrigating Vectis Cannula 219
Koch Aspirating Handpiece 225
Koch Stop & Chop Manipulator 58
Kraff Eye Seculum temporal 19
Kraff Lens Holding Forceps 146
Kratz-Barraquer Wire Speculum 13
Kratz Cystotome Cannula 222
Kuglen Iris Hook 51, 60, 65

L

Lacrimal Cannula 208, 209, 210
Lacrimal Dilator 197, 198
Lacrimal Instruments 6, 194
Lacrimal Intubations Set 204
Lacrimal Probe 203, 210
Lacrimal Probe Cannula 210
Lacrimal Sac Chisel 200
Lacrimal Sac Retractor 21, 27, 29
Lambert Chalazion Forceps 138
Lancaster Eye Speculum 18
LASEK Alcohol Well 84
LASIK Cannula 207, 220
LASIK Corneal Marker 79
LASIK Dulaney Corneal Marker 78
LASIK Flap Forceps 113
LASIK Flap Lifter 67
LASIK Flap Protector 84
LASIK Spatula 64
Lehner Capsulorhexis Forceps 92, 103, 104, 112, 117, 125, 126
Lens cutting Scissors 191
Lens Folding Forceps 92, 112
Lens Forceps 91, 110, 111, 144, 146
Lens Holding Forceps 91, 97, 111, 112, 122, 146, 148
Lens Injector 241

Lens Loop 72
Lens Manipulating Hook 53, 61, 62, 65, 66
Lens Manipulator 51, 54, 58, 60, 62, 65
Lens Spatula 59, 66
Lester Conjunctiva Fixation Forceps 130
Lester IOL-Manipulator 52, 53, 61
Lid Clamp 143
Lid Forceps 142
Lid Plate 25
Lid Retractor 24
Lid Speculum 15
Lieberman Aspirating Speculum 14
Lieberman Eye Speculum 14, 18, 19, 20
Ligabue Nucleus Manipulator & Splitter 54
Ligature Scissors 185
Lippmann Micro Cystotome 224
Lister Conjunctiva Forceps 131
Loop 72, 219
Lordan Chalazion Forceps 141
Luer-Lock Adaptor 232

M

Mallet 199
Manipulators 6, 50
Marking Instruments 234
Marking Pen 6, 87
Masket Phaco Spatula 55, 68
McIntyre-Binkhorst Cannula 213
McNeil-Goldman Scleral Fixation Ring 22
McPherson Needle Holder 157
McPherson-Westcott Stitch Scissors 170
Meibomian Gland Compressor 114
Mellinger Eye Speculum 17
Membrane Hook 85
Membrane Scratcher 85
Membrane Spatula 85
Membran-Peeler 85
Mendez Degree Gauge 75
Metal Bole 240
Metzenbaum Scissors 185, 186
Meyerhoefer Chalazions Curette 74
MIC Corydon Capsulorhexis Forceps 115, 116
MIC Jaffe Capsulorhexis Forceps 115
Micro Colibri Forceps 95
Micro Finger 64
Micro Gripping Forceps 192
MICS Capsule Polishing Cannula 219
MICS Diamond Knife 42
MIC Utrata Capsulorhexis Forceps 115
Miscellaneous 6, 238
Monomanual Phaco Safe-Lock I/A Handpiece 230
Moorfields Conjunctiva Forceps 137

Mosquito Forceps 149
Murdoch Eye Speculum 17
Muscle Clamp 144
Muscle Forceps 143
Muscle Hook 30, 31

N

Nagahara Phaco-Chopper 51, 52, 56
Nasal Speculum 197
Needle & Gouge 26
Needle Holder 152, 154, 155, 156, 157, 158, 159, 160, 161, 162, 163, 164, 167
Needle Holder with TC-Dust 156
Needles 6, 34
Neuhann Nucleus Rotator 57
Nevyas Capsulorhexis 92, 113
Nevyas Capsulorhexis Forceps 92, 113
Noyes Iris Scissors 175
Nucleus Cracker 108
Nucleus Cracking Forceps 121
Nucleus Divider 57
Nucleus Loop 72
Nucleus Manipulator 52, 54, 60, 64
Nucleus Manipulator & Splitter 54
Nucleus Rotator 52, 55, 57, 60, 66, 67
Nucleus Rotator & Dividing Hook 55, 67
Nucleus Rotator (push-pull) & Dividing Hook 55, 67
Nucleus Sustainer 64
Nugent Utility Forceps 133

O

Ocular Landers Lens Forceps 144
Oil Pen 239
Optic Zone Marker 81, 82, 83
Orbital Spatula 70
O-Ring Replacement Set 230
Osher Incision Caliper 236

P

Paper Filters 247
Paper Label/Indicator 248
Partsch Lacrimal Sac Chisel 200
Partsch Mallet 199
Paufigue Suture Forceps 101, 102
Phaco-Chopper 51, 52, 56, 62
Phaco Knife 43
Phaco Pre-Chopper 120
Phaco Splitter 51, 60, 68

Phaco Trays 243, 244
Pickup & Folding Forceps 121
Pigtail Probe 206
Pisacano Nucleus Rotator 52, 60
Plastic Safety Seals 248
Pressure Plate 86
Pterygium Drill 40
PTFE Teflon Filters 247
Ptosis Forceps 141
Punches 6, 34
push-pull 51, 55, 60, 65, 67
Putterman Forceps 143

R

Radial Marker 78
Remky Superior Rectus Forceps 131
Removal Cannula Forceps 108
Replacement burr 40
Replacement tip 51, 52, 53, 54, 55
Retractors 6, 10
Retrieving Hook 205
Retrograde Probe 204
ribbon type 181, 183
Rongeur 195, 196
Rosen Phaco Splitter 51, 60, 68
Rowen Phaco Spatula & Bechert Nucleus Rotator 55, 67
Ruler 236

S

Sac Gouge 201
Sac Retractor 21, 27, 29
Safe - Lock I/A 230, 231
Safe - Lock I/A tips 230
Sauer 15, 102
Schepens Scleral Depressor 23
Schocket Scleral Depressor 23
Scissors 6, 168, 169, 170, 171, 172, 173, 174, 175, 176, 177, 178, 179, 180, 181, 182, 183, 184, 185, 186, 187, 190, 191, 192
Scleral Depressors 6, 10
Scleral Knife 36
Scleral Plug Forceps 144
Serrifine 145
Shaaf Cilia Foreign Body Forceps 138
Shepard Lens Forceps 146

Silicone Finger Mats 250
Silicone Instrument Support 250
Silicone Irrigation Sleeve 231
Simcoe Cortex Extractor 206
Simcoe Irrigating - Aspirating Cannula 228
Sinsky II Lens Manipulating Hook 53, 61, 66
Skleral Fixation Forceps 130
Sleeve Spreader 137
Sleeve tip 231
Snapper Hook 64
Snellen Entropium Forceps 142
Spatula 50, 54, 55, 57, 58, 59, 63, 64, 66, 67, 68, 69, 70, 72, 85
Speculum 11, 12, 13, 14, 15, 16, 17, 18, 19, 20, 21, 197
Speculum adjustable 11
Spencer Ligature Scissors 185
Spoon 6, 71, 72, 217
Spud & Gouge 34
step knife 41
Sterilisation Case 245, 246
Sterilisation Case aluminium 246
Sterilisation cassette 251
Steri Oil Pen 239
Stevenson Lacrimal Sac Retractor 21
Stevens Tentotomy Scissors 182, 186
Stitch Scissors 171, 181, 182
Stop & Chop Manipulator 58
Strabismus Hook 32
Strabismus Scissors 183, 184
Super-Cut 186
Superior Rectus Forceps 131
Suturing Forceps 91, 98, 99, 101, 102, 106, 107, 109, 111, 119, 120, 125
Swiss Model 166
Synechiaie Spatula 66

T

„Teaser“ Hook 31
Temporal Approach 12
Tennant-Maumenee Colibri Forceps 147
Tennant Tying Forceps 91
Tenotomy Scissors 170, 174, 176, 183, 187
Textile Filters 247
Thornton Fixation Ring 76

Thorpe Conjunctiva Fixation Forceps 105, 130
Tip only 192
Tissue Forceps 129
Tissue Retractor 26
titanium 49, 51, 87, 90, 94, 98, 99, 100, 109, 110, 111, 112, 113, 122, 123, 147, 153, 162, 241
Tonometer 237
towel clamp 149, 150
Trabeculectomy Marker 84
Trabeculectomy Punch 48
Trabeculectomy Probe 202
Trepine Blade 46
Trepine Handle 46
Trepines 6, 34
Troutman-Chris Blade Breaker & Holder 167
Troutman-Katzin Corneal Transplant Scissors 176
Troutman Tying Forceps 120
Tying Forceps 91, 93, 98, 103, 104, 105, 108, 110, 118, 120, 122, 124, 125, 126

U

Universal Corneal Scissors 170, 173
Universal Trepine Handle 46
Uribe-Stern Capsulotomy Scissors 178
Urretts-Zavalía Depressor and Fixation 33
Utility Forceps 94, 95, 96, 100, 130, 133
Utility Scissors 185
Utrata Capsulorhexis 92, 103, 105, 112, 115, 118, 122, 123, 124, 125, 126
Utrata Capsulorhexis Forceps 92, 103, 105, 112, 115, 118, 122, 123, 124, 125, 126

V

Vannas Capsulotomy Scissors 177, 178
Vannas Scissors 172
Vera Fixation Hook 28
Viscoelastic Aspirating Cannula

221, 222

Vitrectomy 6, 86, 188, 189, 190

Vitrectomy Forceps 189, 190

Vitrectomy Foreign Body Forceps
190

Vitrectomy Scissors 190

W

Wadsworth Todd Kauter 86

Wash Tray 249

Watzke Sleeve Spreader 137

Westcott Tenotomy Scissors 170,
174, 176, 187

Wies Chalazion Forceps 140

Wilder Lacrimal Dilator 198

Williams Eye Speculum 18

Wire Speculum 11, 12, 13, 16

Y

Y-Hook 54, 58

Z

Zaldivar Implantation Forceps 147

Zone Marker 80, 81, 82, 83